Method/méthode 1 (drawings/dessins 1-10)

DE

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EN

Short stretch bandages for support, relief and compression.

Dressing instruction:
We recommend that bandaging be performed by a trained person.

Indications:

Compression bandages for phlebological and lymphological indications to support muscle pump action and improve the haemodynamics of the limb. For compression bandages on the back, e.g. in venous insufficiency, e.g. varicose ulcer and venous oedema, thrombophlebitis, phlebotrombosis, after vein operations (stripping), after sclerosing obliteration, primary and secondary lymphatic oedemas.

Contraindications:

Incompatibility to fabric, advanced peripheral arterial occlusive disease, decompensated cardiac insufficiency, septic phlebitis, phlegmasia coerulea dolens.

Cautions:

Wearing dermatosis; cutaneous infections. Open wounds must be covered with an appropriate dressing before compression bandaging. Severely compromised skin sensibility and impaired sensitivity of the limb; advanced peripheral neuropathy; conditions in which increased venous return is not desired. The attending physician decides whether compression in general can be applied to patients suffering from the following diseases: rheumatoid arthritis, complex regional pain syndrome, CRPS (M. Sudeck), malignant lymphedema, gangrene.

Longitudinal elasticity:

Compridur®: ca. 50%
JOBST® Comprihaf®: ca. 75%
Comprilan®: ca. 90%

Special notes:

Applies only to **JOBST® Comprihaf®**: initially unroll approximately 20-30 cm of the bandage before applying. This is the only way to guarantee a continuous drop in pressure from distal to proximal.

Compridur® and **Comprilan®**: The attached fixation clips should be not used for fixation while the bandage is applied on a patient. Due to the risk of unintended injuries both clips are solely designed for transportation purpose.

Any serious incident that has occurred in relation to this device should be reported to BSN medical Inc. (BSN medical GmbH within the European Union) and the competent authority of your state.

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JBR13976.05

4605018/44.04

	Artwork ID: JBR13876.05	Review Loop: 2	Printing Colours:  Pantone 281 C
Project Title: C2738_201 - JOBST Canada-Compliance BifG Comprilan	Date of Review Loop: 03.12.2020	Replaces Artwork ID: JBR13876.03	
Contact: ArtworkManagementMedical@essity.com	Technical Drawing for print ID: ---		
Contact Address: BSN medical GmbH Schützenstraße 1-3 22761 Hamburg Germany	Note LOT and expiry-date symbol will be part of variable coding.		
Languages: DEU, GBR, FRA, ESP, ITA, NLD, SWE, NOR, FIN, PRT, DNK			