Komprex®

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Product description and performance characteristics:

Komprex® Foam Rubber Pads: Ready-to-use

Komprex® Foam Rubber Bandages: Ready-to-use and can be cut to size

Komprex® Foam Sheets: can be cut to size

Product composition:

Komprex® Foam Rubber Pads: styrene-butadiene rubber

Komprex® Foam Rubber Bandages: 80% styrene-butadiene rubber

Komprex® Foam Sheets: 80% styrene-butadiene rubber

Intended purpose:

For padding under compression bandages

Indications:

Phlebological and lymphological indications: Varicosis, deep venous insufficiency, thrombophlebitis, phlebothrombosis, post-thrombotic syndrome, Widmer’s stage I to III b chronic venous insufficiency, venous and mixed ulcers, oedema (phleboedema, lymphoedema, lipoedema, post-traumatic, post-operative, during pregnancy), after sclerotherapy and venous surgery. Contusions, sprains, tendon injuries and joint dislocations.

Contraindications:

Advanced peripheral arterial occlusive diseases, congestive heart failure, septic phlebitis, phlegmasia cerulea dolens, sensation disorders of the skin, known allergies and/or hypersensitivity to any of the product components.

Side effects:

In (very) rare cases, skin irritations and/or allergies may occur.

Warnings and precautions:

- Incorrect bandaging can result in skin necrosis and pressure damage to nerves.
- Before applying the bandage, clean/debride any wound that is present and cover with a wound dressing appropriate for the wound status.
- Not intended to be sterilised.

Usage information:

A. Komprex® Foam Rubber Pad:

- a) Wrap Komprex® Pad with gauze
- b) Wound pad
- c) Ulcer

- Apply the rounded side in direction to the body without tension
- Apply the compression bandage using an appropriate wrapping technique

B. Komprex® Foam Rubber Bandage:

- Wrap the limb with a bandage liner

Wrap Komprex® Bandage without pulling tight or overlapping or cut to size

- Apply the compression bandage using an appropriate wrapping technique

C. Komprex® Foam Rubber Sheet:

- Cut out the desired shape
- Angle the edges
- For the remaining application steps, see Komprex® Foam Rubber Pad

Disposal:

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Packaging labelled as recyclable should be taken into the respective national recycling systems.

General instructions:

In case of a serious incident, contact the manufacturer and the responsible health authorities.

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