NÀVOD NA POUŽITIE
使用说明
ΥΠΟΔΕΙΞΕΙΣ ΧΡΗΣΗΣ
ANVISNINGER
GEBRAUCHSANWEISUNG
INSTRUCCIONES DE USO
KÄYTTÖOHJEET
INSTRUCTIONS
FOR USE
GEBRUIKSAANWIJZING
BRUKSANVISNING
ISTRUZIONI PER L'USO
INSTRUÇÕES

Gebrauchsanweisung	2	Nàvod na použitie	11
Instructions for use	3	Napomene za upotrebu	12
Istruzioni per l'uso	4	Kullanım bilgileri	13
Mode d'emploi	5	使用说明	14
Brugsanvisning	6		
Gebbruksaanwijzing	7		
Указания по применению	8		
Wskazówki dotyczące użycia	9		
Upute za upotrebu	10		

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Instructions for use:**Putting on • Taking off • Care****Intended purpose:**

Two-layer compression treatment on the leg to support venous reflux and the associated faster wound healing in *venous ulcers*.

Important information

- Do not wear ULCER X on open wounds. The wounds should be covered with the appropriate wound dressings.
- Always wear the special SIGVARIS rubber gloves when putting on or taking off the stockings. The rubber gloves enable you to don the stockings with relatively little effort and thus preserve the compression stocking fabric.
- The understocking has to be worn for 24 hours. Put the overstocking on after getting up and remove it before going to bed.
- SIGVARIS ULCER X should only be worn as prescribed by your doctor.
- Do not remove the excess threads inside the stocking.
- Protect your TRADITIONAL compression stocking from ointments, oils and grease (including body lotion), from direct sunlight, bright light and considerable heat. The understocking is largely insensitive towards ointments, oils, grease and UV rays.

Putting on the understocking

1. Turn the understocking inside out at the level of the ankle and pull the understocking carefully over the instep and heel like a normal sock. This inversion technique ensures that the wound dressing remains in place.
2. Take hold of the understocking by the cuff and carefully pull it half way up the lower leg.
3. Distribute the stocking fabric evenly over the lower leg by gradually easing the understocking up the leg with flat hands and stroking movements. This allows you to distribute the knitted fabric evenly over the lower leg with very little effort.

Putting on the TRADITIONAL compression stocking

4. Turn the TRADITIONAL compression stocking inside out at the level of the ankle and pull the stocking carefully over the instep and heel like a normal sock.
5. Take hold of the TRADITIONAL compression stocking by the cuff and carefully pull it half way up the lower leg.
6. Distribute the stocking fabric evenly over the lower leg by gradually easing the TRADITIONAL compression stocking up the leg using flat hands and stroking movements.

Taking off the TRADITIONAL compression stocking

7. To doff the TRADITIONAL compression stocking, grip the upper edge and pull the band down over the heel and almost to the toe.
8. Using flat hands and stroking movements, ease the double-layered stocking over the instep and heel, and off towards the toe.

Correct care

- Wash daily by hand or in a laundry bag in a washing machine (delicate wash cycle) at max. 40 degrees (understocking: max. 60 degrees). 
- For hand washing, we recommend SIGVARIS WASHING SOLUTION hand washing solution (do not use fabric softener). 
- Do not tumble-dry. 
- Do not iron. 
- Do not dry-clean. 
- Do not use any chemicals.

Indications:

For the treatment of leg ulcers (*ulcus cruris venosum*) subsequent to wound treatment.

Relative contra-indications:

Pronounced weeping dermatoses, intolerance to compression stocking material, serious sensitivity disorders of the extremities, advanced peripheral neuropathy (e. g. diabetes mellitus), primary chronic polyarthritis.

Absolute contra-indications:

Advanced peripheral arterial occlusive disease, congestive heart failure, septic phlebitis, phlegmasia coerulea dolens.

This quality product by SIGVARIS contains high quality refined natural rubber (latex) that does not come into direct contact with the skin.

The medical benefit of SIGVARIS compression stockings is guaranteed for 6 months with correct use.

All serious incidents that occur when using the product should be reported to the manufacturer and the competent supervisory authority.

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

5

6

Anziehen Kompressionsstrumpf (TRADITIONAL) / Putting on your compression stocking (TRADITIONAL) / Enfiler le bas de compression (TRADITIONAL) / Come indossare la calza compressiva (TRADITIONAL)

Reference:
ULCER X

7

8

Ausziehen Kompressionsstrumpf (TRADITIONAL) / Taking off your compression stocking (TRADITIONAL) / Ôter le bas de compression (TRADITIONAL) / Come sfilare la calza compressiva (TRADITIONAL)

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SIGVARIS AG (DMCC Branch), Unit No: AG-05-C, Silver (AG) Tower, Jumeirah Lakes Towers, Dubai - United Arab Emirates
Phone +971 4 558 6265

SIGVARIS INC., 1119 Highway 74 South - Peachtree City, GA 30269/USA
Phone +1 770 631 1778, certified ISO 9001

SIGVARIS CORP., 4535 Dobrin road, H4R 2L8, Quebec/CDN
Phone +1 514 336 2362

SIGVARIS DO BRASIL, Av. José Benassi, 2323, Distrito Industrial Fazgran, 13213-085-Jundiá, SP/Brazil
Tel. +55 114525 5700

If you look for further distributors please consult: www.sigvaris.com

 SIGVARIS AG, Gröblistrasse 8, 9014 St. Gallen/Switzerland
Made in Switzerland

SIGVARIS and SIGVARIS GROUP are registered trademarks of SIGVARIS AG, CH-9014 St. Gallen/Switzerland, in many countries worldwide. © 2019 Copyright by SIGVARIS AG, CH-9014 St. Gallen/Switzerland

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SIGVARIS AG, Gröblistrasse 8, 9014 St. Gallen, Switzerland
SIGVARIS GmbH, Dr.-Karl-Lenz-Straße 35, 87700 Memmingen, Germany
SIGVARIS S.r.l., Via Galileo Galilei 2/2, 31057 Silea (TV), Italy
SIGVARIS SAS, z.l. Sud d'Andrézieux, Rue B. Thimonier - BP 60223, 42173 Saint-Just Saint Rambert, France
SIGVARIS BRITAIN LTD., 1 Imperial Court, Magellan Close, Andover, SP10 5NT, Great Britain
SIGVARIS AG, ul. Powidzka 50, Gutowo Małe, 62-300 Września, Poland