DE

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**Produktbeschreibung**

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**Instructions for use**

**Product description**

Hydrocoll is a sterile self-adhesive wound dressing consisting of an adhesive hydrocollid layer covered with a polyurethane film. The wound side is covered with a release liner composed of siliconized paper. The primary packaging is a peel pouch made of film and coated paper. Hydrocoll is a non-active medical device. It does not contain any medicinal substances, tissue or cells of human or animal origin or nanomaterials. It is sterilized using irradiation.

**Composition**

Hydrocoll consists of a gas permeable polyurethane cover. The hydrocollid compound itself mainly contains sodium carboxymethylcellulose, a synthetic rubber mixture based on polystyrene polyisoprene copolymers and polyisobutylene.

**Properties and mode of action**

The product absorbs wound exudate. After the exudate has been absorbed, the hydrocollid layer transforms into a gel which provides a moist wound healing environment. The ingress of germs is prevented by the polyurethane film. The film is water-repellent and allows moisture vapor and gas to pass through. Since the gel is produced during the application, the dressing does not adhere to the wound and the healing process is not impaired by the dressing change. Newly formed tissue is not damaged.

Due to the thin flexible material, Hydrocoll molds easily to body contours. The surface of Hydrocoll is dirt-repellent and water-repellent.

**Intended purpose**

Hydrocoll is a sterile, single-use wound dressing for the treatment of medium to lightly exuding wounds such as chronic wounds (leg ulcers, decubitus, etc.) or burns up to grade II, used by healthcare professionals. Accumulated use of the product can exceed 30 days.

Hydrocoll thin is intended for the treatment of lightly exuding wounds.

Hydrocoll is intended for the treatment of human skin by professional users.

**Mode of application**

Select the correct size and shape of Hydrocoll. The dressing should exceed the wound by at least one centimeter to ensure that it stays firmly in place. If needed, cut the dressing to the size of the wound with sterile scissors.

Hydrocoll / Hydrocoll thin: remove the protective paper backing, place Hydrocoll on the wound, adjust the edges to fit and press on.

Hydrocoll sacral: Schutzpapier von der Rückseite entfernen und press the dressing into the anal cleft. The tip of the dressing can point upwards or downwards depending on the position and depth of the ulcer. Press the edges on firmly.

Hydrocoll concave for the heel and elbow: remove the largest part of the protective paper backing, place the middle part of the Hydrocoll dressing over the wound, remove the rest of the protective paper backing, then overlap the wings to fit.

Due to the absorption of exudate a whitish gel is formed that is visible on the top as discoloration and a blister. Hydrocoll should be changed if the blister has almost reached the edge of the dressing. Should exudate leak from the dressing, it must be changed immediately. To change the dressing, lift Hydrocoll at the edges and remove the whole dressing carefully.

It may happen that a layer of hydrocollid gel remains on the wound. It will not impair the healing of the wound and can be rinsed off if necessary.

The wear time of the dressing depends on the condition of the wound and the level of exudate but can remain on the wound for up to 7 days.

**Contraindications**

Hydrocoll is contraindicated in any patients who may be allergic to any of its ingredients. Hydrocoll must not be used on wounds extending into muscle, tendon or bone. It must not be applied to clinically infected wounds or third-degree burns.

**Product disposal**

In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

**Special precautions**

In the absence of available data supporting the use of this dressing on sensitive population groups such as infants, children, pregnant or nursing women, and in the absence of data to the contrary on these population groups, this dressing should be used with caution following a clinician's recommendation.

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

**Incident reporting**

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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