INSTRUCTIONS FOR USE

**BLADDER CATHETERIZATION SYRINGES AND KIT**

**Instruction for use**

English

**Intended use**

Sterile patient coverage and accessories for bladder catheterization procedures.

**Warning**

Sterile. Single use. Before use, check the integrity of the package. Discard any product that shows signs of damage, this could compromise the sterility of the product.

This product must be exclusively used by qualified personnel. Do not reuse. Do not re-sterilize.

It is imperative that the user and/or patient report any serious incidents occurring in relation to the device, to the manufacturer and to the competent authority of the state in which the user and/or patient is established.

**Prefilled syringes: ONLY FOR INDICATED USE (see below); NOT FORMULATED EITHER RECONSTITUTION OR DILUTION OF DRUGS, OR FOR INTRAVENOUS ADMINISTRATION.**

Do not use if the syringe cap is damaged. Do not use if the syringe contents contain particulate matter or precipitates. Do not use if solution is cloudy or colored.

**Indicated uses of pre-filled syringes:**

-10 ml syringe of water or water and glycerin: filling the bladder catheter anchor balloon

-10 ml syringe LS of saline: filling the anchor balloon of the bladder catheter

-10 ml syringe of lubricating gel: lubrication of the bladder catheter

-30 ml syringe of physiological solution: washing of the external genitalia

**Effects or Interactions**

The device contains a substance that is introduced into the human body through an orifice in the body; laboratory tests have been conducted to evaluate its bioavailability effects. These studies demonstrate no interactions or effects in the human body.

A Biological Rationale (BEP) has been drawn up for the physiological solution in which the study confirms the absence of interactions or effects in the human body, as this product is completely equivalent to those commonly used in hospitals.

**Lubrificant Gel Composition**

- Water 81,9 %
- Propylene glycol 8 %
- Propanediol 5 %
- Hydroxyethylcellulose 4 %
- Phenoxyethanol 1 %
- Ethylhexylglycerin 0,1

**Composition Saline Solution**

Sterile saline solution NaCl 0,9 %

100 ml of solution contains:

- Sodium chloride 0.9 g
- Water for injectable preparations q.b. to ml 100

**Instructions for use**

1. Open the sterile package following the correct "PEEL OPEN" method and arrange the material in the correct order inside the sterile field using the aseptic technique.
2. Perform the disinfection of the genital area following the operating procedures provided by the hospital protocol;
3. Perform the bladder catheter insertion technique following the operating procedures outlined by the hospital protocol;
4. Use the pre-filled syringe with lubricating gel to lubricate the catheter;
5. For filling the anchor balloon, use the pre-filled syringe with sterile water or physiological solution or water and glycerin.

**Sterilization**

Sterilized by Ethylene Oxide.

**Storage conditions**

Keep away from light and heat sources.

**Disposal after use**

Follow the current legislation for the medical waste to dispose the device after using.

**Compliance**

- CE marking according to EU Regulation MDR 2017/745, Notified Body 0373
- Non-woven fabric in compliance with the EN13795 standards
- Class IIb
- Product packaging and shelf life validated according to EN 11607 1-2 standards

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