
# evercare®

## inLine

### MONOSET, 1-WAY

<b>de</b>	Monoset, Einweg	<b>lt</b>	Mono rinkinys, 1 kanalas
<b>fr</b>	Monoset, 1 Voie	<b>pl</b>	Monoset, jednodrożny
<b>sv</b>	Blandningsaggregat, 1-vägs	<b>es</b>	Monoset de 1 vía
<b>fi</b>	Infuusiolaite monoset, 1-tie	<b>pt</b>	Monoset de 1 via
<b>da</b>	Mono sæt, 1-vejs	<b>nl</b>	Monoset, 1-weg
<b>no</b>	Monosett 1-vejs	<b>ru</b>	Моносистема
<b>et</b>	Monoset, 1-teeline	<b>ru</b>	односторонняя
<b>lv</b>	Monoset, 1-celja		

REF	Tube length (cm)
<b>10515-01</b>	180
<b>10516-02</b>	180
<b>10517-02</b>	195
<b>10518-02</b>	195

INFUSION SET ISO 8536-4 I S G
IV INFUSION SET, NEUTRAL CODE NO: OM20000

**OneMed Group Oy**  
**Metsäintie 20, FL00320 Helsinki, Finland**

### en

#### Materials Used:

Acrylonitril Butadiene Styrene (ABS), Polyolefin alloy (PVC Free), Polypropylene (PP), Nylon, Polycarbonates (PC), High-density polyethylene (HDPE), Adhesive -Acrylate.

#### Intended Purpose:

Monosets are single-use devices intended for gravity infusion of fluids from a fluid container into the vascular system of

a patient through a vascular access. Monosets can also be used for preparing and mixing drugs.

Intended Users are Healthcare Professionals. There is no restriction to a specific patient population.

Monosets are sterile and intended for single use.

Monosets principle of operation is to, by the means of gravity, conduct intravenous fluids and prepare and mix drugs, by enabling connection to a fluid container and to an infusion device with vascular access.

#### Instructions for use

1. Ensure aseptic technique is followed at all times.

2. Visually inspect the product for any imperfections or surface deterioration before use.

3. For safety reasons, tighten the cap and all disconnectable fittings before use.

4. Move the roller clamp and C-clamp from its position to inspect the tubing for any abnormalities and verify that the air-vent is closed on the side tubing (c).

5. Close the roller clamp and C-clamp and remove the spike cover.

6. Place the fluid container on the tabletop and insert the entire length of the spike (a) into it.

7. Invert the container and hang it on the IV pole or similar.

8. Attached the tube to the tube holder on the roller clamp.

9. Squeeze the drip chamber gentle 3-5 times and allow the solution to fill half of the dripchamber.

10. Open the roller clamp and fill the tubing with the solution until all air is gone from the tubing..

11. If set with integrated stopcocks, While priming stopcock, turn the tap using the directional arrows to desired port position and flush the fluid through the entire stopcock, repeat for all female ports and continue flushing until all air is completely removed.

12. Close roller clamp and remove priming/ cap and connect the male connector to the vascular access device.

13. Insert the entire length of the spike (b) into a glass vial.

14. Follow step 14.1 as below, if the product with needle free connector choose between step 14.1 or 14.2
(14.1) Open the air vent and then C-clamp on the side tubing (c).

(14.2) If the product has needle free connector attached on the side tubing, one can choose to use syringe to add diluent.

15. Allow the fluid to flow from the fluid container to glass vial until half filled.

16. Close the air vent and Shake Gentlyly until the substance is dissolved.

17. Invert the glass vial and rise it above the fluid container level.

18. Open the air vent to allow fluid from glass vial to the fluid container.

19. Close the C-clamp and air vent and leave the empty glass vial on the side tubing.

20. Set the drop rate with the roller clamp as per facility protocol.

21. Inspect and secure the connections continuously throughout the procedure.

22. Recommended to check the drop rate again after about 15-30 minutes and continues every 4-8 hour or according to local protocol.

If the set with integrated stopcock it can be used for intermittent medication or parallel infusion set.

#### Warning:

- The use of this product is restricted to a qualified doctor

or nurse.

- Read the instructions before use. The device should be used in accordance with the instructions.

- Do not use for transfusion of blood or blood components.

- The device is guaranteed sterile and pyrogen-free if the individual package has not been opened or damaged. Do not use if package is damaged or if protective caps are loose or missing.

- Do not over-tighten a luer beyond the friction fit as this may damage both the luer and another connecting device.

- Do not use any instruments for clamping the tubing or luer tightening.

- Re-use may lead to infection or other illness/injury.

#### Indication:

Used for venous Infusion of I.V. fluids, parenteral nutrition, and administration of other drugs.

Maintaining hydration and/or correct dehydration in patients who are unable to take sufficient volume of oral fluids.

#### Contraindication:

Transfusion of blood or blood components.

Not to be used in patients with known hypersensitivity to any of the materials used.

#### Disposal:

According to local regulations.

#### Notice:

Serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

These instructions for use were drawn up in English and have been translated into other languages. In the event of inconsistencies of interpretations, the English text shall be determinative.

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