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## Instructions for the mylife™ AutoLance™

EN

The lancing device is reusable for a single user. It can automatically inject the lancet into and retract it from the fingertip to obtain one or two drops of capillary blood samples for glucose monitoring or other tests.

**Loading and releasing mechanism**  
Simply press the cap on the desired penetration site.

**Depth-adjustable cap**  
Open the cover to replace the lancet.



**Penetration depth window**  
Choose the penetration depth by turning the cap.

**Mechanics holder**  
To insert the disposable lancet.

**Safety button**  
To prevent from accidental releasing.

### Handling of the mylife™ AutoLance™


- Pull off the depth-adjustable cap.
- Insert a new disposable lancet firmly into lancet carrier.
- Twist off and set aside the protective cover of the disposable lancet.
- Replace the depth-adjustable cap.
- Select the penetration depth by turning the cap. 1-3 for soft or thin skin; 4-5 for average skin; 6-7 for thick or calloused skin.
- Place the finger against the lancing device and press the safety button.
- Press your finger smoothly on the cap. The lancing device automatically loads and releases. No priming necessary.

### Alternative site testing (AST)

The AST function can enable you to test the blood glucose on other sites than the fingertips and can benefit to avoid repeated pricking in the fingertips and reduces pain.

- Select the clear AST cap to perform a measurement on an alternative site testing.
- Massage the puncture area of the palm or forearm for a few seconds and immediately after massaging the puncture area, press and release the lancing device with the clear cap against it.
- Continue holding the lancing device against the puncture site, and gradually increase pressure for a few seconds until the blood sample size is sufficient.

- Avoid being pricked by the lancet after removing the protective cover.
- For safety and to prevent cross-contamination, do not share the lancing device with others, and discard used lancets in a puncture-proof or biohazard container.
- Do not disassemble the lancing device other than for handling according to the instructions for use (see figures 1 and 2).
- The lancing device contains small parts which could result in choking if ingested.
- Storage environment: -10 to +60°C.

### Cleaning and disinfecting (Frequency: at least once a week)

Use a clean damp cloth to wipe the outside of the device. To prevent rust of internal spring, do not submerge/wash the device in water/dishwasher. Use a wipe with Isopropanol as the active ingredient to clean any dirt, blood and other body fluids. Take another disinfecting wipe to clean the device thoroughly. Allow the surface to remain wet for 2 minutes.

### Customer service

Guaranteed Times of Operation: 6000 times (5 times / day for around 3 years). If you have any questions/problems with mylife™ products, please contact local Ypsomed customer service. In case of serious incidents caused by the mylife™ product, please call emergency service for help. Please feel free to report to us and the local competent authority.

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**Bionime Corporation**
Manufactured for Ypsomed AG