REF:

13068, 13490, 13491, 13492, 13493, 13494

Spike-G Snap Cap

- de** Spike-G Schnappkappe
- fr** Spike-G Clapet de fermeture
- sv** Spike-G Snäpplock
- fi** Spike-G Napsautettava korkki
- da** Spike-G Snaplåget
- no** Spike-G Klikkhetten
- et** Spike-G Klõps kork
- lv** Spike-G Aizdares vāciņu
- lt** Spike-G Užspaudžiamą dangtelį
- pl** Spike-G Zatraskowym zamknięciem
- nl** Spike-G Snap dop

Spike-B Snap Cap & Fluid Filter

- de** Spike-B Schnappkappe & Flüssigkeitsfilter
- fr** Spike-B Clapet de fermeture & Filtre à liquide
- sv** Spike-B Snäpplock & Vätskefilter
- fi** Spike-B Napsautettava korkki & Nestesuodatin
- da** Spike-B Snaplåget & Væskefilter
- no** Spike-B Klikkhetten & Væskefilter
- et** Spike-B Klõps kork & Vedelikufilter
- lv** Spike-B Aizdares vāciņu & Šķidruma filtra
- lt** Spike-B Užspaudžiamą dangtelį & Škysčių filtrą
- pl** Spike-B Zatraskowym zamknięciem & Filtr cieczy
- nl** Spike-B Snap dop & Vloeistoffilter

Spike-R Chemo, Snap Cap & Fluid Filter

- de** Spike-R Chemo, Schnappkappe & Flüssigkeitsfilter
- fr** Spike-R Chemo, Clapet de fermeture & Filtre à liquide
- sv** Spike-R Chemo, Snäpplock & Vätskefilter
- fi** Spike-R Chemo, Napsautettava korkki & Nestesuodatin
- da** Spike-R Chemo, Snaplåget & Væskefilter
- no** Spike-R Chemo, Klikkhetten & Væskefilter
- et** Spike-R Chemo, Klõps kork & Vedelikufilter
- lv** Spike-R Chemo, Aizdares vāciņu & Šķidruma filtra
- lt** Spike-R Chemo, Užspaudžiamą dangtelį & Škysčių filtrą
- pl** Spike-R Chemo, Zatraskowym zamknięciem & Filtr cieczy
- nl** Spike-R Chemo, Snap dop & Vloeistoffilter

Product	Color	Venting filter 0.1 µm	II. Fluid Filter 5 µm
Spike-G Snap Cap	Green	X	
Spike-B Snap Cap	Blue	X	X
Spike-R Snap Cap	Red	X	X
Spike-G Needle free Connector	Green	X	
Spike-B Needle free Connector	Blue	X	X
Spike-R Chemo, Needle free Connector	Red	X	X

Spike-G Needle free Connector

- de** Spike-G Nadelfreier Konnektor
- fr** Spike-G Connecteur sans aiguille
- sv** Spike-G Nålfri koppling
- fi** Spike-G Neulaton liitin
- da** Spike-G Nålfri konektor
- no** Spike-G Nålfri kobling
- et** Spike-G Nõelavaba konnektor
- lv** Spike-G Savienojums bez adatas
- lt** Spike-G Jungtis be adatos
- pl** Spike-G Złacz bez igły
- nl** Spike-G Naaldloze connector

Spike-B Needle free Connector & Fluid Filter

- de** Spike-B Nadelfreier Konnektor & Flüssigkeitsfilter
- fr** Spike-B Connecteur sans aiguille & Filtre à liquide
- sv** Spike-B Nålfri koppling & Vätskefilter
- fi** Spike-B Neulaton liitin & Nestesuodatin
- da** Spike-B Nålfri konektor & Væskefilter
- no** Spike-B Nålfri kobling & Væskefilter
- et** Spike-B Nõelavaba konnektor & Vedelikufilter
- lv** Spike-B Savienojums bez adatas & Šķidruma filtra
- lt** Spike-B Jungtis be adatos & Škysčių filtrą
- pl** Spike-B Złacz bez igły & Filtr cieczy
- nl** Spike-B Naaldloze connector & Vloeistoffilter

Spike-R Chemo, Needle free Connector & Fluid Filter

- de** Spike-R Chemo, Nadelfreier Konnektor & Flüssigkeitsfilter
- fr** Spike-R Chemo, Connecteur sans aiguille & Filtre à liquide
- sv** Spike-R Chemo, Nålfri koppling & Vätskefilter
- fi** Spike-R Chemo, Neulaton liitin & Nestesuodatin
- da** Spike-R Chemo, Nålfri konektor & Væskefilter
- no** Spike-R Chemo, Nålefri kobling & Væskefilter
- et** Spike-R Chemo, Nõelavaba konnektor & Vedelikufilter
- lv** Spike-R Chemo, Savienojums bez adatas & Šķidruma filtra
- lt** Spike-R Chemo, Jungtis be adatos & Škysčių filtrą
- pl** Spike-R Chemo, Złacz bez igły & Filtr cieczy
- nl** Spike-R Chemo, Naaldloze connector & Vloeistoffilter

POLYMED


PICTURE PARTS
POLY MEDICURE LTD.
 Plot No. 34, Sector 68, IMT, Faridabad-121004 (Haryana), INDIA
 Customer Care No. : +91-129-4287053, M.L. : MFG/MD/2018/000032
 E-mail: customercare@polymedcure.com
 Web: www.polymedcure.com

CE/RECP Obelis S.a

Boulevard Général Wahis 53,
 1030 Brussels, Belgium
 Tel: +32(0)73025954
 E-mail: mail@obelis.net

OneMed Sverige AB
 Tagenevägen 29
 425 37 Hisings Kärra, Sweden

EN
Materials Used:

Acrylonitrile Butadiene Styrene (ABS), Polypropylene (PP), PTFE Membrane, PET (For Fluid Filter), Copolymer ETT (For Chemo Spike), Colorant, Poly Carbonate & Silicon.

Indication:

Needle free access for single vial use and multiple dose withdrawal, maximize the extractable fluid volume from the vial.

Contraindication:

Not to be used in patients with known hypersensitivity to any of the materials used.

Intended Purpose:

The product is used to steady repeated withdrawals and injection. The product is designed to be used to transfer medication from a vial to a syringe and from syringe to vial.

Venting filter 0.1 µm

Protects the user from contamination by aerosols and the medication against contamination by bacteria.

Suited for use of vial with toxic substances.

Fluid Filter 5 µm

Protects against contamination by particles.

Suited for use of vial with powder form medications.

Needle free Connector

If the vial is upside down the valve closes and is drip proof.

Instructions for use

Remove the cap from the vial and swab the with preferred solution and allow it to dry.

Remove the protective cover from the spike.

Hold the body of the Kraken® Spike and insert into the center of the vial/container ensuring the vial/container is on a flat, solid surface. (See figure no. 1)

Ensure the spike is entering vertically and not tilted or at an angle.

Using a sterile alcohol wipe, swab the surface of the valve, let it dry. (For swabbable valve only).

Push firmly using a twisting motion until the syringe enters the spike.

(A) Connect the syringe to the other end of spike. (See figure no. 2 - A For Kraken® Spike with snap fit cap)

(B) Connect the syringe to the other end of spike. (See figure no. 2 - B For Kraken® Spike with Needle free Connector)

Invert the vial and withdraw the desired amount of solution. (See figure no. 3)

Disconnect the syringe from the vial. (See figure no. 4)
 Important: Solution into the container from a filled syringe must only be injected to the vial in upright position.

Warning:

Do not inject fluid while the vial is inverted.

Do not insert the spike at an angle or in a tilted fashion.

Do not hold the vent while inserting the spike into the vial.

Do not hold the spike when the vial is inverted without syringe.

To prevent valve damage do not use needle or blunt cannula to access the Needle free Connector.

The product should be used only by qualified healthcare professionals.

Do not re-sterilize. Discard the set after single vial use.

Reuse and cleaning of product may alter their structural and mechanical properties. It may lead to infection or other illness/injury.

Poly Medicure Limited will not be responsible for any direct incidental or consequential damages resulting from use of product.

Disposal:

According to local regulations.

Notice:

Serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

These instructions for use were drawn up in English and have been translated into other languages. In the event of inconsistencies of interpretations, the English text shall be determinative.

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