

Mölnlycke<sup>®</sup>  
Exufiber<sup>®</sup> Ag+







STERILE EO

CE 2797

Manufacturer

 Mölnlycke Health Care AB  
Gamlestadvägen 3C, Box 13080  
SE-402 52 Göteborg, Sweden

[www.molnlycke.com](http://www.molnlycke.com)

Single use  
Usage unique

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## Exufiber® Ag+

Gelling fibre dressing with silver

### Product description

Exufiber Ag+ is a sterile, soft nonwoven dressing made of highly absorbent and gel-forming polyvinyl alcohol (PVA) fibres coated on both sides with silver sulphate. Exufiber Ag+ has high absorption and retention capacity of wound exudate and maintains a moist wound environment. The structural integrity of Exufiber Ag+ supports removal in one piece.

In contact with liquid, Exufiber Ag+ releases silver ions and inactivates a variety of wound related pathogens such as bacteria, fungi and moulds, as has been shown *in vitro*. By reducing the number of microorganisms, Exufiber Ag+ can prevent and reduce the formation of bacterial biofilm, as shown *in vitro* and *in vivo*. Exufiber Ag+ can be used as an effective barrier against microbial penetration of the dressing.

Product material content: Polyvinyl alcohol (PVA) fibres, hydroxypropyl cellulose and silver sulphate equal to 0.2 mg silver per cm<sup>2</sup>.

### Indications for use

Exufiber Ag+ is intended to be used in the following medium to high exuding wounds by absorbing and retaining exudate, inactivating a variety of wound related pathogens.

- Venous leg ulcers
- Diabetic foot ulcers
- Pressure ulcers

Indicated wear time: up to seven [7] days.

Intended patient population: Exufiber Ag+ can be used in patients with venous leg ulcers, diabetic foot ulcers and pressure ulcers showing signs of clinical infection.

### Contraindication(s)

For external use only.

### Precaution(s)

- Clinicians/ Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
- Do not use Exufiber Ag+ on patients and/or user with a known hypersensitivity to the ingoing materials/ components of the product.
- Frequent or prolonged use of this product may result in permanent discoloration of the skin.
- Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.
- Sterile. Do not use if sterile barrier is damaged or opened prior to use. Do not re-sterilise.

### Instructions for use

Exufiber Ag+ should be used by or under the supervision of a qualified health care professional.

#### Shallow wounds

##### Application:

Cleanse the wound with saline solution or water according to clinical practice.

1. Dry the surrounding skin thoroughly.
2. Select appropriate dressing size. For best clinical outcomes, ensure that the product covers the wound bed and overlaps the wound edges. The dressing will swell, causing contraction of the dressing edges, as it absorbs wound fluid and starts gelling.
3. Apply the dry Exufiber Ag+ dressing to the wound.
4. Fixate with an appropriate secondary dressing. The choice depends on the exudate level.
5. Compression therapy may be used in conjunction with Exufiber Ag+.

#### Deep wounds and cavities

##### Application:

Cleanse the wound with saline solution or water according to clinical practice.

1. Dry the surrounding skin thoroughly.
2. Loosely pack ribbon or pad into the wound to allow room for swelling of the dressing. When using the ribbon dressing, cut to appropriate length leaving a small overhang of 2- 3cm outside the wound for easy retrieval.
3. Cover with an appropriate secondary dressing. The choice depends on the exudate level.
4. Fixate (or cover) with an adhesive dressing or a bandage as appropriate.

#### Dressing change and removal

Both the wound and the dressing should be inspected regularly. Exufiber Ag+ can be left in place for up to 7 days, depending on wound condition or as indicated by clinical practice.

The secondary dressing should be changed based on exudate level or before becoming fully saturated.

1. Remove Exufiber Ag+ by gently cleansing/flushing with saline solution or water according to clinical practice. Any non-gelled material will moisten in contact with the saline.
2. If Exufiber Ag+ dries out and is difficult to remove, it should be moistened with saline solution or water according to local policies and allowed to soak until it lifts easily. It may take several minutes for Exufiber Ag+ to transform into a gel. Remove the dressing by gently cleansing/flushing.

The interaction of Exufiber Ag+ together with cleansing agents, oxidising agents such as hypochlorite solutions/hydrogen peroxide and topical treatments has not been demonstrated.

Prior to commencing radiation therapy, remove Exufiber Ag+ if product is present in the treatment area. A new dressing can be applied following treatment.

Avoid contact with electrodes or conductive gels during electronic measurements, e.g. electrocardiograms (ECG) and electroencephalograms (EEG).

Exufiber Ag+ is intended for short-term use. When silver-containing dressings are used on large area wounds or for prolonged periods of time (i.e. greater than 28 days), raised serum levels of silver may result which may lead to deposition of silver in tissues or organs. The use of silver-containing dressings needs to be reviewed in light of their effectiveness at this point. For longer use, a clinical assessment by a physician is recommended. In case of signs of clinical infection, consult a health care professional. Exufiber Ag+ does not replace the need for systemic or other adequate infection treatment.

#### Special storage conditions and handling conditions

Exufiber Ag+ should be stored in dry conditions below 25°C/77°F and protected from direct sunlight. The PVA fibres used in the product may change colour when it is exposed to light, air and/or heat. The colour change has no influence on product properties when used before expiry date. Disposal should be handled according to local environmental procedures.

#### Other information

For EU-customers, further information on the safety and clinical performance of Exufiber Ag+ please refer to the 'Summary of Safety and Clinical Performance for Exufiber Ag+'. You can access the document on [molnlycke.com](http://molnlycke.com). In the future the document will also be published in the European Database on Medical Devices, EUDAMED, <https://ec.europa.eu/tools/eudamed>.

If any serious incident has occurred in relation to the use of Exufiber Ag+, it should be reported to Mölnlycke Health Care and to your local competent authority.

Exufiber is a registered trademark of Mölnlycke Health Care AB.

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