smiths medical
portex®

BLUselect® Tracheostomy Inner Cannula

en Instructions for use

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sv Bruksanvisning

da Brugsanvisning

no Bruksanvisning

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Rx **CE**
ONLY **2797**

10007583-002 2019-10

BLUselect® Tracheostomy Inner Cannula INSTRUCTIONS FOR USE

These instructions apply to the use of the following Portex® BLUselect® Inner Cannula products:

Replacement Inner Cannula Code	Description	Quantity
101/850/(060-100)	Replacement Smooth Inner Cannula	2
101/850/(060CZ-100CZ)	Replacement Smooth Inner Cannula	2
101/856/(060-100)	Replacement Smooth Inner Cannula	20
101/858/(060-100)	Replacement Smooth Inner Cannula	50
101/851/(060-100)	Replacement Fenestrated Smooth Inner Cannula	2
101/851/(060CZ-100CZ)	Replacement Fenestrated Smooth Inner Cannula	20

NOTE: Product codes identified in this IFU may not be available for sale in all geographies.

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Precautions, before using this product. Failure to properly follow warnings, precautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

IMPORTANT: The list below indicates the part number of Portex® Tracheostomy Inner Cannula that aligns with the size of BLUselect® Tracheostomy Tube to be used.

NOTE 1: Use only BLUselect® Inner Cannula(s) with BLUselect® Tracheostomy Tubes.

NOTE 2: Use of an Inner Cannula reduces the airway diameter by: 1.0 mm in 6, 7, 7.5, 8, 8.5, 9 mm ID tracheostomy tubes; and 1.5 mm in 10 mm ID tracheostomy tubes.

Inner Cannula and/or Trach Tube Code Number Size Identification	Tracheostomy Tube I.D.	Inner Cannula I.D.	Jackson Size (approx.)	Tracheostomy Tube Length 'C'	Color Indication of Trach Tube
101/XXX/060	6.0 mm	5.0 mm	5	64.5 mm	Pink
101/XXX/070	7.0 mm	6.0 mm	-	70.0 mm	Green
101/XXX/075	7.5 mm	6.5 mm	7	73.0 mm	Brown
101/XXX/080	8.0 mm	7.0 mm	8	75.5 mm	Grey
101/XXX/085	8.5 mm	7.5 mm	-	78.0 mm	Purple
101/XXX/090	9.0 mm	8.0 mm	-	81.0 mm	Blue
101/XXX/100	10.0 mm	8.5 mm	10	87.5 mm	Yellow

The primary user of this instructional material is the caregiver in a clinical environment such as the adult or critical care units of a hospital, non-critical care units of a hospital, long term care facilities, and home care.

DESCRIPTION

BLUselect® Tracheostomy Inner Cannula features as follows:

- Intended to assist in the provision of a patent airway in the tracheostomized patient.
- Available in fenestrated (colored red) and non-fenestrated (clear) options.
- Fenestrations aid in ventilator weaning.
- Replacement inner cannula is designed and provided as STERILE, Single Patient Use.
- The BLUselect® Tracheostomy Inner Cannula **must not be re-sterilized** by the end user.
- The manufacturer recommends that the tracheostomy inner cannula usage not exceed twenty-nine (29) days.

Use with the following Smiths Medical Portex® BLUselect® Tracheostomy Tubes part numbers:

101/800/060-100	101/802/060-100	101/810/060-100
101/811/060-100	101/812/060-100	101/813/060-100
101/816/070-090	101/817/070-090	101/818/070-090
101/860/060-100	101/870/060-100	101/875/060-100

MATERIALS

This device is not made with plasticizer Diethylhexylphthalate (DEHP). This device is not made with natural rubber latex. Inner Cannulae are manufactured from Polyethylene.

INTENDED USE

The BLUselect® Inner Cannula is intended to be used with the Smiths Medical Portex® BLUselect® Tracheostomy Tube for patients that require an artificial airway due to trauma or medical condition.

Maximum recommended period of use is 29 days.

INDICATIONS FOR USE

The BLUselect® Inner Cannula is intended to be used with the Smiths Medical Portex® BLUselect® Tracheostomy Tube indicated for airway maintenance of tracheostomy patients.

CONTRAINDICATIONS

There are no known contraindications associated with this product.

WARNINGS

- Do not use a Fenestrated BLUselect® Tracheostomy Tube during manual or mechanical ventilation of a patient unless a Smiths Medical Portex® non-fenestrated inner cannula is in place to prevent gas leakage through the fenestrations in the tracheostomy tube.
- Verify the tracheostomy tube position by bronchoscope visualization and/or chest x-ray, lung auscultation, oximetry and capnography to ensure correct placement.
- Only use the BLUselect® Tracheostomy Inner Cannula with the Smiths Medical BLUselect® Tracheostomy tube or Smiths Medical BLUselect® Suctional® Tracheostomy tube. Product codes for compatible tracheostomy tubes are listed in this IFU.
- Avoid airway blockage by routinely checking, cleaning or replacing the inner cannula at regular intervals.
- Avoid contact with electrosurgery electrodes or laser surgery beams because PVC will produce toxic fumes in air or ignite in an enriched oxygen environment (e.g. anesthesia).
- Do not lubricate the inner cannula as this may occlude the airway patency.
- Do not lubricate the inner cannula as it may cause the inner cannula being retained within the tracheostomy tube.
- Do not fit an inner cannula when using a tracheostomy tube in conjunction with a right-angled connector that incorporates a co-axial fresh gas tube which protrudes into the connector lumen (e.g. Norman elbow). The protruding fresh gas tube could occlude the inner cannula, restricting expiratory flow of gas, resulting in barotrauma/hypoxia.
- Ensure optimal position and patency of the fenestration in the BLUselect® Fenestrated Tracheostomy Tube immediately after insertion, and frequently thereafter, to determine that the correct position of the fenestra has been obtained.
- Do not use a decannulation cap with non-Fenestrated tracheostomy tubes because it may be difficult to remove in a critical or emergency situation, potentially resulting in obstruction of the patient's airway.
- The decannulation cap is indicated for use only with fenestrated tracheostomy tubes.
- Remove the decannulation cap immediately if the patient experiences difficulty in breathing.
- When using a Decannulation Cap in conjunction with a cuffed fenestrated tracheostomy tube the following must be true:
 - The user must ensure that the cuff is completely deflated.
 - Ensure a fenestrated inner cannula is used OR the inner cannula has been removed.
- Failure to do so may lead to partial or complete airway obstruction.
- Do not reinsert kinked or damaged inner cannula into the Tracheostomy tube.
- Do not use a speaking valve (e.g. Smiths Medical Orator™) with non-fenestrated tracheostomy tubes or non-fenestrated inner cannula because it will not allow vocalization and may potentially result in obstruction of the patient's airway if used in coordination with an inflated cuff. A speaking valve is intended for use only with fenestrated tracheostomy tubes and fenestrated inner cannula.
- Patients and caregivers in the home care environment should be carefully instructed by a home health care provider in the proper use and handling of the inner cannula.

PRECAUTIONS

- Adequately humidify air delivered to patients to minimize mucus encrustation of the tracheostomy tube and/or inner cannula lumen.
- Where an inner cannula is used, ensure that it is both the correct diameter and length for the tracheostomy tube in use. Size markings are provided on the tube flange, pack labels and additional patient labels to help identify the correct inner cannula to use.
- Use of a non-fenestrated inner cannula (clear) with a fenestrated tracheostomy tube may result in air leakage between the two components. Adjustment of ventilator settings may be needed to compensate.
- The inner cannula should be routinely checked, cleaned, or replaced at regular intervals to avoid blockage or reduction in lumen of the airway.

- Do not use the cleaning brush to clean the tracheostomy tube.
- Do not use any abrasive cleaning implements to clean the inner cannula. Use only the cleaning brush provided.
- Do not use solutions other than those indicated in the cleaning instructions to clean any part of the tracheostomy tube.
- Do not try to forcibly remove the inner cannula from the tracheostomy tube, both the inner cannula and tracheostomy tube should be removed together and replaced with a new tracheostomy tube and inner cannula.
- Ensure the inner cannula does not become kinked or damaged during cleaning.
- If tracheostomy tubes are used outside the hospital, the patient must be instructed by a healthcare professional in the safe use and handling of the product. Contact Smiths Medical Customer Service for the Guidelines for the Care of Patients with Tracheostomy.

ADVERSE EVENTS

Reported adverse reactions associated with tracheostomy tubes are many and diverse. Consult standard textbooks and medical literature for specific adverse reaction information.

STERILIZATION

The BLUselect[®] Tracheostomy Inner Cannula device is Ethylene Oxide (EO) sterilized and provided to the end user in a STERILE package. The device remains STERILE as long as the package integrity has not been compromised and/or the "use by date" not exceeded. The "use by date" and integrity of the outer tray/lid should be verified prior to use; if the "use by date" has expired or packaging is compromised, the implant must not be used.

HOME CARE USE

Follow the applicable steps of the following inner cannula instructions for use as directed by your healthcare provider for your home care use needs.

WARNING: Patients and caregivers in the home care environment should be carefully instructed by a home healthcare provider in the proper use and handling of the tracheostomy tube and accessory products.

INSTRUCTIONS FOR USE

INSERT NEW REPLACEMENT INNER CANNULA

- Ensure personal protection is used while performing the replacement inner cannula procedure.
- Inspect the replacement inner cannula for kinks or damage before use.
- Following insertion of the BLUselect[®] Tracheostomy Tube and removal of the obturator in accordance with currently accepted medical techniques.
- Insert the inner cannula into the tracheostomy tube. Be sure that the inner cannula is securely inserted. A slight snap (click) should be felt indicating proper positioning.

NOTE:

- In packs where two inner cannulae are provided, one inner cannula can continue to be used while the other inner cannula is cleaned.
- Check the inner cannula prior to insertion/re-insertion and discard if kinked or damaged. Replacement inner cannulae are available (Ensure the correct size inner cannula is used).

REPLACING THE INNER CANNULA

1. When the inner cannula needs to be replaced and/or removed for cleaning disconnect the breathing circuit from the tracheostomy tube connector, if attached.
2. Remove the inner cannula by means of the ring pull on the inner cannula.
3. Select or verify the correct size of inner cannula to be used.
 - When an inner cannula is used, ensure the correct matching diameter and length for the tracheostomy tube in use.
 - Tracheostomy tube size markings and related color indication are located on the tube flange, pack labels and additional patient labels.
 - Inner Cannula size markings are located on the inner cannula and pack labels.
4. Replace soiled inner cannula with clean (or new) inner cannula if required.
5. Be sure that the inner cannula is securely inserted. A slight snap (click) should be felt indicating proper positioning.
6. Re-connect the tracheostomy tube to breathing circuit (if used) and check for security.
7. Clean the soiled inner cannula by means of the recommended cleaning instructions.
8. Check removed inner cannula prior to insertion/re-insertion and discard if kinked or damaged. Replacement inner cannula are available.

CLEANING INSTRUCTIONS

Daily cleaning of the inner cannula is recommended or whenever the inner cannula becomes contaminated by mucus or secretions. This period of time will vary according to individual patient needs.

1. Always remove the inner cannula from the Tracheostomy tube for cleaning.
2. Prepare a mild detergent solution with sterile water following the detergent manufacturer's recommendation.
3. Soak the inner cannula in mild detergent solution for 15 minutes.
4. Agitate the inner cannula in the solution to loosen any secretions. Use the cleaning brush provided to help remove any dried or sticky secretions on the inside and outside surfaces of the inner cannula.
5. After soaking and cleaning, rinse the inner cannula and cleaning brush with fresh sterile saline solution for 30 seconds.
6. Allow to dry naturally.
7. Once the inner cannula is clean and dry, it should be stored clean and dry and free from particulate matter.
8. If secretions remain after cleaning, dispose of the soiled inner cannula and insert a new or spare inner cannula.

NOTE: Refer to the BLUselect[®] Tracheostomy Tube for cleaning instructions.

STORAGE INSTRUCTIONS

The STERILE BLUselect[®] Tracheostomy Inner Cannulae should be stored in their original packaging in a dry environment, away from direct sunlight. **Unit packs must remain in the shelf carton until use to prevent damage to the inner cannula.** Proper handling, storage and utilization are important for the use of the device. All device components should be thoroughly inspected for damage prior to use. If damage is noted that may compromise the proper function of the device do not use and contact customer service or your Smiths Medical representative for a replacement.

DISPOSAL OF DEVICE

The end user is responsible for safe disposal of the tracheostomy device and that the disposal is in accordance to local guidelines for disposal of contaminated waste. Follow universal precautions as specified by the Center for Disease Control and Prevention, USA, or relevant advisory bodies in the country where this device is being used.

WARRANTY

No warranty applies as this is a single-use device for maximum use of 29 days.

CUSTOMER SERVICE

To obtain additional information regarding the BLUselect[®] Tracheostomy Tube, contact your local sales representative or call U.S. Customer Service at +1-800-258-5361 or International Customer Service at +44 (0) 845-850-0445.

COMMONLY USED SYMBOLS FOR MEDICAL DEVICES

NOTE: Refer to the individual package label and the back cover of this leaflet for symbols applicable to the product.

CAUTION: FEDERAL LAW (USA) RESTRICTS THESE DEVICES TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

Portex, BLUselect, Suctionaid, Orator, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners. © 2019 Smiths Medical. All rights reserved.

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en English

Caution

Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. ISO-15223-1 Symbols - Used with Medical Device Labelling. ISO 7000: 0434A

Do not reuse

Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure. ISO-15223-1 Symbols - Used with Medical Device Labelling. ISO 7000: 1051

Catalogue Number

Indicates the manufacturer's catalogue number so that the medical device can be identified. ISO-15223-1 Symbols - Used with Medical Device Labelling. ISO 7000: 2493

Batch Code

Indicates the manufacturer's batch code so that the batch or lot can be identified. ISO-15223-1 Symbols - Used with Medical Device Labelling. ISO 7000: 2492

Date of Manufacture

Indicates when the medical device was manufactured YYYY-MM-DD. ISO-15223-1 Symbols - Used with Medical Device Labelling. ISO 7000: 2497

Use by

Indicates the date after which the medical device is not to be used YYYY-MM-DD ISO-15223-1 Symbols - Used with Medical Device Labelling. ISO 7000: 2607

Manufacturer

Indicates the medical device manufacturer. ISO-15223-1 Symbols - Used with Medical Device Labelling. ISO 7000: 3082

Not made with natural rubber latex

Large X indicates natural rubber latex (NRL) is not used as a material in the medical product or product container. ISO-15223-1 Symbols - Used with Medical Device Labelling Annex B.2. ISO 7000: None

Quantity

Indicates number of pieces in package. Symbol is read and understood by ordinary individuals. ISO 7000: None

Authorized Representative in the European Community

Indicates the Authorized representative in the European Community. ISO-15223-1 Symbols - Used with Medical Device Labelling. ISO 7000: None

Sterilized using ethylene oxide

Indicates a medical device that has been sterilized using ethylene oxide. ISO-15223-1 Symbols - Used with Medical Device Labelling. ISO 7000: 2501

MR Conditional

ASTM: F 2503-05 - Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ISO-7000: None

Temperature Limitation

Indicates the temperature limits to which the medical device can be safely exposed. ISO - 15223-1 Symbols - Used with Medical Device Labelling. ISO 7000: 0632

Do not use if package is damaged

Indicates a medical device that should not be used if the package has been damaged or opened. ISO-15223-1 Symbols - Used with Medical Device Labelling. ISO 7000: 2606

Keep dry

Indicates a medical device that needs to be protected from moisture. ISO-15223-1 Symbols - Used with Medical Device Labelling. ISO 7000: 0626

Keep away from sunlight

Indicates a medical device that needs protection from light sources. ISO-15223-1 Symbols - Used with Medical Device Labelling. ISO 7000: 0624

Prescription Only

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician. 21CFR 801 Code Federal Regulations. ISO 7000: None

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EN
Do not reuse: Medical devices require specific material characteristics to perform as intended. These characteristics have been verified for single use only. Any attempt to re-process the device for subsequent use may adversely affect the integrity of the device or lead to deterioration in performance.
Contents sterile unless package has been opened or damaged.
Do not resterilize.

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Portex, BLUselect, Suctionaid, Orator, designmerket Portex og designmerket Smiths Medical er varemerker som tilhører Smiths Medical. Symboler © angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive eiere. © 2019 Smiths Medical. Med enerett.