Intended purpose

Elastic dressing used for support and fixation of another medical device on a body part or limb of a human being.

Intended users

0151 Tube	Users with arm, hand, leg, and foot injuries. Users with IV catheter inserted in the distal area of the extremities.
0161 Finger-Toe	Users with finger and toe injuries.
0170 Head	Users with head and neck injuries.
0186 Headband	Users with head, neck, and ear injuries.

Indications

0151 Tube	For support and fixation of dressing and/or IV line and tubing.
0161 Finger-Toe	
0170 Head	For support and fixation of dressing.
0186 Headband	

 Use the device according to indications to prevent injury or damage.

Contraindications

- Relative contraindications:
- Skin diseases
 - Contact dermatitis

Care instructions

Warnings and precautions

	0151 Tube	0161 Finger-Toe	0170 Head	0186 Headband
Follow the Instructions for use	X	X	X	X
Apply and remove the device gently to prevent further injury and dressing or catheter detachment	X	-	-	-
Use the device according to indications to prevent injury or damage	X	X	X	X
If any unexpected discomfort or skin irritation occurs, remove the device	X	X	X	X
The device is not intended for direct contact with injured skin or mucous membrane	X	X	X	X
For questions regarding use and correct fitting, consult a healthcare professional	X	X	X	X
By urgent need for removal of the device, open with scissors	X	X	X	X

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Performance characteristics

	0151 Tube	0161 Finger-Toe	0170 Head	0186 Headband
Support and fixation of dressing without tape	X	X	X	X
Support and fixation of IV line and tubing without tape for sizes XS, S and M	X	-	-	-
High and sustained elasticity	X	X	X	X
Breathable, comfortable, and soft material	X	X	X	X
Simple and quick application	X	X	X	X
Pocket for dressing	-	-	-	X
Ready to use pre-cut lengths	X	X	X	X
Simple access to the injection valve on IV catheters for sizes XS, S and M	X	-	-	-
Pre-cut opening for the thumb in sizes XS, S and M	X	-	-	-

Material composition

ART. NO	PRODUCT	MAIN FABRIC COMPOSITION
0151	Tube	70% PA, 21% EL, 9% PES
0161	Finger-Toe	79% PA, 21% EL
0170	Head	60% PA, 18% EL, 16% CO, 6% PES
0186	Headband	50% CO, 22% EL, 18% PA, 10% PES

Raw materials and components used in the device are tested for harmful substances according to STANDARD 100 by OEKO-TEX®, product class II, and relevant Biological Evaluation (ISO 10993). The device and packaging comply with REACH.

 The device may contain softening and antistatic agents up to 15%.

Reporting obligation

If you have any questions or complaints, please contact your supplier directly. In the event of a serious incident that could lead to a significant deterioration of health or to death, contact the manufacturer and the competent authority of the member state. Serious incidents are defined in article 2 No. 65 of Regulation (EU) 2017/745 (MDR).

Lifetime of the device in use

The lifetime of the medical device is determined by natural wear and tear, even when treated correctly and as recommended. In case of signs of deterioration affecting performance and comfort, the device should be disposed of.

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Carefix is a registered trademark of Tytex A/S © 2023.
All rights reserved Tytex A/S, Marsvej 6, DK-7430 Ikast
Tytex A/S has been certified according to ISO 9001, ISO 14001 and ISO 13485
TYTEX A/S • Marsvej 6 • DK-7430 Ikast • www.tytex.com