
PVC Resuscitator Set

for single use

Intended to be left blank.

DE	PVC Beatmungsbeutel Set Gebrauchsanweisung	Seite 4 - 5
EN	PVC Resuscitator Set Instructions for use.....	Page 6 - 7
CS	Sada resuscitátoru PVC Návod k použití.....	Strana 8 - 9
DA	Genoplivelsesapparat af PVC Bruganvisning	Side 10 - 11
ES	Equipo de reanimación de PVC Instrucciones de uso.....	Página 12 - 13
FR	Kit de réanimation en PVC Manuel d'utilisation	Page 14 - 15
IT	Kit pallone rianimatore in PVC Istruzioni d'uso.....	Pagina 16 - 17
RU	Комплект реанимационного аппарата (ПВХ) Инструкция по использованию	стр. 18 - 19
TR	PVC Resüsitatör Seti Kullanım kılavuzu	Sayfa 20 - 21
	Symbol description	Page 22 - 23

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INTENDED USE / INDICATION

Any emergency with limited respiratory function or signs of hypoxia.

This PVC resuscitator is intended to be used for the following patients:

Patient	 Body weight
Infants	< 5 kg
Children	5 - 20 kg
Adults	> 20 kg

Place of use: No restrictions

CONTRAINDICATION

In patients who do not have an empty stomach and with increased risk of aspiration, ventilation must not be employed except in emergency situations.

SAFETY INFORMATION

- Read the instructions for use carefully before use on a patient.

- Intended for use by qualified personnel.
- Prior to use, a functional check must be performed (see section "Functional check"). Damaged products must not be used.
- Do not use the product in a toxic atmosphere.
- Do not use oil, grease or any hydrocarbon-based substance on the product so as to avoid possible flammability during use.
- Do not disassemble the pressure relief valve.
- Use original VBM individual parts only. Using individual parts from other manufacturers may have a negative impact on the function of the medical device. The manufacturer accepts no liability for this.

- For best performance of this product, please use within 5 years from the date of manufacture.

- The product is intended for single use and must not be reused. Reprocessing influences the function of the product. Any reuse entails a potential risk of cross-contamination.

FUNCTIONAL CHECK

The valve functions must be checked as follows to ensure proper functionality of the resuscitator:

A reservoir bag is needed to complete the test procedures described below:

1.0 Intake valve

- a) ▶ Firstly compress the resuscitator with one hand and then close the patient valve of the resuscitator with the other hand. Release your grip on the resuscitator.

Rapid re-expansion of the resuscitator confirms efficient air intake via the intake valve.

- b) ▶ Close the patient valve of the resuscitator and try to compress the resuscitator.

If the resuscitator cannot be compressed with reasonable force, or if resuscitator compression forces the air out between your hand and the neck of the resuscitator, the intake valve efficiently prevents backward leakage of air.

2.1 Patient valve (pressure relief valve in closed position)

- a) ▶ Hold a reservoir bag over the patient valve and use your thumb to press on the reservoir bag connector.

- ▶ Ensure a tight seal between the patient valve and reservoir bag.

- ▶ Compress the resuscitator with your other hand several times. Ensure that the lip valve opens during compression.

Filling of the reservoir bag confirms that the patient valve efficiently directs air to the patient.

- b) ▶ Compress the filled reservoir bag attached to the patient valve and inspect the flutter valve.

Rising of the flutter valve confirms that air is properly released into the atmosphere, rather than returning into the resuscitator.

2.2 Patient valve (pressure relief valve in opened position)

- ▶ Close the patient valve with your thumb and compress the resuscitator several times. The resuscitator can be compressed and the visible and audible opening of the pressure relief valve confirms its proper operation.

3.0 Reservoir valve (in the intake valve)

- a) ▶ Proceed as described above in 2.1a in order to fill the reservoir bag with ambient air.

- ▶ Attach the filled reservoir bag to the intake valve and press on the reservoir bag.

The visible rise of the diaphragm when compressing the reservoir bag confirms that the reservoir bag efficiently vents excessive gas to the atmosphere.

- b) ▶ Proceed as described above in 2.1a in order to fill a reservoir bag with ambient air.

- ▶ Attach the filled reservoir bag to the intake valve.
- ▶ Perform several compression-release cycles on the resuscitator until the reservoir bag is flat and empty. Rapid re-expansion of the resuscitator confirms efficient air intake.

USE

- ▶ Optional: Connect reservoir bag and oxygen tube to the resuscitator.

- ▶ Connect the mask to the patient valve.

- ▶ Connect the oxygen tube to a regulated oxygen source.

- ▶ Adjust the oxygen flow rate. The reservoir bag expands completely during inflow and nearly collapses as the compressed resuscitator refills during exhalation.

- ▶ Prior to connecting to a patient, check the function of the resuscitator and make sure that all connections are in proper position. Observe the intake valve, reservoir bag and patient valve in all phases of ventilation; leakages must not occur.

- ▶ Press the mask tightly to the patient's face so that it forms an airtight seal. If the mask does not fit properly, the ventilation may be affected.

- ▶ Fix the mask by hand.

- ▶ Compress the resuscitator to deliver a breath. Observe the patient's chest wall rise to confirm inspiration.
- ▶ Release the resuscitator to allow patient exhalation. Observe the chest wall fall to confirm exhalation.

SHELF LIFE

The shelf life of the product is 5 years.

Use-by date: see product label.

STORAGE AND TRANSPORT CONDITIONS

Operating temperature:
-18 °C ... +50 °C (-0.4 °F ... 122 °F)
Storage temperature:
-40 °C ... +60 °C (-40 °F ... 140 °F)

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

Used or damaged products must be disposed of in accordance with the applicable national and international regulations.

TECHNICAL DATA

	Infants	Children	Adults
Inspiratory resistance	< 5 cmH ₂ O @ 50 l/min		
Expiratory resistance	< 5 cmH ₂ O @ 50 l/min		
Pressure relief valve	Gas release when pressure ≥ 40 cmH ₂ O		
Dead space	< 6 ml		
Expected delivery volume	up to 125 ml	up to 250 ml	up to 675 ml
Product dimension assembled (L x W x D)	460 x 140 x 120 mm	550 x 140 x 145 mm	630 x 140 x 150 mm
Mass	< 500 grams		

Conversion of pressure units: 1 hPa = 1.01973 cmH₂O = 0.75006 mmHg

SCOPE OF DELIVERY

REF	84-10-399	84-10-395	84-10-299	84-10-295	84-10-199	84-10-195
	Infants		Children		Adults	
Resuscitator	280 ml		450 ml		1500 ml	
Mask	# 1		# 3		# 5	
Reservoir bag	600 ml		2500 ml		2500 ml	
Oxygen tube	200 cm					
Packaging unit	1	6	1	6	1	6

MATERIAL DATA

Item	Material
Patient valve	PC (polycarbonate), silicone
Resuscitator	PVC (polyvinyl chloride)
Intake valve set	PC (polycarbonate), silicone
Reservoir bag	PVC (polyvinyl chloride)
Mask	PVC (polyvinyl chloride)
Oxygen tube	PVC (polyvinyl chloride)

Oxygen flow rate	Percent oxygen delivered		
Resuscitator patient category	Device under test Infants	Device under test Children	Device under test Adults
2 l/min	92 %	69 %	50 %
4 l/min	96 %	95 %	74 %
6 l/min	98 %	98 %	96 %
8 l/min	99 %	99 %	98 %
10 l/min	99 %	99 %	99 %
15 l/min	99 %	99 %	99 S
All with assembled oxygen reservoir bag.			

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