EN	NEUTRAL ELECTRODES
DE	NEUTRALELEKTRODEN
FR	ELECTRODES NEUTRES
IT	ELETTRODI NEUTRI
BG	НЕУТРАЛЕН ЕЛЕКТРОД
CS	NEUTRÁLNÍ ELEKTRODY
DA	NEUTRALELEKTRODER
EL	ΟΥΔΕΤΕΡΑ ΗΛΕΚΤΡΟΔΙΑ
ES	ELECTRODOS NEUTROS
ET	NEUTRAALELEKTROODID
FI	NEUTRAALIELEKTRODIT
HR	NEUTRALNE ELEKTRODE
HU	NEUTRÁLIS ELEKTRÓDÁK
LT	NEUTRALŪS ELEKTRODAI
LV	NEITRĀLIE ELEKTRODI
NL	NEUTRALE ELEKTRODEN
PL	ELEKTRODY NEUTRALNE
PT _{EU}	ELÉTRODOS NEUTROS
PT _{BRASIL}	PLACAS TERRA
RO	ELECTROZI NEUTRI
SK	NEUTRÁLNE ELEKTRÓDY
SL	NEUTRALNE ELEKTRODE
SV	NEUTRALELEKTRODER
KO	중립 전극
NO	NØYTRALE ELEKTRODER
RU	НЕЙТРАЛЬНЫЕ ЭЛЕКТРОДЫ
TH	อิเล็กโทรดที่เป็นกลาง
TR	NÖTRAL ELEKTROTLAR
ZH	中性电极
AR	الإلكترودات المحايدة

NEUTRAL ELECTRODES FOR ADULTS AND CHILDREN

Please read this document and keep it safe. Ensure that any individual using this product is familiar with and understands the information contained in this document.

INTENDED USE

Neutral electrodes are intended as accessories for electrosurgical devices. Neutral electrodes are non-sterile; they are intended only for single use on intact (uninjured) skin and lead energy from the patient to the electrosurgical device.

WARNING

Improper use of neutral electrodes can cause damage to tissue. These instructions for use serve to ensure patient safety. **NOT FOLLOWING THESE INSTRUCTIONS MAY LEAD TO BURNS, PRESSURE NECROSES, OR OTHER SKIN TRAUMA DURING USE.**

- Do not use the neutral electrode if it is damaged, modified, or expired.
- If an electrosurgical unit offers a neutral electrode contact quality monitoring system always use a split electrode. Never deactivate the acoustic alarm of this monitoring system during surgery.
- Use the lowest effective safe power setting to achieve the desired surgical effect.
- Do not operate the electrodes in the immediate presence of flammable gases (including concentrated oxygen) and solvents to avoid possible explosion or fire hazard. Avoid contact between the product and solvent-based liquids.

PRODUCT LIMITATIONS:

- Exceeding the following limitations may overload the neutral electrode with current. This may result in a patient burn despite a fully and correctly applied neutral electrode and an activated contact quality monitoring system.
- These neutral electrodes are designed for use during conventional monopolar electrosurgical procedures.
- During non-conventional procedures (e.g. high current mode) involving high electrical currents, long activation periods, or both (e.g. tissue ablation, tissue vaporization, or procedures involving the application of a conductive rinsing fluid to the surgical site), there is a risk of inflicting burns on the patient. For such procedures, consult the generator and accessory manufacturers' instructions, in particular regarding limitations of activation time. Use additional neutral electrodes when indicated.

Product limitations for neutral electrodes for children

This neutral electrode for children is intended for use only on children weighing 5–15 kg (11–33 lbs). Restrict the duration of activation to a maximum of 60 seconds within each 4 minute period. Only use this neutral electrode for children at power settings of up to 120 W.

		
✓		
5 - 15 kg 11 - 33 lbs		

Product limitations for universal neutral electrodes

These universal neutral electrodes are only intended for use on patients weighing more than 5 kg (11 lbs). Restrict the duration of activation to a maximum of 60 seconds within each 2 minute period. **During conventional procedures on patients weighing more than 15 kg [according to IEC 60601-2-2:2017] only split neutral electrodes must be used.**

		
✓	✓	
> 5 kg > 11 lbs		

Product limitations for neutral electrodes for adults

These neutral electrodes for adults are only intended for use on patients weighing more than 15 kg (33 lbs). Restrict the duration of activation to a maximum of 60 seconds within each 2 minute period. **During conventional procedures on patients weighing more than 15 kg [according to IEC 60601-2-2:2017] only split neutral electrodes must be used.**

		
✓		

COMPATIBILITY

For any compatibility questions, particularly regarding compatibility with a specific electrode contact quality monitoring system, please contact your local distributor. In any case, the split neutral electrodes are compatible with the dynamic contact quality monitoring systems REM™, NESSY®, and ARM™.

CAUTION

- When using on children with a body weight of up to 15 kg (33 lbs): Surgeons should be familiar with the effects of electrosurgery on small patients, and should first consider the use of bipolar electrosurgery, which does not require a neutral electrode!
- Do not use the neutral electrode if the packaging is damaged.
- Do not trim down the neutral electrodes to reduce their size! Do not use additional gel!
- Do not reposition the neutral electrode!
- Before using the electrosurgical device, familiarize yourself with its instructions for use, in particular with respect to power setting restrictions and the maximum uninterrupted activation period, as well as information on the use of neutral electrodes for applications with a high energy output.
- Should the coagulating or cutting effect diminish during surgery, or a higher than normal power setting be required, immediately make sure that the neutral electrode is adequately placed and in full contact with the skin. Inspect all connections of the neutral electrode, i.e. clamp, cable, connector, plus the active accessories and the generator before turning up the power setting!
- If you reposition the patient, make sure the entire surface of the neutral electrode still sticks well to the skin and check all cable connections.
- Do not reuse the neutral electrode! If re-used, the adhesive strength and the electrical properties may be insufficient, which could lead to a patient injury. There is also a risk of cross-infection from one patient to another.
- Store any unused electrodes in the original pouch. Close the pouch by folding the open end once or several times in order to keep the electrodes fresh.
- Use all electrodes within 7 days of opening the pouch.


INSTRUCTIONS FOR USE

These neutral electrodes must only be used by appropriately qualified medical personnel who have been trained in their use based on these instructions for use.

SELECTION AND PREPARATION OF THE PLACEMENT SITE

- Select a convex skin surface that is muscular or well supplied with blood and as close as possible to the surgical site, but if possible at least 15 cm away from it:
 - **on adults** preferably on the upper arm or upper thigh.
 - **on children** between 5 and 15 kg, preferably on the torso or upper thigh.
- Whenever the size of the child patient allows, use a neutral electrode for adults or a universal neutral electrode in order to reduce the likelihood of patient burns.
- Ensure that the patient's weight is not bearing upon the selected site during the operation, and that it is not subject to any other pressure (for example due to compression stockings). Make sure that the selected site is not thermally insulated or subject to additional heat from heat sources during the operation.
- Avoid very hairy skin surfaces and surfaces above metallic implants, scar or adipose tissue, bony protrusions, injection sites, tattoos, and inflamed or damaged areas. Avoid areas in which fluids could collect.
- If the patient has a pacemaker or other active implants, consult a qualified professional to determine whether high frequency surgery is suitable, and if so where to place the neutral electrode and neutral electrode cable.


- For surgeries during which the HF current may flow through parts of the body with a relatively small cross-sectional area, bipolar methods may be more suitable in order to avoid unintentional damage to tissue.
- Monitoring electrodes or other devices that may provide an alternative earthing route for the high-frequency current should be applied as far from the surgical site as possible. It is recommended to use only ECG and other monitoring systems, cables, and wiring equipped with HF current limiting devices (e.g. high frequency filters or chokes). If this is not possible, the neutral electrode must always be applied closer to the surgical site than these electrodes or devices. The use of needle monitoring electrodes is not recommended.
- Remove any hair from the chosen skin area and clean it carefully, e.g. of cosmetics. Carefully dry the area, especially if alcohol or other skin cleaning agents have been used. Avoid using flammable skin prepping agents or disinfectants, e.g. acetone degreasers. **Be aware that failure to shave the area may lead to burns.**
- Avoid skin-to-skin contact (for example between the patient's arms and body) by inserting insulation between the relevant areas, for example dry gauze.
- Remove any metallic jewelry.

NEUTRAL ELECTRODE APPLICATION

- If the electro-surgical device is equipped with a system for monitoring neutral electrode contact quality (such as REM™, NESSY®, ARM™, etc.), always use a split electrode. A contact quality monitoring system may not work in conjunction with an unsplit electrode, meaning that the loss of secure contact between the patient and the neutral electrode will not trigger an audible alarm. Test the functioning of the monitoring system by attempting to use the device without the neutral electrode connected. The device should not be activated and an alarm should sound.
- Open the pouch immediately before use, and remove the neutral electrode. Check the expiry date on the pouch. Do not use electrodes that have passed their expiry date. Store any unused electrodes in the original pouch. Close the pouch by folding over the open end once or several times in order to keep the electrodes fresh. Use all electrodes within 7 days of opening the pouch.
- If applicable, remove one of the traceability labels from the pouch and stick it in the patient's record. In the patient's record, document the position of the electrodes, the skin preparation performed, and the condition of the skin.
- Remove the electrode, as indicated, from its protective liner using the removal tab. Check the electrode and the cable and all connectors for defects, e.g. dried out or missing gel and damaged cable insulation. Do not use defective devices.
- Align the neutral electrode so that one of its long sides is closest to the surgical site. For electro-surgical generators, which are equipped with the system to monitor the neutral electrode contact quality as well as a system to monitor an even current distribution across both conductive surfaces of the neutral electrode, the split electrode must be aligned so that the split is pointing towards the surgical site, otherwise alarms may be triggered.
- Apply the neutral electrode to the prepared area of skin, starting with one corner, and applying even pressure over the entire surface without stretching the skin or the electrode. Prevent the emergence of air bubbles or skin folds beneath the electrode. Gently rub your hand over the electrode to ensure good contact of the entire adhesive surface to the skin.
- Do not wrap the electrode entirely around an extremity. The electrode must not be in contact with or overlap itself.
- Neutral electrodes without cables: Check the reusable neutral electrode cable for defects. Do not use any neutral electrode cables with dirty metallic electrode contacts or other defects such as damaged insulation. Open the clamp of the neutral electrode cable by lifting the locking lever. Fully insert the connection lug of the electrode into the clamp. Close the clamp by pushing the locking lever

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- down as far as it will go. Check that the connection lug is fully inserted into the clamp and will not come into contact with the patient's skin. The clamp must not end up underneath the patient.
- Position the neutral electrode cable so that the cable is not pulling the neutral electrode away from the patient's skin. The neutral electrode cable should also be positioned to avoid contact with the patient and other cables, and to ensure there are no loops. In particular, do not wrap the neutral electrode cable around one of the patient's extremities or other grounded objects, in order to prevent capacitive coupling and resulting burns.
- The neutral electrode cable must be pointing away from the surgical site at an angle of 90° or 180°.
- Check that the neutral electrode is securely adhered to the skin across its entire surface, and that the clamp is securely attached to the connection lug of the electrode. Make sure that the clamp is not exerting any unnecessary pressure on the patient's skin.
- Connect the neutral electrode cable with the electro-surgical device in accordance with the instructions for use for the generator.
- Before using the electro-surgical device, familiarize yourself with its instructions for use, in particular with respect to power setting restrictions and the maximum uninterrupted activation period, as well as information on the use of neutral electrodes for applications with a high energy output.
- Never deactivate the acoustic alarm of the contact quality monitoring system for neutral electrodes during surgery.
- If the patient is moved, ensure that the entire surface of the neutral electrode remains securely adhered to the skin, and check all cable connections (clamp, cable, connector).

REMOVING AND DISPOSING OF THE NEUTRAL ELECTRODE

- After use, carefully remove the neutral electrode using one hand whilst supporting the skin beneath it with the other. Use the removal tab to slowly detach the electrode (not using the neutral electrode cable). Ripping off the electrode or peeling it off quickly may damage the skin. Exercise particular caution if the skin is exceptionally sensitive, especially with children but also with elderly patients, diabetics, or patients on a prolonged course of certain medications known to cause dermatological side effects.
- To open the cable clamp, lift the locking lever.
- Dispose of the product in accordance with local regulations or the guidelines of the hospital/health care facility.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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