evercare® inLine

GRAVITY SET

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en

Materials Used:

Acrylonitril Butadiene Styrene (ABS), Polyolefin alloy (PVC Free), Polypropylene (PP), Nylon/ Polytetrafluoroethylene (PTFE), Polycarbonates (PC), High-density polyethylene (HDPE), Adhesive -Acrylate.

Intended Purpose:

Gravity sets are single-use devices intended for gravity infusion of

fluids from a fluid container into the vascular system of a patient through a vascular access.

Intended Users are Healthcare Professionals. There is no restriction to a specific patient population.

Gravity Sets are sterile and intended for single use.

Gravity Sets principle of operation is to conduct intravenous fluids through gravity method, by enabling connection to a fluid container and to an infusion device with vascular access.

Instructions for Use:

- Ensure aseptic technique is followed at all times
- Visually inspect the product for any imperfections or surface deterioration before use.
- For safety reasons, tighten the caps and all disconnectable fittings before use.
- Move the roller clamp from its position to inspect the tubing for any abnormalities and verify that the air-vent is closed in vented devices.

- Close the roller clamp and remove the spike cover.

- Place the fluid container on the tabletop and insert the entire length of the spike into it.

- Invert the container and hang it on the IV pole or similar.

- Attached the tube to the tube holder on the roller clamp.

- Squeeze the drip chamber gentle 3-5 times and allow the solution to fill half of the drip chamber.

- Open the roller clamp and fill the tubing with the solution until all air is gone from the tubing.

- Close roller clamp and remove priming cap and connect the male connector to the vascular access device.

- If venting is required open the air vent now.

- Set the drop rate with the roller clamp as per facility protocol.

- Inspect and secure the connections continuously throughout the procedure.

- Recommended to check the drip rate again after about 15-30 minutes and continue every 4-8 hour or according to local protocol.

- Device with integrated needle free connector can be used for intermittent medication.

Warning:

- The use of this product is restricted to a qualified doctor or nurse.

- Read the instructions before use. The device should be used in accordance with the instructions.

- The device is guaranteed sterile and pyrogen-free if the individual package has not been opened or damaged. Do not use if package is damaged or if protective caps are loose or missing.

- Do not over-tighten a luer beyond the friction fit as this may damage both the luer and another connecting device.

- Do not use any instruments for clamping the tubing or luer tightening.

- Re-use may lead to infection or other illness/injury.

Indication:

Used for venous Infusion of I.V. fluids, parenteral nutrition, and administration of other drugs.

Maintaining hydration and/or correct dehydration in patients who are unable to take sufficient volume of oral fluids.

Contraindication:

Not to be used in patients with known hypersensitivity to any of the materials used.

Transfusion of blood or blood components.

Disposal:

According to local regulations.

Notice:

Serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

These instructions for use were drawn up in English and have been translated into other languages. In the event of inconsistencies of interpretations, the English text shall be determinative.

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