

Tennis Elbow Gel/Air Support




Figure 1




Figure 2



Figure 3

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

Intended for single patient use.



WARNING: If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.

APPLICATION INSTRUCTIONS

1. Select gel pad or air pad and insert under the elastic band on the support.

NOTE: For cold therapy, freeze gel pad for at least one hour before application.

2. Pull strap end through the plastic link to form a loose circle (figure 1).

3. Slide the support onto the affected arm and position 2-3" below elbow.

4. Rotate support to place pad insert on muscle just below the point of pain.

5. Pull strap and secure at desired fit and level of compression (figure 2).

6. Adjust strap as needed for comfort

NOTE: Please refer to sew-in label on device for care instructions.

NOTICE: Every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort with this product, however, there is no guarantee that injury will be prevented with its use.

**NOT MADE WITH
natural rubber latex**