

Drain-Guard

Intended Use

Drain-Guard is intended to anchor the percutaneous drainage catheter.

Contraindications

Do not use on patients with known tape or adhesive allergies.

Patient population

All groups dedicated to have the percutaneous drainage catheter anchored.

Product description

The Drain-Guard is intended to anchor percutaneous drainage catheters that are protruding from the patient's skin in an angle approximately perpendicular to the patient's skin surface. The device consists of a skin-side surface with a skin friendly adhesive combined with a specially designed foam block that safely anchors catheters from 5 to 16 French. The soft foam support block can prevent kinking of the catheter. A transparent hydrocolloid cover allows continuous monitoring of the puncture site and provides an efficient protection against contamination.

Precautions

- For single use only. Do NOT reuse, reprocess or re-sterilize this device. Any reprocessing may impede the functions of the device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device results in risk of device malfunction.
- Read the manufacturer's instruction prior to use.
- Carefully inspect the packaging for any damage prior to use. Do NOT use if the packaging is damaged. Do NOT use past expiration date. DO NOT use if the device is exposed to non-sterile surface before procedure.
- Remove the Drain-Guard and replace with a new Drain-Guard fixation device if needed.
- While removing the Drain-Guard fixation device use special care not to move or accidentally pull out catheter.
- In case of an appearance of skin lesions remove the device.

Instruction for Use

In order to achieve effective adherence of the Drain-Guard, device must be applied to a clean and dry skin surface.

Application instructions are as follows:

1. Take the Drain-Guard out of the pouch. (Figure 1)
2. Open front slit for catheter insertion. (Figure 2a)
3. Place Drain-Guard so catheter is positioned almost at bottom of slit. (Figure 2b)
4. Remove two small pieces of siliconized paper (a and b) by pulling the tab and press lightly to affix adhesive to the skin. (Figure 3)
5. Remove large siliconized paper (c) by pulling the tab and press lightly on product to affix adhesive to the skin. (Figure 3)
6. Remove siliconized paper on top of Drain-Guard. (Figure 4)
7. Place catheter on the exposed adhesive area. (Figure 5)
8. To apply the top cover pull tab of siliconized paper. (Figure 6a)
9. Position top cover with the siliconized paper. (Figure 6b)

10. Press the top cover (a) to secure good adhesion to the top of the carrier and press the tongue-shaped hydrocolloid (b) around the puncture site. (Figure 7)
11. Drain-Guard application is finished. (Figure 8)

In order to remove Drain-Guard follow the steps:

1. Peel off tongue-shaped hydrocolloid around puncture site. (Figure 9a)
2. Peel off lid material along pre-cut lines. (Figure 9a)
3. Lift catheter (Figure 9b) away from adhesive top surface and cautiously remove Drain-Guard from skin surface. (Figure 9c)
4. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy. (Figure 9d)

In case of any unforeseen medical condition please contact healthcare professional immediately.

Reporting

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical ApS and the competent regulatory authority of the country in which the user and/or patient is established.

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