	WORK INSTRUCTION FM-DVSOP2002-01 Rev. 001 Eff. Date: 25-OCT-2016 CO No.: CO-10029546		Doc. Number AS10002537-002	Revision 100
	WI Type: ARTWORK SPECIFICATION		Effective Date: 19 APR 2017	Page 1 of 19
Title: INSTRUCTIONS FOR USE - BLUE LINE ULTRA® TRACHEOSTOMY INNER CANNULA, SINGLE PATIENT USE (FORM FILL SEAL - SOFT PACK ONLY) (10002537-002)				

DOCUMENT VERIFICATION

Prior to using this document, the user is responsible for verifying that the revision and effective date are current.

REVISION HISTORY

Owning Site: MASH		Additional Sites Using Document: MCZH		
Rev.	Effective Date	Change Order No.	Changes Made to Document	
100	19 APR 2017	CO-10037072	PWR: 5406 This artwork replaces CCN. 10002537-001 (LLS10002537-001 Issue 100). Updates include the following changes; <ul style="list-style-type: none"> Warning 3 updated: 'us' has been amended to 'use' in the English text only. 	

LANGUAGES

en	fr	de	it	es	pt	nl	sv	da	no
English	French	German	Italian	Spanish	Portuguese	Dutch	Swedish	Danish	Norwegian
Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
fi	el	ja	cs	sk	hu	pl	sl	et	lv
Finnish	Greek	Japanese	Czech	Slovakian	Hungarian	Polish	Slovenian	Estonian	Latvian
Y	Y	N	Y	N	Y	Y	N	Y	Y
lt	tr	bg	ro	ru	hr	zh	ko	ar	
Lithuanian	Turkish	Bulgarian	Romanian	Russian	Croatian	Chinese	Korean	Arabic	
N	Y	Y	N	N	N	N	N	N	

ARTWORK DESIGN NOTES

CCN: 10002537-002

Size: Internal pages: 210mm (W) x 148mm (H). **Finished Folded Size:** 105mm (W) x 148mm (H).

External pages: 315mm (W) x 148mm (H). **Finished Folded Size:** 105mm (W) x 148mm (H).

IFU Format: Printed Instruction Booklet, 6 pages external cover and 60 pages internal. (See folding diagram on page 2).

A) For detailed IFU specifications refer to item PRS.

B) The layouts comply with the requirements of DVSOP8005 REV. 000.


WI Type: **ARTWORK SPECIFICATION**

Effective Date: 19 APR 2017


Page 2 of 19

Title: **INSTRUCTIONS FOR USE - BLUE LINE ULTRA® TRACHEOSTOMY INNER CANNULA, SINGLE PATIENT USE (FORM FILL SEAL - SOFT PACK ONLY) (10002537-002)**

FOLDING DIAGRAM




Internal Pages

Back Cover

Front Cover

Internal Flap

FINISHED FOLDED SIZE A6

○ = Stitch locations.
Internal pages stitched twice to cover.

	WORK INSTRUCTION FM-DVSOP2002-01 Rev. 001 Eff. Date: 25-OCT-2016 CO No.: CO-10029546		Doc. Number AS10002537-002	Revision 100
	WI Type: ARTWORK SPECIFICATION		Effective Date: 19 APR 2017	Page 3 of 19
Title: INSTRUCTIONS FOR USE - BLUE LINE ULTRA® TRACHEOSTOMY INNER CANNULA, SINGLE PATIENT USE (FORM FILL SEAL - SOFT PACK ONLY) (10002537-002)				

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66

Blue Line Ultra® Tracheostomy Inner Cannula, Single Patient Use

en INSTRUCTIONS FOR USE	sv BRUKSANVISNING	pl INSTRUKCJA STOSOWANIA
fr MODE D'EMPLOI	da BRUGSANVISNING	et KASUTUSJUHISEID
de GEBRAUCHSANLEITUNG	no BRUKSANVISNING	lv LIETOŠANAS INSTRUKCIU
it ISTRUZIONI PER L'USO	fi KÄYTTÖOHJEET	tr KULLANIM TALIMATI
es INSTRUCCIONES DE USO	el ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	bg УКАЗАНИЯ ЗА УПОТРЕБА
pt INSTRUÇÕES DE UTILIZAÇÃO	cs NÁVOD K POUŽITÍ	
nl GEBRUIKSAANWIJZING	hu HASZNÁLATI UTASÍTÁS	


Manufacturer:
Smiths Medical Czech Republic a.s.
 Olomoucká 306, 753 01 Hranice
 Czech Republic
 Tel: +44 (0)1233 722100

www.smiths-medical.com

CCN: 10002537-002 03/17

Page intentionally blank

English

en

Portex® Blue Line Ultra® Tracheostomy Inner Cannula, Single Patient Use

INSTRUCTIONS FOR USE

These instructions are intended for use with the following Portex® products:
 Catalogue No: 100/850/060-100 Replacement smooth inner cannula x 2
 Internal REF: 100/850/060CZ-100CZ

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

Please note that the code number of the Portex® Single Patient Use Inner Cannula indicates the size of Portex® Tracheostomy tube with which it should be used e.g.

Inner Cannula Code No.	For use with Tracheostomy Tube	Inner Cannula I.D.	Jackson Size (approx.)	Tracheostomy Tube Length 'C'
100/850/060	6.0 mm	5.0 mm	5	64.5 mm
100/850/070	7.0 mm	6.0 mm	-	70.0 mm
100/850/075	7.5 mm	6.5 mm	7	73.0 mm
100/850/080	8.0 mm	7.0 mm	8	75.5 mm
100/850/085	8.5 mm	7.5 mm	-	78.0 mm
100/850/090	9.0 mm	8.0 mm	-	81.0 mm
100/850/100	10.0 mm	8.5 mm	10	87.5 mm

NOTE:

- Use only Blue Line Ultra® Inner Cannula with Blue Line Ultra® Tracheostomy Tubes.
- The use of an Inner Cannula reduces the airway diameter by 1.0mm (excluding size 10.0mm tracheostomy tube which reduces the airway diameter by 1.5mm).

Description

A range of sterile single patient use replacement inner cannula intended to assist in the provision of a patent airway in the tracheostomised patient, available in both fenestrated (coloured red) and unfenestrated (clear) versions for use with the following Portex® Tracheostomy tubes.

100/800/060-100	100/811/060-100	100/831/070-090	100/860/060-100
100/802/060-100	100/812/060-100	100/832/070-090	100/870/060-100
100/810/060-100	100/813/060-100	100/833/070-090	

Instructions for use

- Following insertion of the tracheostomy tube and removal of the obturator, insert the inner cannula fully until a click is felt indicating correct positioning.

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- If the inner cannula needs to be removed for cleaning disconnect from breathing circuit if attached.
- Remove the inner cannula by means of the ring pull on the inner cannula.
- Replace the soiled inner cannula with a clean inner cannula if required ensuring the inner cannula is correctly positioned by means of the click.
- Re-connect the tracheostomy tube to breathing circuit (if used) and check for security.
- Clean the soiled inner cannula by means of the recommended cleaning instructions.
- Check the inner cannula prior to insertion/re-insertion and discard if kinked or damaged. Replacement inner cannulae are available (100/850 or 100/851).

Cleaning instructions

Daily cleaning of the inner cannula is recommended or whenever the inner cannula becomes contaminated by mucous or secretions. This period of time will vary according to individual patient needs. Always remove the inner cannula from the tracheostomy tube for cleaning.

- Soak the inner cannula in sterile saline or mild detergent solution.
- Agitate the inner cannula in the solution to loosen any secretions. Using the cleaning brush provided with all Portex® arc shaped tracheostomy tubes (except 100/800 and 100/802) to help remove any dried or sticky secretions.
- After soaking and cleaning, rinse the inner cannula and cleaning brush with fresh sterile saline solution.
- Allow to dry naturally. Once clean and dry the inner cannula should be stored clean and dry and free from particulate matter.

Cleaning instructions – Tracheostomy Tube

The tracheostomy tube should only be cleaned with a sterile saline solution.

Precautions

- Unit packs must remain in the shelf carton until use to prevent damage to the inner cannula.
- The security of all breathing system connectors should be checked when established and frequently thereafter. Disconnection may be facilitated with the use of a disconnection wedge.
- Patients should be adequately humidified to minimise encrustation of the tracheostomy tube and/or inner cannulae.
- The patency of the tracheostomy tube lumen must be assured by regular suctioning and regular cleaning of the inner cannula. Check routinely and replace as required to maintain a patent airway. Maximum recommended period of use 30 days.
- Care should be taken to ensure that the inner cannulae do not become kinked or damaged during cleaning, and that no kinked or damaged inner cannulae are re-inserted into the tracheostomy tube.
- The position and patency of the fenestration in the Portex® Fenestrated Tracheostomy tube and inner cannulae should be checked immediately after insertion to determine optimal position of fenestra and frequently thereafter to ensure that tissue granulation has not caused an obstruction.
- If Portex® Tracheostomy tubes and inner cannulae are used outside the hospital, the patient must be instructed by a healthcare professional in the safe use and handling of the product.
- Do not use solutions other than those indicated in the cleaning instructions to clean any part of the tracheostomy tube.
- Do not use any abrasive cleaning implements to clean the inner cannula. Use only the cleaning brush provided. Do not use this brush to clean the tracheostomy tube.
- Blue Line Ultra™ Inner Cannula is designed and intended for single patient use. It should not be re-sterilised by any method.
- Dispose of per local procedure.

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WARNINGS

- Do not lubricate the inner cannula as this may occlude the inner cannula and prevent the inner cannula being retained within the tracheostomy tube.
- Contact with electrosurgery electrodes or laser surgery beams must be avoided because PVC will produce toxic fumes in air, or ignite in an enriched oxygen environment (e.g. anaesthesia).
- Never use a Portex® Fenestrated Tracheostomy tube for manual or mechanical ventilation unless a Portex® non fenestrated inner cannula is in place to prevent gas leakage through the fenestrations in the tracheostomy tube.
- If a Portex® Fenestrated Tracheostomy tube is used with a Decannulation cap, either a fenestrated inner cannula must be used or if a non fenestrated inner cannula is present, it should be removed, and in both cases the tracheostomy tube cuff must be deflated. Failure to do so may lead to partial or complete airway obstruction.
- Caution must be exercised when using a Decannulation Cap in conjunction with a cuffed fenestrated tracheostomy tube. The user must ensure that the cuff is completely deflated and that a fenestrated inner cannula is used, or that the inner cannula has been removed. Failure to do so may result in total or partial obstruction of the airway. Observe the patient for signs of respiratory distress and immediately remove the Decannulation Cap if required.
- If it is not possible to remove the inner cannula from the tracheostomy tube, do not try to forcibly remove. Both the inner cannula and tracheostomy tube should be removed together and replaced with a new tracheostomy tube and inner cannula.
- Ensure that the inner cannula is both the appropriate diameter and length for the tracheostomy tube in use. This can be checked by comparing size markings on the tracheostomy tube flange and patient labels with those on the inner cannula unit pack label. Use this inner cannula only with tracheostomy tubes from the Portex® Blue Line Ultra® range.
- Use of an inner cannula of the wrong diameter may cause difficulties in introduction or unnecessarily restrict gas flow. An inner cannula which is too long may protrude excessively from the outer tube, leading to tracheal damage or occlusion. Use of an inner cannula which is too short may lead to build-up of secretions which could cause infection or blockage.
- The use of an Inner Cannula reduces the airway diameter by 1.0mm (excluding size 10.0mm tracheostomy tube which reduces the airway diameter by 1.5mm).
- Do not fit an inner cannula when using a tracheostomy tube in conjunction with a right-angled connector that incorporates a co-axial fresh gas tube which protrudes into the connector lumen (e.g. Norman elbow). The protruding fresh gas tube could occlude the inner cannula, restricting expiratory flow of gas, resulting in barotrauma/hypoxia.
- During and after attachment of the breathing system to the tracheostomy tube connector, avoid application of excessive rotational or linear forces on the tube to prevent accidental disconnection or occlusion.

Portex, Blue Line Ultra, Thermoment, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

© 2015 Smiths Medical. All rights reserved.

	WORK INSTRUCTION FM-DVSOP2002-01 Rev. 001 Eff. Date: 25-OCT-2016 CO No.: CO-10029546	Doc. Number AS10002537-002	Revision 100
	WI Type: ARTWORK SPECIFICATION	Effective Date: 19 APR 2017	Page 6 of 19
Title: INSTRUCTIONS FOR USE - BLUE LINE ULTRA® TRACHEOSTOMY INNER CANNULA, SINGLE PATIENT USE (FORM FILL SEAL - SOFT PACK ONLY) (10002537-002)			

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Portex, Blue Line Ultra, Thermovent, designmerket Portex og designmerket Smiths Medical er varemerker som tilhører Smiths Medical. Symbolet © angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive eiere.