mylife™ Pura® BLOOD GLUCOSE TEST STRIP INSERT


EN

Intended use

The mylife™ Pura® blood glucose test strips are intended for the quantitative measurement of glucose in capillary whole blood from the fingertip, palm or forearm. The mylife™ Pura® blood glucose test strips are intended for *in vitro* (outside the body) diagnostics use for self testing. It is used as an aid in management of diabetes at home. The mylife™ Pura® blood glucose test strips must only be used with mylife™ Pura® / Pura X blood glucose monitoring system and mylife™ Control control solution. The application of other blood glucose test strips or control solutions might lead to incorrect results.

Test procedure

- Wash your hands with warm soapy water and dry thoroughly.
- Take one test strip from the vial. Recap the vial immediately.
- Insert the test strip into the strip port of the meter with the indication symbol facing up.
- A blood drop symbol appears and is accompanied by a further beep (if volume is switched on).
- Place the lancing device against the pad of your fingertip and press the release button.

Sample size example

0.75 µL	1.0 µL	1.5 µL	2.0 µL	3.0 µL
•	•	•	•	•

Please take a minimum of 0.75 µL of blood to perform the test with mylife™ Pura®. A blood sample size above 3.0 µL might contaminate the test strip port and the meter, a sample size below 0.75 µL may give an inaccurate test result.

- Touch and hold the blood drop to the edge of the test strip until you hear a beep (if volume is switched on) and the test strip view window is totally filled with blood. If the view window is not totally filled with blood or the test does not start, please discard the test strip and repeat the test with a new test strip.
- You will see the countdown mode on the screen. After 5 seconds, the test result will appear.
- Remove the test strip from the meter. Please follow the local regulations and discard the used test strip properly.

For more information on how to use your mylife™ Pura® blood glucose meter and mylife™ Softlance lancing device and questions regarding high and low results, see the user manual for the mylife™ Pura® blood glucose monitoring system.

Test result

Blood glucose test results are shown on the meter as mg/dL or mmol/L. // If your blood glucose result is unusually high or low, or if you question your results, repeat the test with a new test strip. You can run a quality control test with the control solution to check your meter and test strip. Quality control section please refer to the user manual for mylife™ Pura® system. // If your test result is below 10 mg/dL (0.6 mmol/L) or the test strip was not completely filled with blood, "Lo" will appear on the screen. Please repeat your test with a new test strip. // If your test result is above 600 mg/dL (33.3 mmol/L), "Hi" will appear on the screen. Please repeat your test with a new test strip. // If you repeat the test and still get "Hi" or "Lo" result, or you are experiencing symptoms that are not consistent with your blood glucose results and you have followed all the instructions in this manual, contact your healthcare professional immediately or follow their instruction for therapy.

Expected values⁽¹⁾

Fasting Blood Glucose

GLUCOSE LEVEL	INDICATION
From 70 to 99 mg/dL (3.9 to 5.5 mmol/L)	Normal fasting glucose
From 100 to 125 mg/dL (5.6 to 6.9 mmol/L)	Pre-diabetes (Impaired fasting glucose)
126 mg/dL (7.0 mmol/L) and above on more than one testing occasion	Diabetes

Precautions

Check the expiration date printed on the strip vial. Do not use expired test strips. // Close the vial cap immediately after taking the test strip out of the vial. // Do not perform a quality control test with expired control solution. // Do not bend or twist the test strip. Damaged test strips may result. // Do not reuse test strips. // Do not reuse lancets. // Discard used lancets result in an incorrect correctly. // If the meters and test strips are exposed to a high temperature difference, please wait 30 minutes before measurement.

Warning

Keep the test strips or vial cap away from children. They may cause a choking hazard. If a test strip or vial cap is swallowed, contact your physician immediately.

Limitations

The test strips are designed for use with capillary whole blood samples. Do not use serum or plasma samples. // Caution is advised in the interpretation of glucose values below 50 mg/dL (2.8 mmol/L) or above 250 mg/dL (13.9 mmol/L). Consult a physician as soon as possible or follow their instruction therapy for these situations, if values in this range are obtained. // Test results may be falsely low if the patient is severely dehydrated, excessive water loss, in shock, or in a hyperosmolar state (with or without ketosis). Critically ill patients should not be tested by blood glucose meters. // Do not perform the blood glucose test at temperatures below 10°C (50°F) or above 40°C (104°F) below 10% or above 90% relative humidity. // Hematocrit should be between 30% - 57%. If you do not know your hematocrit, ask your healthcare professional. // High blood fat values may influence the blood glucose test and patients should have baseline glucose values established by a clinical laboratory method prior to the start of home glucose monitoring with their physicians. // These baseline glucose values should be checked periodically. // Incorrect test results may be obtained at high altitude more than about 3,048 meters (10,000 feet) above sea level. // Fluoride should not be used as a preservative when collecting blood glucose samples. // The mylife™ Pura® has not been validated for use on neonates. Therefore, it should not be used for neonates.

NOTE

- Hands and fingers contaminated with sugar from foods or beverages may cause false elevated results.
- Storage of test strips near bleach as well as products containing bleach will affect the results of the test strips.
- Healthcare professionals should evaluate their test results and the measurement procedure at periodic intervals. It is recommended that blood glucose meter results should be compared to a well characterized laboratory measurement (hexokinase or glucose oxidase method) with the same blood sample.

Storage and handling

Store the test strips in the original capped vial at temperatures between 4°C to 30°C (39°F to 86°F) and relative humidity below 90%. // Do not freeze. // Replace the vial cap immediately and close tightly after taking the test strip out of the vial. Do not leave the cap of vial opened to avoid incorrect test results. // When new vial of test strips is first opened please write the opening date on the label. // Use test strips within 12 months after first opening or until the expiration date printed on the label (whichever comes first).

Measurement range

The measurement range of mylife™ Pura® is 10 to 600 mg/dL (0.6 to 33.3 mmol/L).

Accuracy

The accuracy of mylife™ Pura® was demonstrated by comparing plasma equivalent glucose values on the meter with plasma glucose values on a lab instrument.

Trouble shooting and customer service

For more information on error messages and trouble shooting, please refer to the error messages and trouble shooting section of the user manual. If you have any questions or in case of problems with the mylife™ Pura® / Pura X products, please contact your local Ypsomed customer service.

Additional information for healthcare professionals

Detection principle⁽²⁾

The glucose oxidase and potassium ferricyanide in the test strip reacts with the glucose in the sample to produce an electrical current which is proportional to the amount of glucose in the sample. The meter measures the current and converts it to the corresponding glucose concentration.

Precision

The precision was evaluated including (i) venous whole blood sample (ii) 3 levels glucose control solution in period of 10 days, by 10 meters and 3 batches of test strips.

(i) Venous whole blood sample: (ii) Control solution:

Glucose levels	P-01	P-02	P-03	P-04	P-05	CS-L	CS-N	CS-H
(1) Total test numbers (n)	300	300	300	300	300	300	300	300
(2) Mean mg/dL (mmol/L)	47.3(2.6)	105.6(5.9)	139.5(7.7)	234.9(13.0)	374.0(20.8)	46.7(2.6)	98.7(5.5)	289.9(16.1)
(3) SD mg/dL (mmol/L)	1.5(0.08)	2.0(0.11)	2.8(0.15)	4.3(0.24)	5.5(0.31)	1.2(0.06)	2.0(0.11)	4.1(0.23)
(4) CV (%)	3.2%	1.9%	2.0%	1.8%	1.5%	2.5%	2.0%	1.4%

Accuracy

The accuracy of the mylife™ Pura® system was demonstrated by comparing whole blood (plasma equivalent) glucose values on the mylife™ Pura® meter with plasma glucose values on a lab instrument. A total of 106 patients were enrolled. Each patient collected and tested their own blood samples (from the fingertip, palm and forearm) using the mylife™ Pura® system. Then the blood samples were centrifuged immediately after collection to obtain plasma. Analyze the plasma by the lab instrument - YSI 2300. 100% of meter values were within ± 15% of the YSI values at glucose concentrations ≥ 100 mg/dL (5.55 mmol/L) and within ± 15 mg/dL (0.83 mmol/L) at glucose concentrations < 100 mg/dL (5.55 mmol/L). The results and differences between the two methods, mylife™ Pura® system and YSI 2300 (as the reference method) are proved in the tables below.

Table 1: represents samples for glucose concentrations < 100 mg/dL (5.55 mmol/L).

Difference range in values between the YSI value and the meter value	The percent (and number) of samples for which the difference between the meter value (Alternative site) and the YSI value was within the difference range shown in the side row.
Within ± 5 mg/dL (0.28 mmol/L)	50.5%(94/186)
Within ± 10 mg/dL (0.56 mmol/L)	91.9%(171/186)
Within ± 15 mg/dL (0.83 mmol/L)	100%(186/186)

Table 2: represents samples for glucose concentrations ≥ 100 mg/dL (5.55 mmol/L).

Difference range in values between the YSI value and the meter value	The percent (and number) of samples for which the difference between the meter value (Alternative site) and the YSI value were within the difference range shown in the side row.
Within ± 5 mg/dL (0.28 mmol/L)	63.8%(287/450)
Within ± 10 mg/dL (0.56 mmol/L)	94.0%(423/450)
Within ± 15 mg/dL (0.83 mmol/L)	100%(450/450)

* Acceptance criteria in ISO 15197:2013 are that 95% of all differences in glucose values should be within ± 15 mg/dL (0.83 mmol/L) at glucose concentrations < 100 mg/dL (5.55 mmol/L), and within ± 15% at glucose concentrations ≥ 100 mg/dL (5.55 mmol/L).

Note: For glucose concentrations < 100 mg/dL (5.55 mmol/L), difference values are expressed in mg/dL (mmol/L), and for glucose concentrations ≥ 100 mg/dL (5.55 mmol/L), difference values are compared in percentage.

Lay User Evaluation

A total of 100 users were enrolled. Each user tested their fingertip blood samples with 3 lots of mylife™ Pura® test strips and the mylife™ Pura® meter. Professionally collected blood samples were centrifuged immediately after collection to obtain plasma and then analyzed by the lab instrument (YSI 2300 analyzer). 100% of the mylife™ Pura® values were within ± 15% of the YSI values at glucose concentrations ≥ 100 mg/dL (5.55 mmol/L) and within ± 15 mg/dL (0.83 mmol/L) at glucose concentrations < 100 mg/dL (5.55 mmol/L).

Interferences

The following compounds may interfere with the glucose measurement at the concentrations listed:

- Ascorbic acid ≥ 6 mg/dL (0.34 mmol/L)
- Glutathione reduced ≥ 70 mg/dL (2.28 mmol/L)
- Uric Acid ≥ 16 mg/dL (0.95 mmol/L)

Reagents

Each test strip contains the following reagents: Glucose Oxidase (GOD) 18.8% Non-reactive Ingredients 43.5%

References

- 1) Diabetes Information - American Association for Clinical Chemistry (AACC)[Electronic Version] Retrieved Aug 21, 2015 from www.labtestsonline.org/understanding/analytes/glucose/test.html
- 2) In vitro Diagnostics in Diabetes: Meeting the Challenge. Clinical Chemistry 45:9, 1596-1601 (1999).

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