en

Materials Used:

Acrylonitril Butadiene Styrene (ABS), Polyolefin alloy (PVC Free), Polypropylene (PP), Nylon, Polycarbonates (PC), High-density polyethylene (HDPE), Adhesive -Acrylate.

Intended Purpose:

Transfusion sets are single-use devices intended to administer blood from a container into the patient through a needle or catheter inserted into a vein.

Intended Users:

Healthcare professionals.

There is no restriction to a specific patient population.

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
NEUTRAL CODE NO: OM20000

OneMed Group Oy
Metsäntie 20, FI-00320 Helsinki, Finland

- Re-use may lead to infection or other illness/injury.

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Instructions for Use:

1. Ensure aseptic technique is followed at all times.

2. Visually inspect the product for any imperfections or surface deterioration before use.

3. For safety reasons, tighten the caps and all disconnectable fittings before use.

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fins prévues :

• Re-use may lead to infection or other illness/injury.

Indication:

Transfusion of blood and blood component.
Maintaining blood requirement during surgery and excessive blood loss.

Contraindication:

Not to be used in patients with known hypersensitivity to any of the materials used.

Disposal:

According to local regulations.

Notice:

Serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Instructions for Use:

1. Ensure aseptic technique is followed at all times.
2. Visually inspect the product for any imperfections or surface deterioration before use.
3. For safety reasons, tighten the caps and all disconnectable fittings before use.

4. Move the roller clamp from its position to inspect the tubing for any abnormalities and verify that the air-vent is closed in vented devices.

5. Close the roller clamp and remove the spike cover. If Transfusion set is Y-type close C-clamp on non-used tube.

6. Place the blood/ blood component container horizontal on the tabletop and insert 3/4 of its length of the spike horizontally into it using 1/4 turn motion.

7. Invert the container and hang it on the IV pole or similar.

8. Attached the tube to the tube holder on the roller clamp.

9. Squeeze the drip chamber gentle 3-5 times and allow the solution to fill half of the drip chamber.

10. Open the roller clamp and fill the tubing with the blood/ blood component until all air is gone from the tubing.

11. Close roller clamp and remove priming cap and connect the male connector to the vascular access device.

12. Set the drip rate with the roller clamp as per facility protocol.

13. Inspect and secure the connections continuously throughout the procedure.

14. Recommended to check the drip rate at regular interval.

Warning:

- The use of this product is restricted to a qualified doctor or nurse.
- Read the instructions before use. The device should be used in accordance with the instructions.
- The device is guaranteed sterile and pyrogen-free if the individual package has not been opened or damaged. Do not use if package is damaged or if protective caps are loose or missing.
- Do not over-tighten a luer beyond the friction fit as this may damage both the luer and another connecting device.
- Do not use any instruments for clamping the tubing or luer tightening.
- Den Behälter umdrehen und an der IV-Stange oder

Indication:

Transfusion of blood and blood component.
Maintaining blood requirement during surgery and excessive blood loss.

Contraindication:

Not to be used in patients with known hypersensitivity to any of the materials used.

Disposal:

According to local regulations.

Notice:

Serious incident that has occurred in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Instructions for Use:

1. Ensure aseptic technique is followed at all times.
2. Visually inspect the product for any imperfections or surface deterioration before use.
3. For safety reasons, tighten the caps and all disconnectable fittings before use.

4. Move the roller clamp from its position to inspect the tubing for any abnormalities and verify that the air-vent is closed in vented devices.

5. Close the roller clamp and remove the spike cover. If Transfusion set is Y-type close C-clamp on non-used tube.

6. Place the blood/ blood component container horizontal on the tabletop and insert 3/4 of its length of the spike horizontally into it using 1/4 turn motion.

7. Invert the container and hang it on the IV pole or similar.

8. Attached the tube to the tube holder on the roller clamp.

9. Squeeze the drip chamber gentle 3-5 times and allow the solution to fill half of the drip chamber.

10. Open the roller clamp and fill the tubing with the blood/ blood component until all air is gone from the tubing.

11. Close roller clamp and remove priming cap and connect the male connector to the vascular access device.

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