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Cuticell® Contact

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Soft silicone wound contact layer

Product description

Cuticell® Contact is a primary wound contact layer consisting of an elastic and transparent polyurethane film coated with soft silicone. Cuticell® Contact is designed to adhere to the surrounding skin but not to the moist wound bed and to allow exudate to pass through into an outer absorbent dressing. Cuticell® Contact is highly conformable with two-way stretch.

Indications

Cuticell® Contact is indicated for acute and chronic wounds such as: ulcers of the lower legs and decubitus, surgical wounds, blisters, second degree burns, skin graft donor sites, lacerations, abscesses, skin tears. Cuticell® Contact may be used as a wound contact layer in conjunction with Negative Pressure Wound Therapy (NPWT).

Contraindications

Do not use on patients with a known sensitivity to any of the product components.

Warnings and precautions

Cuticell® Contact is packaged for single use. Do not re-use or re-sterilize as this is a risk of transmission of body fluids or contaminated tissue between patients. Do not use if the pouch is open or damaged. Discard open or unused material of the device. The wound should be inspected for signs of infection and -in case- treated according to clinical practice, untersucht und bei Bedarf gemäß den Regeln Mühscher Praxis behandelt werden.

Directions for use

- Clean the wound according to good clinical practice.
- Make sure the surrounding skin is dry.
- Choose a dressing size that covers the wound and overlaps the wound edges for at least 2 cm.
- Open the pouch and take out Cuticell® Contact. Ensure sterility.
- Remove one release liner (1).
- Apply Cuticell® Contact with the sticky side towards the wound (2), remove second release liner (3) and smooth dressing in place.
- Cover with a secondary dressing, appropriate for the specific amount of exudate.
- Frequency of dressing change depends on the condition of the wound. It might be sufficient to change only the secondary dressing.
- Product can be left in place for up to 14 days, depending on the condition of the wound and taking local clinical protocols into consideration.
- To change Cuticell® Contact, carefully remove the dressing, peel off the adhesive side of the product and discard to application of a new dressing.
- Clean the wound according to good clinical practice prior to application of a new dressing.
- Product may be used in conjunction with topical treatments, like antiseptics or hydrogels.

Storage

Cuticell® Contact should be stored in dry conditions. Do not expose to direct sunlight.

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BSN medical GmbH