

curea medical GmbH
Münsterstraße 61-65
48565 Steinfurt
info@curea-medical.de
+49 36071 9009500

REV 02.1/2020-08

curea clean/curea clean BTBS

ENGLISH

1) Description

curea clean/curea clean BTBS is a non-sterile polyacrylate-based wound dressing and class 1 medical device [MD].

2) Indication

Absorption of wound exudate as a secondary absorbent pad.

3) Contraindication

- Known hypersensitivity to the components, e.g. Acrylic acid.
- Do not use on mucous membranes, eyes or exposed bones or tendons, as these parts of the body may be permanently damaged.
- Heavily bleeding wounds because the product is not a hemostat.

4) Application

Make sure to cover the wound with the white side of the dressing. Apply the product under medical supervision and according to the indication. We recommend to change the dressing at least after 7 days and a frequent (best daily) check of the periwound area for signs of infection and of the dressing for depletion of capacity.

5) Residual risk

The application of the products is generally very low risk when used as intended and according established therapeutic standards, yet there are few residual risks:

- The swelling of the wound dressing can cause fixation bandages to cause tension bubbles.
- Saturated wound dressings can lead to maceration (disintegration) of the skin.
- Saturated wound dressings can cause circulatory disorder due to pressure and weight.
- Infections could be overseen due to long application times.
- Wound dressings can absorb larger amounts of blood, depending on their size, and cause critical health status.

6) Side effects/obligation to report

There are no known physiological side effects. The use of the products is not limited to any specific population group. If any side effects occur in connection with the use of this product, please inform us informally via our mailing address or by phone +49 36071 9009 500 or by email to info@curea-medical.de

If the side effects were so serious that one of the following consequences occurred directly or indirectly, could have occurred or could occur in the future:

- death of a patient, user or other person,
- a temporary or permanent serious deterioration in the health of a patient, user or other person,

- a serious threat to public health, report this immediately to the relevant regulatory authority and the manufacturer.

7) Warnings

- Do not tear or cut the wound dressing as core material may leak out and contaminate the patient. If core material leaks, clean the wound with sterile wound irrigation. There is no danger to the patient or user due to spilled sodium polyacrylate.
- The wound dressings are intended for single use. Do not reuse to avoid potential infection risks.

8) Storage

Store in a dry place and protect from high humidity. Protect from direct sun light.
Regard temperature limits.

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