
SPIRO®

REF 710 CE 0402

REF 710R CE 0402

REF 711 CE

UNCONTROLLED

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en Instructions for Use

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


no Bruksanvisning

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sv Bruksanvisning

FOGLESS
International

CE LOT REF Rx Only

| | | |
|----|--------------------------|----|
| da | Brugervejledning | 2 |
| de | Gebrauchsanweisung | 4 |
| en | Instructions for Use | 6 |
| es | Instrucciones de uso | 8 |
| fi | Käyttöohjeet | 10 |
| fr | Mode d'emploi | 12 |
| it | Istruzioni per l'uso | 14 |
| nl | Gebruiksaanwijzing | 16 |
| no | Bruksanvisning | 18 |
| pl | Instrukcja obsługi | 20 |
| pt | Instruções de utilização | 22 |
| sv | Bruksanvisning | 24 |

SPIRO TALEVENTIL **REF** 710 + 710R + 711

Kun en-patients brug

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SPIRO SPEAKING VALVE **REF** 710 + 710R + 711

single patient use only

The SPEAKING VALVE is a one-way valve which opens ONLY FOR INHALATION. The exhalation passages MUST BE OPEN. Ensure that the patient can tolerate an uncuffed cannula. Measurement of the tracheal pressure at passive exhalation with the Speaking Valve in place is advised (< 5 cm H₂O is recommended).

N.B. THE SPEAKING VALVE MUST NEVER BE USED IN COMBINATION WITH A CUFFED CANNULA WITHOUT FENESTRATION, AS THE PATIENT WILL BE UNABLE TO EXHALE!

CONTRAINDICATIONS: This device is not for use with patients who have respiratory infection, extensive phlegm, stenosis, vocal cord paralysis or acute airway obstructions. This device is not for use with inline ventilators, endotracheal tubes or other artificial airways.

The COVER (5) is provided with a tapered section which is seated in the valve HOUSING (3) during use. The further pressed in, the tighter the cover is secured. Avoid pressing the cover too hard, as it reduces the valve's life and makes it difficult to remove.

DAILY CHECKS: Check valve operation daily and inspect for damage. Replace valve with new unit in the case of damage or defective operation. **NEVER DISASSEMBLE THE MEMBRANE!**

CLEAN the valve daily. Do not use sharp tools. Use the fingers or **PLASTIC TWEEZERS** (P/N 705, sold separately) to remove/re-insert the filter. The valve housing and cover shall be air dried after washing with luke-warm water and mild fragrance-free detergent. Rinse the valve cover in both directions to remove all cleaning liquids. Other detergents or solvents must not be used. Always check the operation again after cleaning and check that the air passages are free. Do not wash and re-use the filter!

The valve should always be used with **FILTER (4)** (P/N 740, sold separately). **THE FILTER SHOULD BE REPLACED** daily or more frequently if required. Ensure that the filter is seated evenly against the bottom of the valve housing, and that air flow through the filter is not obstructed by phlegm or other substances.

A HEAT SHIELD/SHOWER PROTECTOR (6) (P/N 730, sold separately) should be used as required with the inhalation opening facing downward.

An ADAPTER (2) (P/N 707, 708, sold separately) is only required for connection to a cannula (1) without an ISO 15 connector (taper).

SUPPLEMENTAL OXYGEN can connect via the Ø5.5 mm nipple on the valve housing. The WHITE O₂ NIPPLE CAP (7) (P/N 710/710R only) shall be fitted to the Speaking Valve and removed only when oxygen is being supplied.

The SPEAKING VALVE should be replaced with a new unit after 29 days of use, or more frequently if necessary. The valve must not be used while the patient is sleeping or otherwise when the patient's consciousness is impaired. Use a SPIRO HUMIDIFIER (P/N 702, sold separately) in such cases.

The Speaking Valve must be removed before inhalation of medicine.

RESISTANCE TO FLOW for SPIRO Speaking Valve is < 2 cm H₂O*, and with the heat shield/shower protection on < 3 cm H₂O*.

All products should be stored in a dark and dry environment and at normal indoor temperature.

For the STORAGE PERIOD OF SPIRO PRODUCTS: See year and month next to the 'HOURGLASS' symbol on the product label. The products shall be disposed of in accordance with the applicable national regulations.

ALWAYS CONSULT YOUR DOCTOR, SUPPLIER OR MANUFACTURER IF YOU HAVE ANY QUESTIONS.

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Fogless International AB