Instructions for Use

## Cannulated Bone Tunnel Plug

EN English

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TR Türkçe

Smith & Nephew, Inc.  
Endoscopy  
150 Minuteman Road  
Andover, MA 01810  
USA

[www.smith-nephew.com](http://www.smith-nephew.com)  
[ifu.smith-nephew.com](mailto:ifu.smith-nephew.com)

+1 978 749 1000 Telephone  
+1 978 749 1108 Fax  
+1 800 343 5717 U.S. Customer Service

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©Trademark of Smith & Nephew, registered U.S. Patent & Trademark Office. This product may be covered by one or more U.S. patents. See [smith-nephew.com/patents](http://smith-nephew.com/patents) for details.  
©1997, 2024 Smith & Nephew, Inc. All rights reserved.

01/2024 1060440 Rev. K

# Cannulated Bone Tunnel Plug

## Instructions for Use

### Device Description

The Smith & Nephew Cannulated Bone Tunnel Plug is designed for maintaining fluid and joint distention during insertion, and extraction of instruments.

The three-rib, tapered construction accommodates tunnel diameters of 7 mm to 12 mm.

The internal diameter tapers from 7 mm to 5 mm. The internal diameter may be customized, as instructed below, to accommodate the insertion of a range of instrument sizes.

### Performance Characteristics

The Smith & Nephew Cannulated Bone Tunnel Plug has a double seal design capable of preventing fluid leakage.

### How Supplied

The Smith & Nephew Cannulated Bone Tunnel Plug is supplied sterile, for single use only.

### Intended Users

The Cannulated Bone Tunnel Plug is intended for use by health care professionals in accordance with these instructions for use. The use environment is a professional healthcare facility.

### Intended Use

The Cannulated Bone Tunnel Plug is intended for arthroscopic surgical use for maintaining fluid and joint distension during insertion and extraction of instruments.

### Indications for Use

The Cannulated Bone Tunnel Plug is indicated for arthroscopic surgical use for maintaining fluid and joint distension during insertion and extraction of instruments.

### Contraindications

- Known hypersensitivity to the material. Where sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to use.
- Insufficient quantity or quality of bone.
- Insufficient blood supply and previous infections which may tend to retard healing.
- Active infection.

- Conditions which tended to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.

### Target Patient Groups

The patient population for the Cannulated Bone Tunnel Plug are patients in whom arthroscopic surgery is determined necessary by the prescribing healthcare professional.

### Warnings

- **Do not use if the product sterilization barrier or its packaging is compromised.**
- **The Cannulated Bone Tunnel Plug is supplied sterile and intended for SINGLE USE ONLY. DO NOT clean, resterilize, or reuse the Cannulated Bone Tunnel Plug as this may result in product malfunction, failure, or patient injury and may expose the patient to the risk of transmitting infectious diseases. Do not use after the expiration date.**
- **It is the healthcare professional's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.**
- **Read these instructions completely prior to use.**

### Precautions

**R<sub>x</sub> only** U.S. Federal law restricts this device to sale by or on the order of a physician.


- Prior to use, inspect the device to ensure it is not damaged. Do not use a damaged device.
- After use, this device may be a potential biohazard and should be handled in accordance with accepted medical practice and applicable local and national requirements.
- As with any surgical instrument, careful attention should be paid to assure that excessive force is not placed on the instrument. Excessive forces applied to the instrument can result in failure.

### Adverse Event Reporting

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

complaints@smith-nephew.com

## Instructions for Use

(Do not cut beyond this point)

Figure 1.

### **To Customize Inner Diameter**

- Referring to Figure 1, cut off the tip of the Cannulated Bone Tunnel Plug at the length which corresponds to the internal diameter required for passing instruments. Surgical scissors or scalpel may be used to cut the device.
- Do not cut beyond the third rib.

### **To Insert**

- Insert the Cannulated Bone Tunnel Plug into bone tunnel until a snug fit is accomplished.

## Safe Disposal

After use, this single use device may be a potential biohazard and should be handled in accordance with accepted medical practice and applicable local and national requirements.

## For Further Information

If further information, including warranty, on this product is needed, please contact Smith & Nephew Customer Service at +1 800 343 5717 in the U.S., or an authorized representative.

## Symbols Glossary

	Manufacturer		CE Mark		Consult instructions for use
	Date of Manufacture		Do not use if package is damaged. Consult instructions for use.		Prescription use only
	Use-by date		Sterilized using irradiation		Medical Device
	Catalog Number		Do not re-sterilize		Authorized Representative in the European Community
	Batch Code		Do not re-use		European Importer
	Unique Device Identifier		Single sterile barrier system		
	Caution				

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