

evercare[®] inLine Urological set, Y-type, PVC-free

REF: 10555

evercare[®] inLine Urology sets are used to conduct fluids through gravity method, by enabling connection to a fluid container to carry out washing and/or cleaning of urethra or bladder through the irrigation tube connected to a final medical device directly in contact with patient body. PVC-free Urology sets can be used for irrigation and for cystoscopy, T.U.R and other arthroscopic surgical procedures. Y-type make it possibility to have a continuous flow and use the same set for continuous irrigation according to internal procedure or guidelines..

Key Features

- ▶ Large bore set with long tubing
- ▶ C-Clamp to stop the flow of the infusion
- ▶ Transparent, Soft & larg drip chamber
- ▶ For easy squeezing and filling
- ▶ For safe high flow rate
- ▶ To minimize the risk of air pockets within the tubing
- ▶ White Roller clamp
- ▶ Three different connectors for all kind of applications
- ▶ 15 cm silicone tubing
- ▶ Standard cone
- ▶ Luer lock adaptor (bi-packed)
- ▶ A secondary sterile packaging with a primary protective packaging inside & can be used to open within a sterile field
- ▶ Y-type with 2 spikes
- ▶ 2-way
- ▶ 2 C-clamp



REF	DESCRIPTION	SIZE - Tube Inner/outer, L(length)	PRIMING VOLUME(ml) Approx	Pack Count dispenser/carton
10555	Urological set, Y-type, PVC-free	Ø 4.8/6.8 mm, L 215cm	NA	5 / 20

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-01
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class IIa
Compliance with product standards:	EN 8536-4, EN ISO 80369-7
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Spike cover	PP, transparent
Spike	ABS
C-clamps	PP
Y-connector	Polyolefin Alloy, transparent, PVC-free
Tube	Polyolefin Alloy, transparent, kink resistant, PVC-free
Drip Chamber cap	ABS
Drip Chamber	Polyolefin Alloy, transparent, PVC-free
Roller clamp	ABS, roller wheel white
Stepped connector	Polyolefin alloy, transparent, PVC-free
Tube Silicone	Silicone
Luer lock adaptor	ABS

Warnings & Recommendations

Warnings

The use of this product is restricted to a qualified doctor or nurse.
Read the instructions before use. The device should be used in accordance with the instructions.
Do not use for IV infusions or blood transfusion.
The device is guaranteed sterile and pyrogen-free if the individual package has not been opened or damaged. Do not use if package is damaged or if protective caps are loose or missing.
Do not over-tighten a luer beyond the friction fit as this may damage both the luer and another connecting device.
Do not use any instruments for clamping the tubing or luer tightening.
Re-use may lead to infection or other illness/injury.

Storage Recommendation

Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.

Disposal Recommendation

Dispose in accordance with local regulations.

Shelf Life

5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PE film	354 x 200
Dispenser Pack	PP	305 x 460
Transport Carton	TCF quality, Carton	400 x 295 x 287
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
10555	06438129117516	N/A	06438129517514

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.