

evercare® inLine Gravity set, Vented, 3-way stopcock glued, PVC-free

REF: 10577-01, 10578-01

evercare[®] inLine Gravity Sets are sterile, single use products that are used to conduct intravenous fluids into the vascular system of a patient by connecting to a source of infusion in one end and to an infusion device with vascular access in the other end. The fluid is transferred by the means of gravity and the flow rate is regulated by a roller clamp component. Suitable for glass bottles, glass vials and plastic bags or bottles.

Key Features

- ▶ Transparent soft drip chamber for easy squeezing and filling
- ▶ 15-micron particle filter at bottom of the drip chamber
- Hydrophobic filter prevents leakage while priming
- ▶ 3-way stopcock
- Lipid resistant
- ▶ 45° click feature
- Can withstand up to 2 bar
- ▶ Glued
- ▶ Rotating luer lock



REF	DESCRIPTION	SIZE - Tube Inner/outer, L(length)	PRIMING VOLUME(ml) Approx	Pack Count dispenser/carton
10577-01	Gravity set, Vented, 3-way stopcock glued, PVC-free	Ø 3/4.1 mm, L 195cm	12.9	20 / 80
10578-01	Gravity set, Vented, 3-way stopcock glued, PVC-free	Ø 3/4.1 mm, L 335cm	22	15 / 60

Sales Unit: Transport Carton

Page 1 2024-06-13

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:

The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-01

Compliance with legal requirements: Medical Device Regulation (EU) 2017/745, MDR Class IIa

Compliance with product standards: EN 8536-8, EN ISO 80369-7

Quality standards: EN ISO 13485

Label information and packaging compliance: EN ISO 15223-1, EN ISO 20417

Material

Raw Material & Characteristics
PP, transparent
ABS, transparent blue
LDPE/HDPE, blue
PP + hydrophobic filter
Polyolefin alloy, PVC-free
PP/nylon, pore size 15 μm
Polyolefin alloy, PVC-free
ABS, roller body blue, roller wheel white
ABS, transparent
Transparent channel PC, lipid resistant, white handle HDPE
HDPE, white
ABS & PC, transparent
HDPE, light blue

Page 2 2024-06-13

Warnings & Recommendations

TI CILL			1. 6.	
The use of this	product is	restricted to a	a dualified	doctor or nurse.

Read the instructions before use. The device should be used in accordance with the instructions.

Do not use for transfusion of blood or blood components.

The device is guaranteed sterile and pyrogen-free if the individual package has not been opened or damaged. Do not use if package is damaged or if protective caps are loose or missing.

Do not over-tighten a luer beyond the friction fit as this may damage both the luer and another

connecting device.

Do not use any instruments for clamping the tubing or luer tightening.

Re-use may lead to infection or other illness/injury.

Storage Recommendation	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PE film	235 x 100
Dispenser Pack	PP	305 x 460 x 0
Transport Carton	TCF quality, Carton	400 x 295 x 287
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
10577-01	06438129122169	N/A	06438129522167
10578-01	06438129122176	N/A	06438129522174

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.