equine Progesterone VET test kit

English

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of equine progesterone in li-hep plasma or serum on a solo or CUBE-VET analyser.

Order information	Indication	Kit size
Order number: VT0250	equine Prog. test kit	6 tests
Order number: VT2500	equine Prog. control kit	1 x 2 ml (dec. level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

During early pregnancy, progesterone is produced in the equine ovary by the corpus luteum (CL). Its concentration remains elevated and peaks between 60 and 120 days of gestation. Circulating progesterone is used diagnostically to evaluate luteal function during early pregnancy. When the circulation progesterone concentration reaches a distinct level, this can be considered adequate to maintain early stage pregnancy. There are a number of reasons for monitoring during pregnancy, such as uterine infections, history of pregnancy loss and luteal insufficiency.

Method

Homogeneous immunoturbidimetric test.

Measurement Range

1.4 – 13.0 ng/ml

Sample Material

Use 40 µl li-hep plasma or serum.

- Use fresh sample material and perform the measurement within 2 hours.
- Storage up to 24 hours: Separate serum/plasma from the cells and store it in a plain tube at 2 – 8 °C.
- Storage up to 2 weeks: Separate serum/plasma from the cells and store it in a plain tube at -20 °C

Test Kit

ERS cuvette filled with antibody reagent
ERS cap filled with buffer reagent

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 – 8 °C. Opening the vacuum package may limit the reagent stability to three months (stored at 2 – 8 °C) from the date of opening. DO NOT FREEZE!

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Decision Limits

Horse: > 2 ng/ml (gestation)

It is recommended that each laboratory establishes its own decision limits.

Quality Control

For internal quality control the EUROLyser equine Progesterone control kit is recommended. Order number: VT2500

Precision

Reproducibility within-run: Control; N = 20; mean = 5.4 ng/ml; CV = 9.69%;
Precision accepted: CV < 12%

Correlation

Equine sample correlation: N = 50
 y (EUROLyser eProg) = 1.0055x (Siemens Progesterone) - 0.0119; R^2 = 0.9522;

Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 15% of initial values.

Haemoglobin	525 mg/dl	Cholesterol	620 mg/dl
Human albumin	12 g/dl	Rheumatoid factor	1080 IU/ml
Bilirubin (conjugated)	72 mg/dl	Triglycerides	835 mg/dl
Bilirubin (unconjugated)	72 mg/dl	Uric acid	30 mg/dl

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


English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

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2. Sample preparation

- 2.1 Aspirate 40 µl sample material from centrifuged sample tube

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3. Sample processing

- 3.1 Dispense 40 µl sample INTO LIQUID in ERS cuvette
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic sample processing by pressing the  start button on the solo laboratory photometer or by closing the door of the CUBE-VET laboratory photometer

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