



Breathing Filters, HMEs and HMEFs



Use of Breathing Filters

Breathing filters provide an effective barrier that prevent cross contamination between patients, respiratory breathing systems, equipment and the clinical environment. Their use is widely recognised as beneficial and is recommended by a number of Anaesthetic Associations¹.

The threat to patients


Patients who require an artificial airway have their natural physiological protection bypassed. This will increase the risk of cross contamination between patients and healthcare equipment. The cross contamination of patients via an anaesthetic system has been reported, and documented areas of concern regarding infection includes *Hepatitis C*, *Mycobacterium tuberculosis*, blood in sputum and the *SARS* virus.

Critically ill patients are commonly at risk of infection and particularly from Ventilator Associated Pneumonia (V.A.P.). This nosocomial infection increases morbidity and potential mortality as well as the cost of treating the patient. The strategic use of an efficient breathing filter will provide an effective barrier between patients, breathing systems and ventilatory equipment.


Proven efficiency

The Intersurgical range of breathing filters has been designed for the protection of the patient, breathing system and equipment. They have been independently tested and proven to be highly efficient in preventing the passage of bacteria and viruses. Clinically relevant testing is carried out on all products using *Bacillus subtilis* (1.0µm x 0.7µm) and Ø174 bacteriophage, additional testing includes *Mycobacterium tuberculosis* (0.3µm x 1.0µm), *Hepatitis C* (0.03µm) and *MS-2 coliphage* (0.02µm). These tests provide you with clinically relevant information to allow evidence-based decisions to be made on the most appropriate product to meet your clinical requirements.

Potential infectious viruses [Particle sizes µ microns]




Coliphage T1 [0.017µ]
MS-2 coliphage [0.02µ]
Hepatitis C [0.03µ]
Adenovirus [0.07µ]




HIV [0.11µ]
Cytomegalovirus (CMV) [0.1µ]
Orthomyxovirus [0.1µ]

Potential infectious bacteria [Particle sizes µ microns]



Mycobacterium tuberculosis
[0.3µ x 1.0µ smallest size]
Serratia marcescens [0.45µ]
Pseudomonas aeruginosa [0.5µ]
Brevundimonas diminuta [0.3µ]



Staphylococcus aureus [1.0µ]
Bacillus subtilis [1.0µ x 0.7µ]

Essential requirements

The Intersurgical range of breathing filters offers a choice of electrostatic and pleated mechanical filters with a range of patient connections, providing a choice of products to meet various clinical situations.

- Independently validated Filtration Efficiency²
- Proven filtration against *Mycobacterium tuberculosis* and *Hepatitis C*
- Proven efficiency not affected by anaesthetic agent
- Safe inert material
- Option of patient connections – conveniently packed and ready for use
- Lightweight – reducing patient trauma
- Low compressible volume – reducing rebreathing of CO₂
- Low resistance to flow – over 24 hours
- Safe, secure ISO connectors
- Compliance to all relevant international standards

Filtration efficiency

Filtration performance is determined by independent microbiological testing against clinically relevant bacterial and viral challenges. The level of breakthrough of the challenge determines the efficiency. This efficiency is reported as a percentage based upon this break-through.

Number of organisms challenging the filter	Number of organisms passing through the filter	Efficiency of the filter
100,000	1,000	99%
	100	99.9%
	10	99.99%
	1	99.999%

References

1. Association of Anaesthetists of Great Britain and Ireland 1996. Danish Society of Anaesthetists 1998. French Society of Anaesthetists 1998.
2. All filters are independently validated for filtration efficiency at the Health Protection Agency, Porton Down, Salisbury, Wiltshire, U.K and Nelson Laboratories Inc, USA. All quoted performance figures are mean values.

Filta-Guard™ range

High-efficiency

A dedicated breathing filter designed for use in breathing systems in anaesthesia and intensive care, the high-efficiency Filta-Guard protects the patient, hospital personnel and equipment from potential microbial contamination. The Flow diffuser improves performance and optimises resistance to flow.



Code	1944000	1944003	1944011
Box Qty.	70	70	20
Luer lock port		✓	✓
Filtration efficiency	>99.999%	>99.999%	>99.999%
Resistance to flow at 30L/min	1.0cm H ₂ O	1.0cm H ₂ O	1.1cm H ₂ O
Resistance to flow at 60L/min	2.3cm H ₂ O	2.48cm H ₂ O	2.8cm H ₂ O
Compressible volume	67ml	67ml	67ml*
Weight	40g	41g	54g
Connections	22F-22M/15F	22F-22M/15F	22F-22M/15F
Minimum tidal volume	>200ml	>200ml	>200ml
Accessories			Superset™ catheter mount** 3504000

Inter-Guard™ range

Sterile

The Inter-Guard range of sterile breathing filters is designed for use in breathing systems in the operating theatre and the intensive care unit for the protection of the patient, breathing system and equipment. The perfect combination between size and performance.



Code	1344007S S	1344000S S	1344711S S
Box Qty.	50	50	50
Luer lock port		✓	
Filtration efficiency	>99.998%	>99.998%	>99.998%
Resistance to flow at 30L/min	0.8cm H ₂ O	0.8cm H ₂ O	1.1cm H ₂ O
Resistance to flow at 60L/min	2.0cm H ₂ O	2.0cm H ₂ O	2.7cm H ₂ O
Compressible volume	41ml	42ml	42ml*
Weight	22g	23g	36g
Connectors	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F
Minimum tidal volume	>150ml	>150ml	>150ml
Accessories			Superset™ catheter mount** 3504000

S Sterile

* Filter only

** Refer to information on catheter mount volumes in Patient Connections.

Clear-Guard™ range

Medium efficiency

The Clear-Guard range of breathing filters includes a number of options, all designed for use in breathing and anaesthetic systems for the protection of the patient, hospital personnel and the equipment from potential microbial contamination. Designed with a rounded ergonomic polypropylene housing, the Clear-Guard 3 range represents our most cost-effective filter option, and is available with an integral 90° elbow, reducing the need for an additional catheter mount or separate patient elbow.

The Clear-Guard Midi low volume filter provides a further option with minimum deadspace; ideal for use in anaesthesia.

Clear-Guard 3



Code	1544007	1544000	1544011	1544012
Box Qty.	150	150	75	75
Luer lock port		✓	✓	✓
Filtration efficiency	>99.99%	>99.99%	>99.99%	>99.99%
Resistance to flow at 30L/min	0.9cm H ₂ O	0.9cm H ₂ O	0.9cm H ₂ O	1.2cm H ₂ O
Resistance to flow at 60L/min	1.9cm H ₂ O	1.9cm H ₂ O	2.5cm H ₂ O	3.3cm H ₂ O
Compressible volume	60ml	60ml	60ml*	60ml*
Weight	27g	28g	42g	43g
Connectors	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F
Minimum tidal volume	>200ml	>200ml	>200ml	>200ml
Accessories			Superset™ catheter mount** 3504000	Superset™ catheter mount with fixed elbow** 3502000



Code	1544197 NEW	1544351	1545000
Box Qty.	50	70	75
Luer lock port	✓	✓	✓
Filtration efficiency	>99.99%	>99.99%	>99.99%
Resistance to flow at 30L/min	1.0cm H ₂ O	1.2cm H ₂ O	0.8cm H ₂ O
Resistance to flow at 60L/min	3.4cm H ₂ O	3.6cm H ₂ O	2.1cm H ₂ O
Compressible volume	60ml*	60ml*	75ml
Weight	34g	45g	34g
Connectors	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F
Minimum tidal volume	>200ml	>200ml	>225ml
Accessories	Fixed elbow 1997000	Flexible catheter mount with fixed elbow** 3501000	

* Filter only

** Refer to information on catheter mount volumes in Patient Connections.

Clear-Guard™ range

Clear-Guard Midi



Code	1644000	1644007 <i>NEW</i>	1644011	1644012
Box Qty.	100	100	75	75
Luer lock port	✓	✓	✓	✓
Filtration efficiency	>99.9%	>99.9%	>99.9%	>99.9%
Resistance to flow at 30L/min	0.7cm H ₂ O	0.7cm H ₂ O	0.9cm H ₂ O	1.0cm H ₂ O
Resistance to flow at 60L/min	1.8cm H ₂ O	1.8cm H ₂ O	2.2cm H ₂ O	2.8cm H ₂ O
Compressible volume	34ml	34ml	34ml*	34ml*
Weight	19g	19g	33g	33g
Connectors	22F-22M/15F	22F-22M/15F	22F-22M/15F	22F-22M/15F
Minimum tidal volume	>100ml	>100ml	>100ml	>100ml
Accessories			Superset™ catheter mount** 3504000	Superset™ catheter mount with fixed elbow** 3502000



Clear-Guard Midi

Code	1644137
Box Qty.	50
Luer lock port	
Filtration efficiency	>99.9%
Resistance to flow at 30L/min	1.0cm H ₂ O
Resistance to flow at 60L/min	3.0cm H ₂ O
Compressible volume	34ml*
Weight	25g
Connectors	22F - mouthpiece
Minimum tidal volume	>100ml
Accessories	mouthpiece

* Filter only

** Refer to information on catheter mount volumes in Patient Connections.

Hydro-Guard™ range

Pleated membrane filter

A versatile, low-volume breathing filter with a pleated mechanical membrane for use in anaesthesia as an HMEF (moisture return: 23mg H₂O/L at VT 500ml), or in ITU as a filter only. Available as a sterile option if required.



Code	1745000	1745011	1745012	1745197
Box Qty.	40	75	75	50
Luer lock port	✓	✓	✓	✓
Filtration efficiency	>99.999%	>99.999%	>99.999%	>99.999%
Resistance to flow at 30L/min	1.3cm H ₂ O	1.4cm H ₂ O	1.7cm H ₂ O	1.6cm H ₂ O
Resistance to flow at 60L/min	2.9cm H ₂ O	3.4cm H ₂ O	4.4cm H ₂ O	3.9cm H ₂ O
Compressible volume	63ml	63ml*	63ml*	63ml*
Weight	30g	44g	45g	37g
Connectors	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F
Minimum tidal volume	>200ml	>200ml	>200ml	>200ml
Accessories		Superset™ catheter mount** 3504000	Superset™ catheter mount with fixed elbow** 3502000	fixed elbow** 1997000



Code	1745711	1745180 NEW	1744000S	1744011S	1744580S
Box Qty.	75	50	50	50	50
Luer lock port	✓	✓	✓	✓	✓
Filtration efficiency	>99.999%	>99.999%	>99.999%	>99.999%	>99.999%
Resistance to flow at 30L/min	1.7cm H ₂ O	1.5cm H ₂ O	1.7cm H ₂ O	1.6cm H ₂ O	1.73cm H ₂ O
Resistance to flow at 60L/min	4.4cm H ₂ O	3.4cm H ₂ O	3.6cm H ₂ O	3.2cm H ₂ O	4.0cm H ₂ O
Compressible volume	63ml*	63ml*	63ml	63ml*	63ml*
Weight	50g	59g	36g	50g	62g
Connectors	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F
Minimum tidal volume	>200ml	>200ml	>200ml	>200ml	>200ml
Accessories	Superset™ catheter mount with detachable elbow** 3504000 1997000	Smoothbore catheter mount** 5180000		Superset™ catheter mount** 3504000	Smoothbore catheter mount with double swivel elbow and double flip top cap** 5180001S

S Sterile

* Filter only

** Refer to information on catheter mount volumes in Patient Connections.

Flo-Guard

Low-resistance breathing filter for CPAP and bilevel

Flo-Guard provides a combination of filtration performance with low-resistance, ideal for use in both the hospital and the home, where high flow rates may be used, including CPAP, bilevel and cough-assist applications.

Electrostatic filter media

Providing an excellent filtration efficiency

Large surface area

To reduce resistance to flow

Clear housing

For good visibility

Conical shape

To aid airflow



Code	1690000
Box Qty.	50
Filtration efficiency	>99.99%
Resistance at 30L/min	0.4 cm H ₂ O
Resistance at 60L/min	0.8cm H ₂ O
Resistance at 90L/min	1.4cm H ₂ O
Resistance at 120L/min	2.0cm H ₂ O
Resistance at 150Lmin	2.7cm H ₂ O
Resistance at 180L/min	3.4cm H ₂ O
Resistance at 210L/min	4.3cm H ₂ O
Resistance at 240L/min	5.1cm H ₂ O
Compressible volume	80ml
Weight	27.8g
Connectors	22F-22M



Air-Guard™

For use in respiratory systems and oxygen concentrators


The Intersurgical Air-Guard is a hydrophobic pleated mechanical filter for the protection of oxygen concentrators and other respiratory equipment.

Providing an excellent level (>99.9999%) of protection against bacterial and viral challenges, the Air-Guard's product performance has been tested, validated and verified at independent microbiological laboratories.


The Air Guard pleated mechanical filter is validated for 24-hours use but this can be extended when used in non-humidified gas flow. The product can be used for extended periods of time for the protection of oxygen concentrators, upon the discretion of the clinician.

Please note: This product is contraindicated for use at the patient connection end of the breathing system.

Respiratory system protection




Oxygen concentrator protection



Air-Guard Clear

Code	1790000 (S*)
Box Qty.	50 (40*)
Filtration efficiency	>99.9999%
Resistance to flow at 30L/min	0.9cm H ₂ O
Resistance to flow at 60L/min	1.9cm H ₂ O
Compressible volume	120ml
Weight	56g
Connectors	22F-22M/15F
Minimum Tidal Volume	>360ml



 Sterile option available

(S*) Add an S to the seven digit code number for the sterile version of this product eg. 5000000S (sterile box quantity is shown in brackets).

Intersurgical Pulmo-Protect™ lung function filter


Forced expiratory and inspiratory manoeuvres are used to assess a patient's lung function. This assessment helps diagnose the extent of any disease process in the patient's lung, as a pre and postoperative assessment, during smoking cessation and occupational health programmes.

During these tests patients may generate peak flows as high as 12 L/sec (720 LPM) and expel infective droplets, which could contaminate the pulmonary function test equipment. As a result cross infection between patients is a real risk¹.

Pulmo-Protect combines low resistance with a high bacterial and viral efficiency filter designed specifically to protect the patient and equipment during lung function tests.

Pulmo-Protect provides:

- Protection of pulmonary function test equipment
- High filtration properties reduces risk of contamination between patients
- Low resistance performance ensures efficacy of results and complies with ATS/ERS recommendations²
- Low functional volume
- Microbiological filtration effectiveness has been independently tested and validated to provide >99.99% efficiency against bacteria and viruses³
- Flexible, comfortable, disposable mouth piece helps with patient compliance and improve the effectiveness of the test
- Comfortable single patient nose clip prevents patient-to-patient touch contamination
- Device compatible range of colour coded filters available individually or as a complete pulmonary function test kit



Pulmo-Protect options

For use with the following test devices	Pulmo-Protect lung function filter	Filter, mouthpiece and nose clip	Filter, flexible mouthpiece and nose clip
JAEGER®, MasterScreen, SensorMedics®, Vmax™, Micro Medical®, Chest and Microgard®	1691000 (Qty 50) ○	1691050 (Qty 30) ○	1691010 (Qty 30) ○
Medisoft, BodyBox, HypAir Compact + and SpiroAir	1691001 (Qty 50) ●	1691051 (Qty 30) ●	1691011 (Qty 30) ●
Fukuda Denshi®, SP-350 and Fudac-77	1691002 (Qty 50) ●	1691052 (Qty 30) ●	1691012 (Qty 30) ●

Accessories

Code	Description	Box Qty.
1691020	Flexible mouthpiece	10
1691040	Nose clip	35

References:

References: 1. Journal of Respiratory Medicine 2005;09:015 An audit into the efficacy of single use bacterial/viral filters for the prevention of equipment contamination during lung function assessment. 2. European Respiratory Journal 2005; 26: 319–338 Standardisation of spirometry M.R. Miller, J. Hankinson, V. Brusasco, F. Burgos, et al. 3. Nelson Labs 771942B.1 4. Nelson Labs 771943B.1

JAEGER is a registered trademark of CareFusion Germany 234 GmbH. SensorMedics is a registered trademark of SensorMedics Corporation. Vmax is a trademark of SensorMedics Corporation. Micro Medical is a registered trademark of CareFusion U.K. 232 Limited. Microgard is a registered trademark of Becton Dickinson and Company. Fukuda Denshi is a registered trademark of Fukuda Denshi Co., Ltd.


Heat and Moisture Exchangers (HMEs)

In normal respiration the anatomy of the upper airway helps to warm and humidify inspired air, and to retain the warmth and moisture contained in expired air. During inspiration, even cold or dry air is typically heated to 37 °C and, fully saturated, contains 44mg H₂O per litre. In mechanical ventilation or anaesthesia, the patient's upper airway may be bypassed by the introduction of a tracheal tube. As a result the patient's lungs may be confronted with cold dry inspired gas.

Prolonged exposure to dry ventilatory gases can lead to:


- Localised inflammation of the trachea
- A reduction in ciliary function
- Retention and thickening of secretions
- Lowering of patient temperature
- Reduction in cardiopulmonary function
- Increased risk of tracheostomy tube occlusion
- Extended duration and cost of care

Respiratory epithelium adversely affected by heat and moisture loss



Hydro-Trach™ T range

A heat and moisture exchanger designed for use on tracheostomised patients, the Hydro-Trach T is an ideal product for prolonged use with spontaneously breathing patients.



▶ Video available at www.intersurgical.com

Average Fi O₂ at variable O₂ flow rates

Code	1873000 (S*)	1874000 (S*)
Box Qty.	25 (100*)	40 (30*)
Moisture loss	13.2mg H ₂ O/L	13.2mg H ₂ O/L
Calculated moisture return	26mg H ₂ O/L	26mg H ₂ O/L
Resistance at 30L/min	0.2cm H ₂ O	0.2cm H ₂ O
Resistance at 60L/min	0.7cm H ₂ O	0.7cm H ₂ O
Compressible volume	19ml	19ml
Weight	8g	8g
Connectors	15F	15F
Minimum tidal volume	>60ml	>60ml
Accessories		1.8m oxygen tube


Oxygen (L/min)	Fi O ₂ at 15 BPM
1	26.5%
2	32.2%
3	35.0%
4	38.2%
5	41.8%
6	44.9%
7	47.9%
8	50.6%

Sterile option available

(S*) Add an S to the seven digit code number for the sterile version of this product

Hydro-Therm™ HME range

A dedicated range of Heat and Moisture Exchangers designed to replicate the functions of the body's upper airway by conserving expired heat and moisture and returning them to the patient during inhalation. The Hydro-Therm is a small volume, lightweight device which is clinically suitable over a wide range of patient sizes.



Code	1850000	1855000
Box Qty.	20	20
Luer lock port		✓
Moisture loss	7.8mg H ₂ O/L	7.8mg H ₂ O/L
Calculated moisture return	30.8mg H ₂ O/L	30.8mg H ₂ O/L
Resistance at 30L/min	0.5cm H ₂ O	0.5cm H ₂ O
Resistance at 60L/min	1.6cm H ₂ O	1.6cm H ₂ O
Compressible volume	20ml	20ml
Weight	11.8g	11.8g
Connectors	15M-22M/15F	15M-22M/15F
Minimum tidal volume	>60ml	>60ml

Hydro-Therm™ 3 HME range

A dedicated range of Heat and Moisture Exchangers designed to replicate the functions of the body's upper airway by conserving expired heat and moisture and returning them to the patient during inhalation. The Hydro-Therm 3 is a large-volume HME with rounded housing for use in anaesthesia and intensive care.



Code	1560000	1560011
Box Qty.	150	100
Luer lock port	✓	✓
Moisture loss	6.7mg H ₂ O/L	6.7mg H ₂ O/L
Calculated moisture return	31.6mg H ₂ O/L	31.6mg H ₂ O/L
Resistance at 30L/min	0.2cm H ₂ O	0.3cm H ₂ O
Resistance at 60L/min	0.8cm H ₂ O	1.3cm H ₂ O
Compressible volume	59ml	59ml*
Weight	31g	45g
Connectors	22F/15M-22M/15F	22F/15M-22F/15M
Minimum tidal volume	>200ml	>200ml
Accessories		Superset™ catheter mount** 3504000

* HME only

** Refer to information on catheter mount volumes in Patient Connections.

The Heat and Moisture Exchanging Filters (HMEFs)

The range of Heat and Moisture Exchanging Filters (HMEFs) combines the filtration efficiency of dedicated breathing filters with optimum moisture return provided by the addition of an HME element. Designed for use at the patient connection.

Filta-Therm™ range

High efficiency

A high-efficiency HMEF offered with a variety of patient connections. Filta-Therm Plus (1941001) provides the optimum solution for intensive care with improved HME performance and high filtration efficiency.



Code	1941001	1941000	1942000	1941011
Box Qty.	70	70	70	20
Luer lock port	✓	✓		✓
Moisture loss	8.3mg H ₂ O/L	9.3mg H ₂ O/L	9.3mg H ₂ O/L	9.3mg H ₂ O/L
Calculated moisture return	30.3mg H ₂ O/L	29.5mg H ₂ O/L	29.5mg H ₂ O/L	29.5mg H ₂ O/L
Filtration efficiency	>99.999%	>99.999%	>99.999%	>99.999%
Resistance at 30L/min	1.1cm H ₂ O	1.0cm H ₂ O	1.1cm H ₂ O	1.1cm H ₂ O
Resistance at 60L/min	2.6cm H ₂ O	2.4cm H ₂ O	2.5cm H ₂ O	2.8cm H ₂ O
Compressible volume	65ml	65ml	66ml	65ml*
Weight	44g	43g	42g	57g
Connectors	22F-22M/15F	22F-22M/15F	22F-22M/15F	22F-22M/15F
Minimum tidal volume	>200ml	>200ml	>200ml	>200ml
Accessories				Superset™ catheter mount** 3504000

Code	1941197	1941351	1906000	2015000
Box Qty.	50	20	70	20
Luer lock port	✓	✓	✓	✓
Moisture loss	9.3mg H ₂ O/L	9.3mg H ₂ O/L	9.9mg H ₂ O/L	9.9mg H ₂ O/L
Calculated moisture return	29.5mg H ₂ O/L	29.5mg H ₂ O/L	28.9mg H ₂ O/L	28.9mg H ₂ O/L
Filtration efficiency	>99.999%	>99.999%	>99.999%	>99.999%
Resistance at 30L/min	1.4cm H ₂ O	1.4cm H ₂ O	1.2cm H ₂ O	1.2cm H ₂ O (HMEF only)
Resistance at 60L/min	3.8cm H ₂ O	4.0cm H ₂ O	2.4cm H ₂ O	2.4cm H ₂ O (HMEF only)
Compressible volume	65ml*	65ml*	70ml	70ml (HMEF only)
Weight	49g	60g	45g	45g (HMEF only)
Connectors	22F-22M/15F	22F-22M/15F	22M-22M-22M/15F	22F-22F-22M/15F
Minimum tidal volume	>200ml	>200ml	>210ml	>210ml
Accessories	Fixed elbow** 1997000	Flexible catheter mount with fixed elbow** 3501000		Breathing system 1.6m

* HMEF only

** Refer to information on catheter mount volumes in Patient Connections.

Inter-Therm™ range

Sterile

The Inter-Therm range of sterile HMEFs is designed for use in breathing systems in the operating theatre and intensive care unit.

The Inter-Therm includes a wound paper HME media, providing excellent humidification and low resistance properties.

Inter-Therm Mini Angled offers an easy to use option with an integral 90° elbow for paediatrics, reducing the need for an additional catheter mount or separate patient elbow.

Turbo webs

Evenly distributes the air flow to ensure all of the filter media is utilised

Wound paper HME element

Providing high moisture return





HME on the patient side

Ensuring low resistance to flow during use

Self retaining port cap

Improving patient safety



Code	1341000S 	1341007S 	1341197S 	1341351S 
Box Qty.	50	50	50	50
Luer lock port	✓		✓	✓
Moisture loss	6mg H ₂ O/L	6mg H ₂ O/L	6mg H ₂ O/L	6mg H ₂ O/L
Calculated moisture return	32.3mg H ₂ O/L	32.3mg H ₂ O/L	32.3mg H ₂ O/L	32.3mg H ₂ O/L
Filtration efficiency	>99.998%	>99.998%	>99.998%	>99.998%
Resistance at 30L/min	1.6cm H ₂ O	1.6cm H ₂ O	1.5cm H ₂ O	1.7cm H ₂ O
Resistance at 60L/min	2.7cm H ₂ O	2.7cm H ₂ O	4.0cm H ₂ O	4.6cm H ₂ O
Compressible volume	57ml	57ml	57ml*	57ml*
Weight	31g	30g	38g	47g
Connectors	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F
Minimum tidal volume	>180ml	>180ml	>180ml	>180ml
Accessories			Elbow** 1997000	Flexible catheter mount with fixed elbow** 3501000

 Sterile

* HMEF only

** Refer to information on catheter mount volumes in Patient Connections.

Inter-Therm™ range



Code	1341211S 	1341974S 	1341011S 	1341012S 
Box Qty.	50	50	50	50
Luer lock port	✓	✓	✓	✓
Moisture loss	6mg H ₂ O/L	6mg H ₂ O/L	6mg H ₂ O/L	6mg H ₂ O/L
Calculated moisture return	32.3mg H ₂ O/L	32.3mg H ₂ O/L	32.3mg H ₂ O/L	32.3mg H ₂ O/L
Filtration efficiency	>99.998%	>99.998%	>99.998%	>99.998%
Resistance at 30L/min	1.4cm H ₂ O	1.6cm H ₂ O	1.4cm H ₂ O	1.6cm H ₂ O
Resistance at 60L/min	3.5cm H ₂ O	4.4cm H ₂ O	3.5cm H ₂ O	4.4cm H ₂ O
Compressible volume	57ml*	57ml*	57ml*	57ml*
Weight	49g	51g	45g	46g
Connectors	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F
Minimum tidal volume	>180ml	>180ml	>180ml	>180ml
Accessories	Superset™ catheter mount and safety cap** 3504000	Superset™ catheter mount with detachable fixed elbow** 3504000, 1997000	Superset™ catheter mount** 3504000	Superset™ catheter mount with fixed elbow** 3502000



Inter-Therm™ Mini






Inter-Therm™ Mini



Inter-Therm™ Mini Angled



Code	1341580S 	1331000S 	1331197S  <i>NEW</i>	1332000S 
Box Qty.	50	50	50	50
Luer lock port	✓	✓	✓	✓
Moisture loss	6mg H ₂ O/L	8.4mg H ₂ O/L	8.4mg H ₂ O/L	8.5mg H ₂ O/L
Calculated moisture return	32.3mg H ₂ O/L	30.2mg H ₂ O/L	30.2mg H ₂ O/L	30.1mg H ₂ O/L
Filtration efficiency	>99.998%	>99.99%	>99.99%	>99.99%
Resistance at 30L/min	1.6cm H ₂ O	2.2cm H ₂ O	2.6cm H ₂ O	2.2cm H ₂ O
Resistance at 60L/min	3.9cm H ₂ O	N/A	N/A	N/A
Compressible volume	57ml*	28ml	28ml*	29ml
Weight	58g	20g	26g	19g
Connectors	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F	15M-22M/15F
Minimum tidal volume	>180ml	>90ml	>90ml	>90ml
Accessories	Smoothbore catheter mount with double swivel elbow and double flip top cap** 5180000		Elbow** 1997000	

 Sterile

* HMEF only

** Refer to information on catheter mount volumes in Patient Connections.

Clear-Therm™ range

Clear-Therm 3 and Clear-Therm Angled

HMEF for use in anaesthesia and intensive care with the option of an integral 90° elbow, reducing the need for an additional catheter mount or separate patient elbow.



Code	1541000	1541011	1541012	1541197
Box Qty.	150	75	75	50
Luer lock port	✓	✓	✓	✓
Moisture loss	7.8mg H ₂ O/L	7.8mg H ₂ O/L	7.8mg H ₂ O/L	7.8mg H ₂ O/L
Calculated moisture return	30.8mg H ₂ O/L	30.8mg H ₂ O/L	30.8mg H ₂ O/L	30.8mg H ₂ O/L
Filtration efficiency	>99.99%	>99.99%	>99.99%	>99.99%
Resistance at 30L/min	0.8cm H ₂ O	1.0cm H ₂ O	1.2cm H ₂ O	1.2cm H ₂ O
Resistance at 60L/min	2.1cm H ₂ O	2.6cm H ₂ O	3.4cm H ₂ O	3.3cm H ₂ O
Compressible volume	60ml	60ml*	60ml*	60ml*
Weight	29g	44g	45g	36g
Connectors	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F
Minimum tidal volume	>200ml	>200ml	>200ml	>200ml
Accessories		Superset™ catheter mount** 3504000	Superset™ catheter mount with fixed elbow** 3502000	Fixed elbow** 1997000



Code	1541351	1541974	1542000
Box Qty.	70	40	75
Luer lock port	✓	✓	✓
Moisture loss	7.8mg H ₂ O/L	7.8mg H ₂ O/L	7.4mg H ₂ O/L
Calculated moisture return	30.8mg H ₂ O/L	30.8mg H ₂ O/L	31.1mg H ₂ O/L
Filtration efficiency	>99.99%	>99.99%	>99.99%
Resistance at 30L/min	1.3cm H ₂ O	1.2cm H ₂ O	0.9cm H ₂ O
Resistance at 60L/min	3.7cm H ₂ O	3.5cm H ₂ O	2.3cm H ₂ O
Compressible volume	60ml*	60ml*	74ml
Weight	47g	50g	38g
Connectors	22F/15M-22M/15F	22F/15M-22M/15F	22F/15M-22M/15F
Minimum tidal volume	>200ml	>200ml	>225ml
Accessories	Flexible catheter mount with fixed elbow** 3501000	Superset™ catheter mount with detachable fixed elbow** 3504000 + 1997000	

* HMEF only

** Refer to information on catheter mount volumes in Patient Connections.

Clear-Therm™ range

Clear-Therm Midi

A low volume option for minimising dead space in anaesthesia.



Code	1641000	1641197	1641012
Box Qty.	100	50	75
Luer lock port	✓	✓	✓
Moisture loss	13.3mg H ₂ O/L	13.3mg H ₂ O/L	13.3mg H ₂ O/L
Calculated moisture return	25.9mg H ₂ O/L	25.9mg H ₂ O/L	25.9mg H ₂ O/L
Filtration efficiency	>99.9%	>99.9%	>99.9%
Resistance at 30L/min	0.9cm H ₂ O	1.2cm H ₂ O	1.2cm H ₂ O
Resistance at 60L/min	2.1cm H ₂ O	3.1cm H ₂ O	3.2cm H ₂ O
Compressible volume	34ml	34ml*	34ml*
Weight	19g	26g	35g
Connectors	22F-22M/15F	22F-22M/15F	22F-22M/15F
Minimum tidal volume	>100ml	>100ml	>100ml
Accessories		Fixed elbow** 1997000	Superset™ catheter mount with fixed elbow** 3504000 + 1997000

Clear-Therm Mini and Micro

For paediatrics and neonates.



Clear-Therm Mini

Clear-Therm Mini

Clear-Therm Micro

Clear-Therm Micro

Code	1831000	1831197	1441000	1441197
Box Qty.	40	40	20	40
Luer lock port	✓	✓	✓	✓
Moisture loss	6.8mg H ₂ O/L	6.8mg H ₂ O/L	12.3mg H ₂ O/L	12.3mg H ₂ O/L
Calculated moisture return	31.7mg H ₂ O/L	31.7mg H ₂ O/L	26.8mg H ₂ O/L	26.8mg H ₂ O/L
Filtration efficiency	>99.99%	>99.99%	>99.99%	>99.99%
Resistance at 30L/min	1.5cm H ₂ O	2.0cm H ₂ O	N/A	N/A
Resistance at 10L/min	N/A	N/A	1.0cm H ₂ O	0.9cm H ₂ O
Compressible volume	26ml	26ml*	11ml	11ml*
Weight	22g	29g	11g	19g
Connectors	22F/15M-22M/15F	22F/15M-22M/15F	15M-15F	15M -22m/15F
Minimum tidal volume	>90ml	>90ml	>35ml	>35ml
Accessories		Fixed elbow** 1997000		Fixed elbow** 1997000

* HMEF only

** Refer to information on catheter mount volumes in Patient Connections.