

# THE ULTIMATE SOLUTION FOR WOUND CLOSURE




**exofin<sup>®</sup>**

**TOPICAL SKIN ADHESIVE**  
High Viscosity, 2-Octyl Cyanoacrylate


## Comparative Strength Data

### Tensile Strength




ASTM F2258-05\*  
Measurement = PSI

### Lap Shear



ASTM F2255-05\*  
Measurement = PSI

### Wound Closure



ASTM F2458-05\*  
Measurement = Pounds Force

## Average Cure Speed



\*Data on file

Tack Free Time\*



# ABOUT



**[exo·fin]** **exo-** "Outer, outside, exterior, exoskeleton"  
**-fin** "Finished, completed, finality"

*The results of two years' observation and feedback from healthcare professionals worldwide, clearly showed that no single topical skin adhesive contained all the attributes that mattered most to them. From this valuable information, the renowned research and development facility at Chemence Medical created a topical skin adhesive that delivers on all the requirements asked for. We are excited to introduce the next generation of medical adhesives: exofin® Topical Skin Adhesive.*

**JEFF ROBERSON**

President, Chemence Medical, Inc.

## The exofin® Advantage



### 14-Day Microbial Barrier

- ✓ Exofin® topical skin adhesive provides a 14 day microbial barrier\*
- ✓ The longest in the industry

## PRODUCT FEATURES



### FAST

Average set time is less than 30 seconds, which is 2 x faster than competitive products\*



### STRONG

2-Octyl formulation is proven to be the strongest and most flexible topical skin adhesive available\*\*

The high viscosity formula and proprietary production process gives exofin® a bonding strength between 19% to 25% greater than competitive 2-Octyl adhesives\*



### FLEXIBLE

High purity 2-Octyl technology ensures flexibility on skin and prevents cracking.

Added flexibility allows the adhesive to remain intact for up to 5 -10 days.



### SAFE

Provides a microbial barrier against clinically relevant pathogens, including MRSA †

100% glass free packaging ensures zero risk of sharps related injury to medical professionals or patients.



### INNOVATIVE

Soft and flexible applicator glides across uneven surfaces, reaching difficult to access areas and ensuring that the adhesive creates a perfect seal over the wound or incision.

Transparent Tip allows users to see the adhesive prior to expression, therefore eliminating the "surprise" of too much glue being dispensed all at once.

# GLASS-FREE PACKAGING WITH INNOVATIVE APPLICATOR



## ANGLED TIP

Angled tip moves easily along body contours

Clear applicator allows user to see adhesive before it is dispensed



## LONG CONNECTOR PIECE

Raised grooves allow for easy grip with gloved fingers

## ALUMINUM TUBE

Glass-free packaging for safety and convenience

Simple dispensing

Available in 1ml and 0.5ml sizes



## EASY TO USE

Application is simple. Tighten the applicator tip and begin to express adhesive.

Controlled application tip ensures the same reliable result each time for both multiple or single coats.



## VISCOUS

Proprietary formulation and production process gives a high viscosity.

Reduced risk of adhesive migration (flowing) away from the incision site!



## PATIENT FRIENDLY

Waterproof barrier allows patient to take a shower almost immediately after the medical procedure.

Less trauma to the wound site.

No need for hospital revisit for suture removal nor anaesthesia in emergency or outpatient applications.



## COST EFFECTIVE


Saves time and reduces costs of wound closure procedures.

1 ml tube – 43% more adhesive than Dermabond® Advanced and 25% more than the largest comparable adhesive unit available.

0.5 ml tube – 39% more adhesive than Dermabond® Mini.



## Applications



exofin® Topical Skin Adhesive has an increased viscosity to reduce the risk of unintended placement of the adhesive during application due to migration of the liquid adhesive from the wound site.

In vitro studies have shown that exofin® acts as a barrier to microbial penetration as long as the adhesive film remains intact.

exofin® is intended for topical skin applications only - and to hold closed easily approximated skin edges of wounds - from both surgical incisions, including incisions from minimally invasive surgery, and simple, thoroughly cleansed, trauma-induced lacerations.

exofin® Topical Skin Adhesive may be used in conjunction with, but not in place of, deep dermal sutures.



**GENERAL  
OPERATIONS**



**OBSTETRICS &  
GYNAECOLOGY  
PROCEDURES**



**ORTHOPAEDIC  
SURGERIES**




**CARDIOVASCULAR  
SURGERIES**



**C-SECTIONS &  
COSMETIC  
SURGERY**

## Exofin® in action, examples:



**LAPAROSCOPY**



**C-SECTION**



**ABDOMINOPLASTY**



**KNEE  
ARTHROSCOPY**



**MASTECTOMY**



exofin®  
TOPICAL SKIN ADHESIVE

# DIRECTIONS

## Directions for Use



- 1** The incision or trauma site must be clean and dry before applying adhesive. Ensure that the wound edges are easily approximated.
- 2** Screw the Applicator Tip onto the threaded tube in a clockwise direction until the applicator is tight and seated on the tube. The applicator will puncture the foil membrane on the tube once seated and adhesive will be allowed to flow.



- 3** When holding the tube horizontal and away from the patient, gently squeeze the tube from the bottom, moving the adhesive upward until it becomes visible through the applicator device.



- 4** While approximating skin edges, apply by squeezing the tube gently and continuously as you spread the adhesive back and forth over the wound site.  
  
Continue to hold the wound edges in approximation until the adhesive becomes tacky. Typically less than 20 seconds. A smooth and even coat of adhesive is desirable.



- 5** Once adhesive has completely polymerised (non-tacky) you may cover the site with a secondary bandage. If a secondary bandage is used, DO NOT apply to the adhesive area until it is completely tack free. Applying while tacky may result in adhesive removal when removing secondary bandage and possibly dehiscence of the wound site.
- 6** Discard adhesive device according to normal protocol after use.



**exofin® Micro Topical Skin Adhesive 0.5 ml**  
**EX91510 (single box of 10 tubes)**  
**EX91510-01 (6 boxes, 10 tubes each)**



**exofin® Topical Skin Adhesive 1 ml**  
**EX91011 (single box of 10 tubes)**  
**EX91011-01 (6 boxes, 10 tubes each)**



[chemencemedical-intl.com](http://chemencemedical-intl.com)

Chemence Ltd.  
13 Princeswood Road, Corby  
Northamptonshire, NN17 4XD, UK  
**+ 44 (0)1536 402600**

Chemence Medical, Inc.  
200 Technology Drive  
Alpharetta, GA, 30005, USA  
**+1 844-633-4583**

Chemence Sp. z o.o.  
Aleja Jana Pawła II 22  
00-133 Warsaw, Poland  
**+1 770 658 4302**

All Contents © Copyright 2024 Chemence Medical, Inc. All rights reserved.  
No modification or use, or further reproduction or copying, of any of the content is authorized without the express consent of Chemence Medical, Inc.