# Medtronic

## Septal Button

Bouton septal / Bottone per setto nasale /  
Septumknopf / Botón septal / Septale knop /  
Septumknop / Väliseinänappi / Septumknapp /  
Botão septal / Διαφραγματικό Πλώμα /  
Wkładka przegrodowa / Septální knoflík /  
Septalis záródugó / Septumknapp /  
Septal Dügme

---

## Product Information and Instructions

Informations et instructions concernant le produit /  
Informazioni sul prodotto ed istruzioni per l'uso /  
Produktinformation und Gebrauchsanweisung /  
Instrucciones e información sobre el producto /  
Productinformatie en instructies / Produktoplysninger  
og vejledning / tuotetta koskevat tiedot ja käyttöohjeet /  
Produktinformation och instruktioner / Informação  
do produto e instruções / Πληροφορίες προϊόντος  
και οδηγίες / Informacje o produkcie i instrukcje /  
Údaje o výrobku a pokyny / Termékinformáció  
és használati utasítás / Produktinformasjon  
og instruksjoner / Ürün Bilgileri ve Talimatlar

**DEVICE DESCRIPTION AND INTENDED USE**

The Medtronic Septal Button is a soft silicone button which may be used for septal perforation closure, particularly for those patients who are poor candidates for surgical correction.

**INDICATIONS FOR USE**

Closure with a septal button is indicated when the nature of the perforation (size and location) and septal anatomy are appropriate, and /or when surgical closure is contraindicated. Generally, anterior septal perforations that are 2 cm or less in maximum diameter are amenable for closure with a button.

**CONTRAINDICATIONS**

Septal deformities are relative contraindications for button placement, as the flange on the convex side of the deformity may increase nasal obstruction.

**INSTRUCTIONS FOR USE**

Proper surgical procedures and techniques are the responsibility of the medical professional. Each physician must evaluate the appropriateness of this procedure based on his/her own medical training and experience and the needs of the patient. The following procedure is provided by Manning M. Goldsmith, M.D. for informational purposes only.

**INSERTION OF THE ONE-PIECE SEPTAL BUTTON**


After adequate cocaineization of the nose, 1% lidocaine with 1:100,000 epinephrine is infiltrated about the margins of the septal perforation. The flanges of the button are trimmed to an approximate size that will permit 3 - 4 mm of overlap with respect to the margins of the perforation. The edges of one flange are then sewn together in a rosette fashion and the suture is tied and left long. The sutured flange is directed into the nasal cavity and the suture is passed through the perforation and retrieved from the contralateral side. Gentle traction is applied to the suture and the rosette is pulled through the perforation while a hemostat simultaneously pushes it through from the contralateral side. The suture is then cut and the flange resumes its original configuration. Minor adjustments of the size of the flanges may be made. Care must be taken so that the flanges of the button do not abut the junction of the septum with the upper lateral cartilage, as this is a particularly annoying sensation for the patient.

### INSERTION OF THE TWO-PIECE SEPTAL BUTTON

After cocainization of the nose, the flanges of the button are trimmed to an approximate size that will permit 3 - 4 mm of overlap with respect to the margins of the perforation and grasped with a hemostat as shown in Figure 1. The male connecting piece of the button is then introduced through the perforation and retrieved from the contralateral side with a hemostat (Figure 2). A nasal speculum or hemostat is used to stabilize the female flange of the button as the male connecting piece is pulled through the female connector (Figure 3). The button is adjusted to the desired width by pulling the connecting flanges of the male piece sequentially through the female connector. The excess male connecting piece is severed, leaving the most proximal connecting flange intact to secure the two pieces of the button (Figure 4). Care must be taken so that the flanges of the button do not abut the junction of the septum with the upper lateral cartilage, as this is a particularly annoying sensation for the patient. Should removal be necessary, the connecting flange of the male piece is cut with a scalpel and each piece of the button is delivered through separate nostrils.

### STERILITY

The Septal Button is provided STERILE and is intended for SINGLE PATIENT USE ONLY. The package should be carefully inspected for punctures, tears, security of the seal or any other evidence of the sealed pouch having been compromised prior to placement of the contents in the sterile field.

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**Medtronic Xomed, Inc.
Jacksonville, FL 32216-0980 U.S.A.**