Eurotrol

Your Global Reference Point for Quality Control

HemoTrol® Duo

Low

REF	AN01637A01
-----	------------

Normal

REF	AN01637A02
-----	------------

High

REF	AN01637A03
-----	------------

EN
Intended use

HemoTrol® Duo is an assayed quality control material for professional use to verify the performance characteristics of the HemoCue® Hb 301 and the HemoCue® Hb 801 System. HemoTrol® Duo is intended for the quantitative determination of hemoglobin.

Summary and Principle
 HemoTrol® Duo is a quality control material containing a known hemoglobin concentration in three clinically relevant ranges. When measured like a patient sample, the regular use of the quality control material can provide proof for proper functioning of the HemoCue® Hb 301 and the HemoCue® Hb 801 System. Quality control materials should be used in accordance with local, state, and/or federal regulations or accreditation requirements.

Warnings and precautions
 1. For in vitro diagnostic use.
 2. Only use when the quality control has a reddish color.
 3. Do not freeze.
 4. Only mix by hand.
 5. Caution: HemoTrol® Duo contains material of bovine origin. Bovine based materials do not carry Hepatitis B Surface Antigen (HbsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus (HIV-1/HIV-2). The product is free from transmissible spongiform encephalopathies (TSEs). However, no test method can offer complete assurance that products derived from bovine source material are free from infectious agents and may transmit infectious agents.

6. Follow good laboratory practice while using the product.
 7. Disposition of all waste material should be in accordance with the local guidelines.

Storage and Stability
 For the shelf life at 2 to 8 °C (36 – 46 °F). see the expiration date on the vial label. After opening, the product is stable at 2 to 30 °C (36 - 86 °F) for 31 days when handled according to procedure.

Procedure
 For optimal performance the product should be treated as a patient sample and instructions for use of the HemoCue® Hb 301 or HemoCue® Hb 801 System should be followed.
 1. Equilibrate the vial for 15 minutes to a temperature of 15 to 30 °C (59 – 86 °F).
 2. Gently mix the sample by inverting the vial 8 – 10 times before sampling.
 3. Hold the vial in an upright position and tap the cap to remove liquid from the dropper tip (Figure 1).
 4. Open the cap carefully without squeezing the vial (Figure 2).
 5. Discard the first drop of the quality control liquid.
 6. Apply a drop on a hydrophobic surface and use it immediately for measurement as indicated in the analyzer user manual.
 7. After sampling, remove any liquid from the cap with a clean tissue and recap properly.
 8. Return to the indicated storage conditions after use. The entire procedure has to be repeated upon reuse of the quality control during the open-vial time period.

Limitations
 1. The quality control material should not be used past open-vial time period or expiration date.
 2. The quality control material is not intended for use as calibrator / standard.
 3. Incomplete mixing of the quality control material prior to use invalidates the measured sample and any following measurements made with the remaining quality control material.
 4. The ranges of the product are intended only as guidelines and laboratories should determine the ranges based on their own test system and tolerance limits.

Expected values
 For expected values see the batch code specific table. The means and ranges were derived from replicate analyses on the HemoCue® Hb 301 and HemoCue® Hb 801 System and are specific for that batch code of the product.

EN
Intended use

HemoTrol® Duo is an assayed quality control material for professional use to verify the performance characteristics of the HemoCue® Hb 301 and the HemoCue® Hb 801 System. HemoTrol® Duo is intended for the quantitative determination of hemoglobin.

Summary and Principle

HemoTrol® Duo is a quality control material containing a known hemoglobin concentration in three clinically relevant ranges. When measured like a patient sample, the regular use of the quality control material can provide proof for proper functioning of the HemoCue® Hb 301 and the HemoCue® Hb 801 System. Quality control materials should be used in accordance with local, state, and/or federal regulations or accreditation requirements.

Contents

Purified stroma-free bovine hemolysate and stabilizers.

Warnings and precautions

- For in vitro diagnostic use.
- Only use when the quality control has a reddish color.
- Do not freeze.
- Only mix by hand.

5. Caution: HemoTrol® Duo contains material of bovine origin. Bovine based materials do not carry Hepatitis B Surface Antigen (HbsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus (HIV-1/HIV-2). The product is free from transmissible spongiform encephalopathies (TSEs). However, no test method can offer complete assurance that products derived from bovine source material are free from infectious agents and may transmit infectious agents.

- Follow good laboratory practice while using the product.
- Disposition of all waste material should be in accordance with the local guidelines.

Storage and Stability

For the shelf life at 2 to 8 °C (36 – 46 °F). see the expiration date on the vial label. After opening, the product is stable at 2 to 30 °C (36 - 86 °F) for 31 days when handled according to procedure.

Procedure

For optimal performance the product should be treated as a patient sample and instructions for use of the HemoCue® Hb 301 or HemoCue® Hb 801 System should be followed.

- Equilibrate the vial for 15 minutes to a temperature of 15 to 30 °C (59 – 86 °F).
- Gently mix the sample by inverting the vial 8 – 10 times before sampling.
- Hold the vial in an upright position and tap the cap to remove liquid from the dropper tip (Figure 1).
- Open the cap carefully without squeezing the vial (Figure 2).
- Discard the first drop of the quality control liquid.
- Apply a drop on a hydrophobic surface and use it immediately for measurement as indicated in the analyzer user manual.
- After sampling, remove any liquid from the cap with a clean tissue and recap properly.
- Return to the indicated storage conditions after use. The entire procedure has to be repeated upon reuse of the quality control during the open-vial time period.

Limitations

- The quality control material should not be used past open-vial time period or expiration date.
- The quality control material is not intended for use as calibrator / standard.
- Incomplete mixing of the quality control material prior to use invalidates the measured sample and any following measurements made with the remaining quality control material.
- The ranges of the product are intended only as guidelines and laboratories should determine the ranges based on their own test system and tolerance limits.

Expected values

For expected values see the batch code specific table. The means and ranges were derived from replicate analyses on the HemoCue® Hb 301 and HemoCue® Hb 801 System and are specific for that batch code of the product.

For further information or technical assistance, please contact:

<div></div> <div>Manufacturer</div>						
<div>Eurotrol B.V.</div> <div>Keplerlaan 20, 6716 BS Ede, the Netherlands</div> <div>Phone: +31 318 695 777</div> <div>E-mail: office@eurotrol.com</div> <div>www.eurotrol.com</div>						
<div>Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the member state where the user and/or the patient is established.</div>						
<div>Revision history</div> <table border="1"> <tbody><tr> <th>Revision</th> <th>Effective Date</th> <th>Change to previous revision</th> </tr> <tr> <td>01</td> <td>2019-11-22</td> <td>First revision</td> </tr> </tbody></table>	Revision	Effective Date	Change to previous revision	01	2019-11-22	First revision
Revision	Effective Date	Change to previous revision				
01	2019-11-22	First revision				
<div>All trademarks are the property of their respective owners. © 2019, Eurotrol B.V.</div>						

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