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<b>Important</b>	EN
Before <b>first use</b> , every <b>subsequent</b> use and before returning for <b>repair</b> , the instruments must be cleaned and sterilised according to our processing instructions.	

**General Short Notice**

The instruments are only to be used according to their intended use in the medical fields, and only by the accordingly trained and qualified personnel. For specific applications and operational use, the consultant doctor, purchaser or user are responsible for having the appropriate training and information and the sufficient experience for handling the instruments. **nopa** instruments Medizintechnik GmbH, as manufacturer and seller of the devices, shall not be liable for any incidental or consequential damages caused by the improper use, handling or the improper treatment, sterilisation and maintenance. If the instruments are repaired by companies or persons that have not been authorised by **nopa** instruments for repair, the guarantee is null and void. Before each use the instruments must be examined for breaks, cracks, deformations, damage and functional efficiency. Parts to be especially carefully checked are blades, tips, locks, latches and all moving parts. Worn, corroded, deformed, porous or otherwise damaged instruments must be rejected.

**Storage**

Instruments should be kept in a clean, dry environment and individually in their packaging or in a protective container with individual compartments. Protect the areas such as tips, blades, etc. with the corresponding tubes, protective caps, gauze or cloth. Be especially careful that there are no chemicals in the immediate vicinity of the storage area.

**Materials used**

- >>> Stainless steel
- >>> Pure titanium
- >>> Titanium base alloy
- >>> Light metals aluminium

**Steel instruments**

The stainless steels used for instrument manufacturing (non-corrosive, stainless) form protective layers due to their specific passive layers. These steels are only conditionally resistant to attack by chloride ions and aggressive media and liquids!!

**Titanium instruments**

The instruments used that are made from titanium or from titanium base alloy are to be handled in the same way as the stainless steel instruments. There are no special instructions to observe.

**Aluminium instruments**

Only non-alkaline, neutral detergent and demineralised water is to be used (pH 7), as otherwise damage may be caused to anodised surfaces. Just a few cycles of alkaline cleaning processes cause colour fading and stains, in particular to colour anodised instruments. The following described, validated cleaning procedure has not been validated for medical devices made of aluminium. In addition to the efforts made by the manufacturer when selecting the right materials and during

their careful processing, the user must perform the professional ongoing maintenance and proper treatment of the surgical instruments.

**We recommend the following approved methods for reprocessing our reusable surgical instruments.**

This applies to all used materials as well as possible accessories being suitable for reprocessing, with the exception of instruments made of aluminium which have to be cleaned with non alkaline detergents.

**Preparation for disinfection and cleaning**

The instruments should be disinfected and cleaned immediately after use. Dismountable instruments must be disassembled for cleaning. Impurities should not dry on the objects, as this renders the disinfection and purification processes difficult.Cleaning: (recommendation for machine cleaning)

**Manual pre-cleaning:**

Clean the instruments immediately or no later than two hours after application and remove any visible impurities with a soft brush. Cavities must also be flushed (flushing gun, etc.). For this, use a suitable combination of cleaning and disinfecting solution. Then briefly rinse the instruments under clear water.

**Mechanical cleaning:**

Instruments with lumen must be flushed during machine cleaning.

**1. Purging**

Cold water without any additives to remove coarse dirt and foaming substances.

**2. Cleaning**

Perform cleaning at 55°C ± 2°C for at least 5 minutes. For the mechanical cleaning of thermo-stable and thermolabile instruments we recommend the alkaline cleaner Neodisher® MediClean forte; 0.5% strength in the machine. If there are increased concentrations of chloride in the water pitting corrosion and stress corrosion cracking may occur to the instruments. Such corrosion can be minimised by using alkaline cleaner or demineralised water.

**3. Neutralisation**

By adding a neutralising agent to the acid-base, the rinsing of alkaline detergent residues is facilitated. Even when using neutral cleaners, it is recommended, in the case of unfavourable water quality e.g. high salinity, to use a neutraliser to prevent deposit formation. It is advised that neutralisation with Neodisher® Z, is performed in warm water (about 40°C) 0.1% concentration.

**4. Rinsing**

Deionised water without additives

**5. Thermal disinfection/final rinsing**

Thermal disinfection is performed at 92°C ± 2°C for at least 5 min (A0 value> 3000).

**6. Dehydration**

Sufficient drying must be ensured either by the cleaning and disinfection device or by other appropriate measures. Drying at 55-60°C is performed for about 30 minutes. Should residual moisture still be present, post-drying is carried out in a drying oven at 60°C. However, drying time is not dependent on the load or the flushings.

**7. Sterilisation**

STERILISER: Steam autoclave with fractioned pre-vacuum: Temperature: 134°C, with a holding time of at least > 5 up to a maximum of 20 minutes and subsequent drying. At least three pre-vacuum cycles.

**Instrument sterilisation**

Sterilisation accessories, such as packaging material/storage resources, must be matched to both the instruments and the sterilisation method applied.

To avoid damage (joint stress cracks, deterioration of resilience) to instruments with notch lock, due to stress occurring in heating and cooling during the steriilisation process, these instruments must only be locked onto the first notch.

**ATTENTION:**

In surgical instruments that are chrome-plated but whose chrome protection has been damaged, corrosion occurs from the locations of these black points outwards. For this reason, such instru-ments should not be sterilised together with stainless instruments, as rust is transferred to the stainless instruments. This is not observable on the chromed instruments as rust cannot build up on chrome. Hence, given that the stainless instruments are not covered by this layer, we do not recommend that surgical instruments without a chrome-plated surface be sterilised in the same steriliser as the rust-free instruments.

**Steam sterilisation**

When using steam sterilisers for sterilising surgical instruments, it must be ensured that sterilisation steam without contamination is used. Steam for sterilisation purposes must comply with DIN 58946 Part 7 and DIN EN ISO 17665-1. If these requirements for steam quality are not met, impurities such as oil, chemicals, metal chips or rust pollution and/or consequential damage can be caused to the instruments. Damage caused by extraneous rust mainly have their roots in steam supply systems, which are not designed for stainless steel. In order to avoid detriment the manufacturer's instruction manual for steam sterilisers must be observed. For the avoidance of detriment to the instruments DIN 58946, Part 6 „Operation of large sterilisers“, paragraph 5, is to be observed. Sterilisation must comply respectively with DIN 58952 and DIN 58953.

Following sterilisation, the instruments must be stored dry until use. The drying of the instruments and the sterile inner envelope is achieved after cooling to room temperature (23 + 2°C). Excessive condensate quantities during sterilisation is avoided by noting the maximum recommen- ded weight of the loaded sterilisation wire sieves. Drying is facilitated by wrapping the sieves with a towel on the inside of the container or in the outer paper packaging.

**Process approved according to Report No. 1504.1706 and Report No. 1505.0756 with the listed procedures as well as disinfectants and detergents.**  
**During cleaning, disinfection and sterilisation the following sources in particular have to be observed:**

- >>> DIN EN ISO 17664 manufacturer information for the reproces sing of re-sterilisable instruments.
- >>> EN 285 Sterilisation - steam sterilisers - large-sterilisers.
- >>> DIN EN ISO 17665 Sterilisation of medical devices - validation and routine control for sterilisation with moist heat.
- >>> DIN EN 556-1 Sterilisation of medical devices - requirements for medical devices being sterilized inside the final packaging which are marked as isterileí ñ Part 1: Requirements for medical devices that have been sterilized inside the final packaging.
- >>> DIN 58946-7 Sterilisation; steam sterilisation; structural require- ments for large sterilisers, properly performed instrument preparation.
- >>> AK instrument preparation: http://www.a-k-i.org
- >>> Recommendations for the validation and routine control of sterilisation processes with moist heat for medical devices, recommendation of the DGKH [Deutsche Gesellschaft für Krankenhaushygiene] (German Society for Hospital Hygiene). http://www.dgkh.de
- >>> Hygiene requirements in the preparation of medical devices, recommendation of the Commission for hospital hygiene and infection prevention at the Robert Koch Institute (RKI) and the „Requirements for hygiene in the processing of medical devices“ of the Federal Institute for Drugs and Medical Devices (BfArM). http://www.rki.de

**Limiting the preparation**

Frequent re-preparation has little effect on these surgical instruments. The end of product life cycle is usually determined by wear and damage from use. Perform a professional disposal or a recycling system of the surgical instruments after end of product life.

**Care**

„Care“ includes the application of instrument oil (physiologically acceptable paraffin oil according to DAB 8, Ph. Eur. or USP) or instrument milk (emulsion of hydrocarbon, white oil, in water) onto the surfaces (particularly on the moving parts/joints) of the surgical instruments. Basically, surgical instruments must undergo adequate care, and namely before functional testing. The care agent must ensure that, even with the accumulated effect of their constant use, the „bonding“ of the joint parts is ruled out.

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