Seronorm™

CRP Liquid L-2

REF 202305

LOT 1311495

2016-12

12 x 1 mL

IVD

CE

Collaborating laboratories

The analytical data of Seronorm™ CRP Liquid have been elaborated in collaboration with the following independent laboratories:

- Department of Clinical Chemistry, Aker og Bærum Hospital, Bærum, Norway
- First Medical Laboratory, Oslo, Norway
- The Laboratory of SERO AS, Billingstadi, Norway
- Institute of Clinical Biochemistry, Oslo University Hospital, Oslo, Norway
- Central Laboratory of Vestfold Hospital, Tønsberg, Norway
- In-vicon, München, Germany
- Department of Medical Biochemistry, St. Olav Hospital HF, Trondheim, Norway
- Department of Clinical Chemistry, Selander Hospital Kristiansand, Norway
- Primary Health Care, Uppsala, Sweden

EN

Storage and stability

This product will be stable until the expiry date when stored unopened at 2-8 °C.

Once opened, CRP is stable for 20 days when stored tightly capped in the original vial at 2-8 °C, and for 7 days at 15-25 °C.

Limitations

All stability data require that bacterial contamination is avoided. Increased turbidity may indicate bacterial growth.

Assignment of values

In line with the requirements of the European IVD Directive 98/79/EC, the reference values are based on reference methods or transfered values in accordance with certain reference materials. At the time of establishment of this package insert, this reference method or material constituted the top of the calibration hierarchy for the given components. These values are presented with shaded background.

The mean analytical values and acceptable ranges given in the table below are derived from replicate analysis obtained through collaboration with an unbiassed number of independent laboratories and they are specific for this lot. The uncertainty of the analytical value is presented as a single number, U (with a coverage factor $k = 2$), which takes into account various factors including imprecision of the method applied for the assignment as well as the stability of the given component. Individual laboratory means should fall within the acceptable range. However, variations in the analytical value and acceptable range for each component may vary over time due to differences in reagent and/or calibration lots as well as the instrument and/or method used. It is therefore recommended that each laboratory establishes its own analytical mean and acceptable range and only use the provided values as a guide. If applicable, the RiBÅK range (Richtlinien der Bundesärztekammer zur Qualitätssicherung in medizinischen Laboratorien (2008)) is provided.

The value assignment has been established in accordance with the Essential Requirements (Annex 1) of the IVD Directive 98/79/EC, and the ISO 15177 International standard.

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