ms_06879942001V5.0

Roche CARDIAC D-Dimer

cobas[®]

REF 04877802 190

▽ 10

SYSTEM cobas h 232

Svenska

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Referenser

- 1 Baker WF. Diagnosis of deep venous thrombosis and pulmonary embolism. *Med Clin North Am* 1998;82:459-476.
- 2 Becker DM, Philbrick JT, Bachhuber TL, et al. D-dimer testing and acute venous thromboembolism. A shortcut to accurate diagnosis? *Arch Intern Med* 1996;156:939-946.
- 3 Bernardi E, Prandoni P, Lensing AW, et al. on the behalf of the Multicentre Italian D-dimer Ultrasound Study Investigator Group. D-dimer testing as an adjunct to ultrasonography in patients with clinically suspected deep vein thrombosis: prospective cohort study. *BMJ* 1998;317:1037-1040.
- 4 Bounameaux H, de Moerloose P, Perrier A, et al. Plasma measurement of D-dimer as diagnostic aid in suspected venous thromboembolism: an overview. *Thromb Haemost* 1994;71:1-6.
- 5 Brill-Edwards P, Lee A. D-dimer testing in the diagnosis of acute venous thromboembolism. *Thromb Haemost* 1999;82:688-694.
- 6 Demarmels Biasutti F, Lämmle B. Beitrag des Hämostaselabors bei der Diagnostik der tiefen Venenthrombose. *Therap Umschau* 1996;53:265-271.
- 7 Dempfle CE, Hafner G, Lestin HG, et al. Multizentrische Evaluierung von Tina-quant [a] D-Dimer. *J Lab Med* 1996;20:31-37.
- 8 Dempfle CE, Schraml M, Besenthal I, et al. Multicentre evaluation of a new point-of-care test for the quantitative determination of D-dimer. *Chim Clin Acta*.
- 9 van der Graaf F, van den Borne H, van der Kolk M, et al. Exclusion of deep vein thrombosis with d-dimer testing – Comparison of 13 d-dimer methods in 99 outpatients suspected of deep venous thrombosis using venography as reference standard. *Thromb Haemost* 2000;83:191-198.
- 10 Knecht MF. Stellenwert der D-Dimere in der Diagnostik thromboembolischer Erkrankungen. *medwelt* 1997;48:52-65.
- 11 Lestin HG, Ehrenteich U, Hergert M, et al. Untersuchungen zur Relevanz hämostaseeologischer Marker für die Diagnostik, Therapie und Verlaufskontrolle tiefer Venenthrombosen. *Clin Lab* 1996;42:745-756.
- 12 de Moerloose P, Michiels JJ, Bounameaux H. The place of D-Dimer testing in an integrated approach of patients suspected of pulmonary embolism. *Semin Thromb Hemostasis* 1998;24:409-412.
- 13 Perrier A, Buswell L, Bounameaux H, et al. Cost-effectiveness of noninvasive diagnostic aids in suspected pulmonary embolism. *Arch Intern Med* 1997;157:2309-2316.

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© 2013, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim