„DIAQUICK“ Dengue IgG/IgM Ab Cassette for whole blood, serum and plasma samples

REF

Content

- Z06240**
- 30 cassettes individually packed + disposable pipettes (30 x Ref. No: Z06240B)
 - assay diluent sufficient for 30 tests
 - 1 package insert

For In-vitro diagnostic use only

INTENDED USE

The DIAQUICK Dengue IgG/IgM Ab Cassette is a solid phase immunochromatographic assay for the rapid, qualitative and differential detection of IgG and IgM antibodies to dengue virus in human serum, plasma or whole blood. This test is intended for professional use as an aid in the presumptive diagnosis of primary and secondary dengue infection. This test provides only a preliminary test result. Therefore, more specific alternative diagnosis methods (isolation of the virus, antigen detection in fixed tissues, RT-PCR, serological tests like HA-inhibition test) must be used in order to obtain a confirmation of dengue virus infection.

SUMMARY

Dengue fever is one of the most important mosquito-borne diseases in the world in the terms of morbidity and mortality. Dengue fever virus (serotypes 1 – 4) belongs to the group of flaviviridae, and is transmitted in nature by day-biting *Aedes* mosquitoes. The most important mosquito vector is the species *Aedes aegypti*. Primary Dengue infection, also known as Dengue Fever, is the most common type of dengue illness. It is associated with mild to high fever, headache, muscle pain and skin rash. Secondary infection is known as Dengue Hemorrhagic Fever (DHF) or Dengue Shock Syndrome, and often results in high fever and in many cases in hemorrhagic events and circulatory failure. The fatality rate in patients with Dengue Shock Syndrome can be as high as 44%. Dengue typically presents as a fever of sudden onset with headache, retrobulbar pain, pain in the back and limbs (break-bone fever), lymphadenopathy and maculopapular rash. Patients diagnosed with dengue in endemic areas generally have secondary infection, whereas patients in non-endemic areas are usually diagnosed with primary infection. Specific antibody responses to Dengue virus enable serodiagnosis and differentiation between primary and secondary dengue infections.

The DIAQUICK Dengue IgG/IgM Ab Cassette is a rapid immunochromatographic test using recombinant dengue viral antigens of all four serotypes to detect a specific antibody response.

TEST PRINCIPLE

The DIAQUICK Dengue IgG/IgM Ab Cassette utilizes the principle of immunochromatography. Mouse anti-human IgM and anti-human IgG antibodies are immobilized on the nitrocellulose membrane as two individual test lines (IgM line and IgG line) in the test window of the device. The IgG line in the test window is closer to the sample well and followed by IgM line. As the test sample flows along the membrane of the test device the colored Dengue specific recombinant antigen-colloidal gold conjugate agglutinates with the specific antibodies (IgM and/or IgG) of Dengue virus, if present in the sample. This complex moves further on the membrane to the test region where it is captured by the anti-human IgM and/or anti-human IgG antibodies coated on the membrane, leading to the formation of a colored band, which indicates a positive test results. Absence of this colored band in the test window indicates a negative test result. A built-in control line will always appear in the test window when the test has performed properly, regardless of the presence or absence of anti-Dengue virus antibodies in the specimen.

MATERIALS PROVIDED

The DIAQUICK Dengue IgG/IgM kit contains following items to perform the assay

- 1) Dengue IgG/IgM Cassette in foil pouch
- 2) Sample Buffer
- 3) Instructions for Use
- 4) Disposable pipette

MATERIALS REQUIRED BUT NOT PROVIDED

- 1) Specimen collection container
- 2) Timer

STORAGE AND STABILITY

The sealed pouches in the test kit may be stored between 4-30°C for the duration of the shelf life as indicated on the pouch.

PRECAUTIONS

- This kit is for in vitro diagnostic use only.
- This kit is for PROFESSIONAL use only.
- Read the instructions carefully before performing the test.
- This product does not contain any materials of human origin.
- Do not use kit contents after the expiration date.
- Handle all specimens as potentially infectious.
- Follow standard lab procedures and biosafety guidelines for handling and disposal of potentially infective material. When the assay procedure is completed, dispose specimens after autoclaving them at 121° C for at least 20 min. Alternatively, they can be treated with 0.5% Sodium Hypochlorite for 1-2 hours before disposal.
- Do not pipette reagents by mouth
- Don't smoke or eat while performing assays.
- Wear gloves during the whole procedure.


SPECIMEN COLLECTION AND PREPARATION

1. No special prior preparation of the patient is required before sample collection by approved techniques.
2. The test works best on fresh whole blood / serum / plasma samples. If testing cannot be performed immediately, serum / plasma may be stored at 2-8°C up to 3 days. For long-term storage, serum / plasma specimens can be frozen at -20°C for 3 months or at -70°C for longer period. Blood samples collected with a suitable anticoagulant such as EDTA, Heparin or Oxalate may be stored at 2-8°C for up to 3 days. Blood samples should not be frozen.
3. Repeated freezing and thawing of the specimen should be avoided.
4. Do not use haemolysed, clotted, contaminated, lipemic and viscous/turbid specimen.
5. Specimen containing precipitates or particulate matter must be centrifuged and only the clear supernatant used for testing.
6. Do not inactivate the sample by heating.
7. Shipment of specimens should comply with local regulations for transportation of etiologic agents.

TEST PROCEDURE

1. Bring the kit components to room temperature before testing.
2. Open the pouch and remove the card and pipet. Once opened, the test card must be used immediately.
3. Label the test card with patient identity.
4. Depress the pipet near the middle of the pipet to draw serum, plasma or whole blood sample up to the black 5 µL mark and apply the sample to the "S1" area.
5. Add 2 drops of sample buffer to the well marked as "S".
6. At the end of 20 minutes read the results. A strong positive sample may show results earlier.

Note: Result after 20 minutes may not be accurate

INTERPRETATION OF THE TEST

- 1) **Negative**
Only the control line appears.
- 2) **Positive for IgG/IgM**
Control line and both test lines appear. It indicates the possibility of acute secondary infection.
- 3) **IgM Positive, IgG Negative**
Both control line and the first test line (the higher one) appears. It indicates the possibility of primary infection.
- 4) **IgM Negative, IgG Positive**
Both control line and the second test line (the lower one) appears. It indicates the possibility of secondary infection or past infection.
- 5) **Invalid**
The test result is invalid if a colored band does not form in the control region. The sample must be re-tested, using a new test device.

QUALITY CONTROL

1. The control band is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagents are reactive.
2. Good Laboratory Practice recommends the daily use of control materials to validate the reliability of the device. Control materials which are not provided with this test kit may be commercially available.

LIMITATION OF THE TEST

1. The test is for qualitative detection of anti-Dengue antibodies in human serum, plasma or whole blood samples and does not indicate the quantity of the antibodies.
2. The test is for *in vitro* diagnostic use only.
3. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test but should rather be made after all the clinical findings have been evaluated.

PERFORMANCE CHARACTERISTICS

1) Accuracy

A panel of 60 positive and 40 negative patient sera was tested with a reference rapid Dengue Test, Dengue Duo Cassette (Panbio, Australia). The results are summarized in the following table. The agreement is 100%.

Dengue IgG/IgM Cassette	Dengue Duo Cassette	
	Positive	Negative
IgM Positive	35	0
IgG Positive	12	0
IgM/IgG Pos.	13	0
Negative	0	40

Sensitivity = 60/60 = 100%
 Accuracy = 100/100 = 100%


Specificity = 40/40 = 100%

2) Cross Reactivity

No cross reactivity with bilirubin (10 mg/dL), hemoglobin (18mg/dL) or triglycerides (up to 600 mg/dL) have appeared.

BIBLIOGRAPHY

- 1) Halstead, S.B. (1981), The pathogenesis of Dengue. Amer. J. Epidemiol 114: 632.
- 2) Henchal, E. A. and Putnuk, R. J., The Dengue viruses, Clin. Micro. Rev., Oct. 376 – 396, 1990.
- 3) Advances in Dengue Diagnosis, Maria G. Guzman, Gustavo Kouri. Clinical and Diagnostic Laboratory Immunology, Nov 1996, Vol. 3, No.6, p. 621-627.
- 4) Clinical Evaluation of a rapid immunochromatographic test for the diagnosis of Dengue Virus Infection. Chew Theng Sang, Lim Siew Hoon, Andrea Cuzzubbo, Peter Devine. Clinical and Diagnostic Laboratory Immunology, May 1998, Vol. 5, No. 3 p. 407-409.
- 5) Dengue and Dengue Hemorrhagic Fever, Duane J. Gubler. Clinical Microbiology Reviews, July 1998, Vol. 11, No. 3, p. 480-496.
- 6) Immunoglobulin A-specific Capture Enzyme-Linked Immunosorbent Assay for Diagnosis of Dengue Fever, Antoine Talarmin, Bhety Labeau, Josiane Lelarge, Jean-Louis Sarthou. Journal of clinical Microbiology, May 1998, Vol. 36, No. 5, p.1189-1192.
- 7) Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization 1997.
- 8) Hematological observations as diagnostic markers in dengue hemorrhagic fever – a reappraisal, Sunil Gomber, V.G. Ramachandran, Satish Kumar, K.N. Agarwal, P. Gupta, Piyush Gupta and D.K. Dewan. Indian Pediatrics 2001:38: 477-481.

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- 1) Halstead, S.B. (1981), The pathogenesis of Dengue. Amer. J. Epidemiol 114: 632.
- 2) Henchal, E. A. and Putnuk, R. J., The Dengue viruses, Clin. Micro. Rev., Oct. 376 – 396, 1990.
- 3) Advances in Dengue Diagnosis, María G. Guzman, Gustavo Kouri. Clinical and Diagnostic Laboratory Immunology, Nov 1996, Vol. 3, No.6, p. 621-627.
- 4) Clinical Evaluation of a rapid immunochromatographic test for the diagnosis of Dengue Virus Infection, Chew Theng Sang, Lim Siew Hoon, Andrea Cuzzubbo, Peter Devine. Clinical and Diagnostic Laboratory Immunology, May 1998, Vol. 5, No. 3 p. 407-409.
- 5) Dengue and Dengue Hemorrhagic Fever, Duane J. Gubler. Clinical Microbiology Reviews, July 1998, Vol. 11, No. 3, p. 480-496.
- 6) Immunoglobulin A-specific Capture Enzyme-Linked Immunosorbent Assay for Diagnosis of Dengue Fever, Antoine Talarmin, Bhety Labeau, Josiane Lelarge, Jean-Louis Sarthou. Journal of clinical Microbiology, May 1998, Vol. 36, No. 5, p.1189-1192.
- 7) Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization 1997.
- 8) Hematological observations as diagnostic markers in dengue hemorrhagic fever – a reappraisal, Sunil Gomber, V.G. Ramachandran, Satish Kumar, K.N. Agarwal, P. Gupta, Piyush Gupta and D.K. Dewan. Indian Pediatrics 2001:38: 477-481.