Biosyn™ Glycomer™ 631, Monofilament Synthetic Absorbable Sutures

Description: Biosyn™ synthetic absorbable sutures are prepared from a synthetic polyester, Glycomer™ 631,

composed of glycolide (60%), dioxanone (14%) and trimethylene carbonate (26%).

Characteristics: The advanced extrusion process of the molecule of Glycomer™ 631 gives the suture:

- Excellent initial strength

- Excellent in-vivo strength over the critical wound healing period

Excellent knot security
 Minimal memory

Indications: Biosyn™ synthetic absorbable sutures are indicated for use in general soft tissue approximation and/or

ligation including use in ophthalmic surgery, but not for use in cardiovascular or neurological surgery.

Actions: Biosyn™ synthetic absorbable sutures elicit a minimal acute inflammatory reaction in tissues, followed

by gradual encapsulation of the suture by fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of Biosyn™ sutures occurs by means of hydrolysis, where the suture is broken down to glycolic acid, dioxanoic acid, propane diol and carbon dioxide which are subsequently

absorbed and metabolized by the body.

Tensile Strength: Subcutaneous tissue implantation studies in rats indicate tensile strength for Biosyn™ sutures are

approximately 75% of USP and EP minimum knot strength at two weeks and approximately 40% at

three weeks post-implant.

Absorption: Absorption begins as a loss of tensile strength without appreciable loss of mass. Absorption of

Biosyn™ sutures is essentially complete between 90 and 110 days.

Sterilization: Ethylene oxide.

Packaging: Biosyn™ sutures are available in USP and EP sizes 1 (4 metric) through 6-0 (0.7 metric). They are

available undyed (natural) or violet colored. The sutures are supplied sterile, in pre-cut lengths or ligating reels; non-needled or affixed to various needle types, using both permanent and removable

needle attachment techniques; in one, two or three dozen boxes.



Biosyn™ Glycomer™ 631, Monofilament Synthetic Absorbable Sutures

Taper Point									
	Description	Size	Size (Metric)		Color	Units per box		Phasing-out	Reorder code
1/2 CIRCLE 22 mm		40	1.5	75	Violet	36	Surgalloy™		GM281
	⊙	4-0	1.5	90	Violet	36	Surgalloy ^{the}		GM291
	OV.25	34	2	90	Violet	36	Surgal loy TM		GM282
1/2 CIRCLE 26 mm	· · · · · · · · · · · · · · · · · · ·	40	1.5	5x45	Violet	12	D-Tach TM		GM64M GM321
	V-20	40	1.5	75	Undyed	36	Surgalloy™ Surgalloy Optivis™		GM111
		40	1.5	75	Violet	36	Surgalloy***		GM121
		34	2	5x45	Undyed	12	D-Tach™		GM69M
		34	2	5x45	Violet	12	D-Tach™		GM63M
		3-0	2	75	Undyed	36	Surgalloy™		GM322
		3-0	2	75	Violet	36	Surgalloy TM		GM122
		3-0	2	5x75	Violet	12	D-Tach™		GM67M
		2-0	3	75	Undyed	36			GM323
		2-0	3	75	Violet	36			GM123
		2-0	3	90	Violet	36			GM125
1/2 CIRCLE 26 mm	_	40	1.5	75	Violet	36	Surgalloy™		GM421
	⊙	4-0	1.5	90	Violet	36	Surgalloy TM		GM432
	V-20	30	2	90	Violet	36	Surgalloy™		GM332