



PRODUCT DATASHEET

Particulate Respirator, FFP2 with head straps

REF ZD-002

- ▶ "NR" marked (Non Reusable): the use of this filtering half-face mask is limited to one single shift. It is not reusable.
- ▶ Double elastic straps (adjustable) for an improved fitting time and added safety.
- ▶ Flexible nose-bridge with foam strip for comfort and safety
- ▶ Latex Free, PVC Free, Silicone free
- ▶ Flat-fold respirator



Product Info

FFP2 Particulate Respirator without valve (double head straps).

Product description	REF	Pack count
Particulate Respirator, FFP2 with head straps, colour white	ZD-002	20/600pcs
Sales unit: Carton		

Instruction for use:



Take out the respirator and pull on adjustable buckle with wire or other small gadget



Place respirator over face with the nose clamp sitting upward as shown. Position both head straps behind your head, making sure your lower chin is covered.



Hold the nose clamp to the shape of your nose to ensure it fits your face comfortably.



Adjust the respirator with adjustable buckle to a proper position and ensure there is no leakage.



To test fit: cup both hands over the respirator and inhale vigorously, a negative pressure should be felt inside the respirator. If air leaks at the edge, reposition the mask.

Technical Summary

Regulations:	The product covered by this data sheet does not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Product standards:	Personal Protective Equipment Regulation (EU) 2016/425 Annex 9. EN 149:2001+A1:2009 FFP2 NR The products are manufactured in an ISO 13485:2016 certified factory.
Main material:	First layer (outer): Spunbond polypropylene Second layer (filtration layer): Meltblown polypropylene Third layer (liner): Hot-air Nonwoven Fourth layer (inner): Spunbond polypropylene Head-straps: Spandex Head-strap adjusting buckles: Silicone Nose clamp: Polypropylene and iron
Pack count:	20 pcs/PE bag/box, 30 boxes/carton
Transport carton:	650 * 425 * 320 mm
Storage recommendation:	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight and other intensive light sources.
Disposal recommendation:	Follow local instructions
Shelf Life:	3 years

Manufacturer:
Zhende Medical Co., Ltd.
Gaobu Town, 312035,
Shaoxing City, Zhejiang
Province, China
Distributor:
OneMed Sverige AB
www.onemed.se