

SELEFA® Optimia Nonwoven swabs

REF: 223001, 223002, 223003, 223004

SELEFA® Optimia Nonwoven swabs, 4-ply are produced from soft and highly absorptive nonwoven material. They are primarily used for cleaning and exudate absorption of wounds and abrasions.

Key Features

- ▶ High quality viscose/polyester nonwoven material
- ▶ Low linting and high softness
- ▶ Keeps its shape, even when wet
- ▶ Free of binding agents and optical brighteners
- ▶ Effective in liquid and exudate absorption
- ▶ Less adherent to the wound than gauze devices
- ▶ Available in several sizes
- ▶ For external wound care, treating burns or abrasions, in the first aid area etc.
- ▶ Can be sterilized using EO
- ▶ Optimized packaging: convenient dispensing and better protection from contamination



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
223001	Optimia Nonwoven swabs, 4-ply, non-sterile	5x5 cm	150 / 13500
223002	Optimia Nonwoven swabs, 4-ply, non-sterile	7,5x7,5 cm	150 / 9000
223003	Optimia Nonwoven swabs, 4-ply, non-sterile	10x10 cm	150 / 5400
223004	Optimia Nonwoven swabs, 4-ply, non-sterile	10x20 cm	75 / 900

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices.



REACH regulation: The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-05
EN ISO 10993-10

Compliance with legal requirements: Medical Device Regulation (EU) 2017/745, MDR class I

Compliance with product standards: EN 1644-1, -2, EN ISO 11737-1

Quality standards: EN ISO 13485

Label information and packaging compliance: EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristics
Main material	Nonwoven, 70% viscose and 30% polyester, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Carton, TCF or ECF quality	223001: 230 x 445 223002: 157 x 67 x 53 223003: 232 x 67 x 77 223004: 204 x 92 x 104
Transport Carton	Carton, TCF or ECF quality	223001: 410 x 410 x 350 223002: 490 x 280 x 420 223003: 485 x 405 x 420 223004: 430 x 330 x 570
Packing Tape	PVC-free, Transparent	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
223001	N/A	6438129353976	6438129553970
223002	N/A	6438129353983	6438129553987
223003	N/A	6438129353990	6438129553994
223004	N/A	6438129354003	6438129554007

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.