

SPONGOSTAN™ Absorbable Gelatin Sponge

Gelatin sponge for hemostatic use

SPONGOSTAN™ Absorbable Haemostatic Gelatin Sponge is a sterile, water-insoluble, malleable, porcine gelatin absorbable sponge intended for haemostatic use by applying to a bleeding surface.

When used in appropriate amounts SPONGOSTAN™ Sponge is absorbed completely within 4 to 6 weeks.

When applied to bleeding mucosal regions, it liquefies within 2 to 5 days.



Product	Size	Units per box	Code
SPONGOSTAN™ Film	20 x 7 x 0,05 cm	20	MS0001
SPONGOSTAN™ Standard	7 x 5 x 1 cm	20	MS0002
SPONGOSTAN™ Special	7 x 5 x 0,1 cm	20	MS0003
SPONGOSTAN™ Anal	8 x 3 cm	20	MS0004
SPONGOSTAN™ Dental	1 x 1 x 1 cm	24	MS0005
SPONGOSTAN™ Powder	1 g	6	MS0008

