

# **GLORIOUS ENTERPRISES**

## **PRODUCT SHEET**

**Products Name:** PEAN Artery Forceps 16cm straight



**Glorious Art nr.** 43.0604.16  
**Net weight:** 43 gram (approximate)  
**Material:** Stainless steel AISI 410  
**Sterility:** Non-Sterile (Re-Usable)  
**Exposure Temperature:** 134°C (274°F)

**Exposure Time:** 3 Minutes

**CE Classification:** Class-I

**Storage:** 16°C to 30°C (Store the instruments in a closed cabinet, protected from dust, moisture and temperature fluctuations)

**Shelf life:** 3 years

**life-span:** 20 times of use:- after each use, product need to sterilize for ready to use in next operation or next surgery.

**Disposal:** No special decomposition or disposal required for Surgical Instruments as they don't contain any toxic or hazardous material