

GLORIOUS ENTERPRISES

PRODUCT SHEET

Products Name: PEAN Artery Forceps 14cm curved



Glorious Art nr. 43.0601.14

Net weight: 35 gram (approximate)
Material: Stainless steel AISI 410
Sterility: Non-Sterile (Re-Usable)

Exposure Temperature: 134°C (274°F)

Exposure Time: 3 Minutes

CE Classification: Class-I

Storage: 16°C to 30°C (Store the instruments in a closed cabinet, protected

from dust, moisture and temperature fluctuations)

Shelf life: 3 years

life-span: 20 times of use:- after each use, product need to sterilize for

ready to use in next operation or next surgery.

Disposal: No special decomposition or disposal required for Surgical Instruments

as they don't contain any toxic or hazardous material



P.O. Box 63, Naverland 11, DK 2600 Glostrup, Denmark Tel: +45 43 62 26 93, Fax: +45 43 62 26 11, Email: glorious@glorious-enterprises.dk