

Quality time for better care

"Quality time for better care" is the brand promise of Terumo Medical Care Solutions

**HOSPITAL CARE**

# Terumo Syringe

3-part Syringe  
Product Information



Every single detail is the result of Terumo's long time experience in this field which is appreciated by users throughout the world.

- Siliconisation for optimal plunger movement
- Firm and definite end-stop
- Clear scale markings to enable precise dosage
- Transparent barrel for visualisation
- Finger grips support easy handling

Description	Raw Materials
<b>1</b> Barrel	Polypropylene
<b>2</b> Plunger	Polypropylene
<b>3</b> Gasket	Thermoplastic elastomer (TPE) Polyisoprene rubber (PI) MDSS-10000
Lubricant	Silicone oil

Does not contain components made of PVC or natural rubber latex.



### Sterilization

With radiation (e-beam) (for ,SS'-codes) or ethylene oxide (for ,MDSS'-codes) to a Sterility Assurance Level of at least  $10^{-6}$



### Storage

- Shelf life of 5 years
- The batch number and expiry date are printed on the paper (for ,MDSS'-codes) or embossed in the film of the blister pack (for ,SS'-codes)
- Do not store at extreme temperature and humidity

**Disposable syringes  
for single use**

# Product Information Terumo Syringe



## Product Range

Product Code	Nominal Capacity	Tip	Replacement for
MDSS01SE	1 ml	Concentric Luer Slip	SS+01T1
SS*02LE1	2.5 ml	Concentric Luer Lock	-
SS*02SE1	2.5 ml	Concentric Luer Slip	-
MDSS03SE	3 ml	Concentric Luer Slip	SS*02SE1, SS+03S
MDSS03LE	3 ml	Concentric Luer Lock	SS*03L, SS+03L1
MDSS05SE	5 ml	Concentric Luer Slip	SS*05SE1, SS+05S
MDSS05LE	5 ml	Concentric Luer Lock	-
SS*05LE1	5 ml	Concentric Luer Lock	-
SS*05SE1	5 ml	Concentric Luer Slip	-
MDSS10SE	10 ml	Concentric Luer Slip	SS+10ES1
MDSS10ESE	10 ml	Eccentric Luer Slip	SS+10FS1
MDSS10LE	10 ml	Concentric Luer Lock	-
SS*10LE1	10 ml	Concentric Luer Lock	-
MDSS20ESE	20 ml	Eccentric Luer Lock	SS+20ES1
SS+20L1	20 ml	Concentric Luer Slip	-
SS+30L1MP	30 ml	Concentric Luer Lock	-
SS+30ES1	30 ml	Eccentric Luer Slip	-
SS+50L1	50 ml	Concentric Luer Lock	-
SS+50ES1	50 ml	Eccentric Luer Slip	-
SS+50C1	50 ml	Catheter Tip	-



STERILE EO for MDSS-codes  
STERILE R for SS-codes

## Precautions

- Do not use if package is damaged
- Do not reuse
- Dispose of sharps after single use according to facility protocol to avoid risk of infection

Always consult the product label for a complete overview of warnings, cautions and/or precautions prior to actual use.

Choose from a range of volumes from 1 ml to 50 ml  
Syringes are available in all common tip types.

## Packaging

Unit Box	Shipping Carton	
	Quantity	Gross Weight
100	3000	16.60 kg
100	2400	6.00 kg
100	2400	6.00 kg
100	1800	10.90 kg
100	1800	10.90 kg
100	1800	14.90 kg
100	1800	13.15 kg
100	1600	6.25 kg
100	1600	6.25 kg
100	1200	14.60 kg
100	1200	14.60 kg
100	1200	13.03 kg
100	800	15.00 kg
50	400	14.60 kg
50	600	10.15 kg
25	200	6.32 kg
25	200	6.20 kg
25	100	5.10 kg
25	100	5.10 kg
25	100	5.40 kg

STERILE EO for MDSS-codes  
STERILE R for SS-codes

## Regulatory Status

CE Certification –  
Class IIa medical device

**Notified body**  
TÜV Rheinland 0197

**Manufacturer**  
for SS-codes TERUMO (Philippines)  
CORPORATION, 124 East Main Avenue,  
Laguna Technopark, Binan, Philippines  
for SS-codes TERUMO CORPORATION,  
2-44-1 Hatagaya, Shibuya-Ku,  
Toyko 151-0072, Japan  
for MDSS-codes Yangzhou Medline  
Industry Co., Ltd. (China)

## EC representative

for SS-codes TERUMO EUROPE N.V. –  
Leuven, Belgium  
for MDSS-codes Shanghai International  
Corp. GmbH (Europe) –  
Hamburg, Germany

## Distributor

for MDSS-codes TERUMO EUROPE N.V.,  
Interleuvenlaan 40, 3001 Leuven, Belgium



**TERUMO EUROPE N.V.**

Terumo Corporation Terumo Europe N.V.  
+81 3 6742 8500 +32 16 38 12 11

www.terumo-europe.com

All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners. Published by Terumo Europe N.V.  
Please note that some products may not be available in all regions, check with your local Terumo contact.