


# neria™

## Intended use:

The product is intended for the subcutaneous infusion of medication such as immunoglobulins for the treatment of Primary Immune Deficiency, deferoxamine mesylate for Thalassaemia, apomorphine for Parkinson's Disease and morphine for pain management therapy. The neria™ infusion set has been tested to be compatible with the above-mentioned drugs approved for subcutaneous infusion.



## Product range:

Product code	Description	Needle gauge	Needle length	Tubing length	EAN	Pcs per inner box	Pcs per outer box
<b>neria™</b>							
78-060-2731	neria™ infusion set	G27	10 mm	60 cm	05705244000513	10	200
78-060-2738	neria™ infusion set	G27	8 mm	60 cm	05705244000193	10	200
78-060-2736	neria™ infusion set	G27	6 mm	60 cm	05705244000476	10	200
78-080-2737	neria™ infusion set	G27	12 mm	80 cm	05705244004191	10	200
78-080-2731	neria™ infusion set	G27	10 mm	80 cm	05705244001077	10	200
78-080-2738	neria™ infusion set	G27	8 mm	80 cm	05705244000575	10	200
78-080-2736	neria™ infusion set	G27	6 mm	80 cm	05705244000551	10	200
78-110-2731	neria™ infusion set	G27	10 mm	110 cm	05705244000117	10	200
78-110-2738	neria™ infusion set	G27	8 mm	110 cm	05705244000292	10	200
78-110-2736	neria™ infusion set	G27	6 mm	110 cm	05705244000599	10	200
78-060-2931	neria™ infusion set	G29	10 mm	60 cm	05705244000216	10	200
78-060-2938	neria™ infusion set	G29	8 mm	60 cm	05705244000131	10	200
78-060-2936	neria™ infusion set	G29	6 mm	60 cm	05705244000612	10	200
78-080-2931	neria™ infusion set	G29	10 mm	80 cm	05705244000711	10	200
78-080-2938	neria™ infusion set	G29	8 mm	80 cm	05705244000674	10	200
78-080-2936	neria™ infusion set	G29	6 mm	80 cm	05705244000636	10	200
78-110-2931	neria™ infusion set	G29	10 mm	110 cm	05705244000230	10	200
78-110-2938	neria™ infusion set	G29	8 mm	110 cm	05705244000155	10	200
78-110-2936	neria™ infusion set	G29	6 mm	110 cm	05705244000759	10	200

**Product specifications:**

Insertion angle . . . . .90 degrees  
 Needle length . . . . . 6 mm, 8 mm, 10 mm and 12 mm  
 Needle gauge . . . . .27 gauge and 29 gauge  
 Tubing length in cm . . . . .60, 80 and 110 cm  
 Tubing inner diameter . . . . . 0,385 mm  
 Tubing outer diameter . . . . . 1,50 mm

**Material specifications:**

PVC free . . . . . Yes  
 Made with natural rubber latex . . . . . No  
 Raw materials or derivation of animal origin . . . . . No  
 The product contains phthalates . . . . . No  
 Introducer & connector needles may contain . . ≤ 10% nickel

**Certifications:**

CE-marking . . . . . Yes  
 Product classification . . . . . Class IIB (EU)  
 Quality system certified to . . . . . ISO 13485  
 Tariff code . . . . . HS9018390000  
 GMDN . . . . . 35833






**Sterilization:**

Sterile . . . . . Yes  
 Sterilized by . . . . . Ethylene oxide  
 Shelf life . . . . . 5 years  
 Usage time (drug depended) . . . . . Maximum 48 hours

**Drug compatibility/stability tested with:**

Immunoglobulins  
 Apomorphine  
 Morphine (hydromorphone, morphine sulphate & morphine chloride)  
 Desferal

**Storage and handling:**

	<b>AVOID DIRECT SUNLIGHT</b> - Keep the product away from direct sunlight and any humid environment
	<b>KEEP DRY</b> - Store in a dry place at room temperature
	<b>DO NOT RE-USE</b>
	<b>DO NOT USE</b> - If package is damaged
	<b>CAUTION</b> - See Instruction for Use