

Product Description

Item No.:	1131897
Item:	Alcohol swabs
Material:	Non-woven material saturated with 70% Isopropyl alcohol
Sealing:	Aluminium foil paper pouches
Size:	Folded: appr. 3 x 3 cm Unfolded: appr. 6 x 3 cm
Intended Use:	For skin cleansing as well as for many jobs in the laboratory. For single use only!
Instruction for use:	Swipe over the skin in circular movements. CAUTION: Only for external use: in the event of skin irritation or rash, discontinue use. Avoid contact with eyes; if contacted, immediately flood the eye with plenty of water. If ingested, DO NOT induce vomiting, obtain medical attention immediately.
LOT-No:	The LOT-No. and expiry date EXP in format DD (day) MM (month) YY (year) are printed on each pouch, on each box and on the outer cartons
Shelf life:	5 years from date of production
Packaging:	100 pcs. per box 2,000 pcs. per carton



© www.herenz.de

www.herenz.de

Phone +49-040-739-20-40
Fax +49-040-730-20-439
Email hhh@herenz.de
ILN 40 49434 00000 8
PD_Alkoholtupfer-e.docx

Managing Directors

Dieter Cyll, Antje Herenz-Cyll
Comm. Reg. Hamburg HRB 48884
Tax No. 44/789/02888
V.A.T. No. DE 118 665 646

Postbank Hamburg

Sort Code 200 100 20
Account No. 6846 46 200
IBAN DE97200100200684646200
Swift-Code PBNKDEFF

Deutsche Bank AG

Sort Code 200 700 00
Account No. 55 00 251
IBAN DE82200700000550025100
Swift-Code DEUTDEHH

Hamburger Sparkasse

Sort Code 200 505 50
Account No. 1085 243 002
IBAN DE87200505501085243002
Swift-Code HASPDEHH

08.05.2020 - rev. 2