

BD Nexiva™ Closed IV Catheter System Product configurations and specifications

BD Nexiva™ Closed IV Catheter System – Single Port with BD MaxZero™ Needle-free Connector



Catalog no.	Color	Gauge	Catheter length (in)	Catheter ID in (mm)	Catheter OD in (mm)	Extension tube ID (mm)	Gravity flow rate (mL/min)	Maximum power injector flow rate at 300 PSI (mL/sec)	Packaging
383551	● Yellow	24	0.75/ 19	0.0210 (0.53)	0.0280 (0.71)	1.22	18	Not for use with power injectors	20/box 80/case
383552	● Blue	22	1.00/ 25	0.0265 (0.67)	0.0355 (0.90)	1.22	33	3.0	20/box 80/case
383553	● Blue	22	1.75/ 45	0.0265 (0.67)	0.0355 (0.90)	1.22	30	3.0	20/box 80/case
383556	● Pink	20	1.00/ 25	0.0325 (0.83)	0.0435 (1.10)	1.65	61	4.0-5.5	20/box 80/case
383557	● Pink	20	1.25/ 32	0.0325 (0.83)	0.0435 (1.10)	1.65	58	4.0-5.5	20/box 80/case
383558	● Pink	20	1.75/ 45	0.0325 (0.83)	0.0435 (1.10)	1.65	51	4.0-5.5	20/box 80/case
383559	● Green	18	1.25/ 32	0.0385 (0.98)	0.0515 (1.31)	1.65	84	4.0-7.0	20/box 80/case
383560	● Green	18	1.75/ 45	0.0385 (0.98)	0.0515 (1.31)	1.65	79	4.0-7.0	20/box 80/case

BD Nexiva™ Closed IV Catheter System – Single Port without Needle-free Connector



Catalog no.	Color	Gauge	Catheter length (in)	Catheter ID in (mm)	Catheter OD in (mm)	Extension tube ID (mm)	Gravity flow rate (mL/min)	Maximum power injector flow rate at 300 PSI (mL/sec)	Packaging
383511	● Yellow	24	0.75/ 19	0.0210 (0.53)	0.0280 (0.71)	1.22	18	Not for use with power injectors	20/box 80/case
383512	● Blue	22	1.00/ 25	0.0265 (0.67)	0.0355 (0.90)	1.22	33	3.0	20/box 80/case
383513	● Blue	22	1.75/ 45	0.0265 (0.67)	0.0355 (0.90)	1.22	30	3.0	20/box 80/case
383516	● Pink	20	1.00/ 25	0.0325 (0.83)	0.0435 (1.10)	1.65	61	4.0-5.5	20/box 80/case
383517	● Pink	20	1.25/ 32	0.0325 (0.83)	0.0435 (1.10)	1.65	58	4.0-5.5	20/box 80/case
383518	● Pink	20	1.75/ 45	0.0325 (0.83)	0.0435 (1.10)	1.65	51	4.0-5.5	20/box 80/case
383519	● Green	18	1.25/ 32	0.0385 (0.98)	0.0515 (1.31)	1.65	84	4.0-7.0	20/box 80/case
383520	● Green	18	1.75/ 45	0.0385 (0.98)	0.0515 (1.31)	1.65	79	4.0-7.0	20/box 80/case



BD Nexiva™ Closed IV Catheter System – Dual Port with BD Q-Syte™ Needle-free Connector



Catalog no.	Color	Gauge	Catheter length (in)	Catheter ID in (mm)	Catheter OD in (mm)	Extension tube ID (mm)	Gravity flow rate (mL/min)	Maximum power injector flow rate at 300 PSI (mL/sec)	Packaging
383531	● Yellow	24	0.75/ 19	0.0210 (0.53)	0.0280 (0.71)	1.22	18	Not for use with power injectors	20/box 80/case
383532	● Blue	22	1.00/ 25	0.0265 (0.67)	0.0355 (0.90)	1.22	33	3.0	20/box 80/case
383536	● Pink	20	1.00/ 25	0.0325 (0.83)	0.0435 (1.10)	1.65	61	4.0-5.5	20/box 80/case
383537	● Pink	20	1.25/ 32	0.0325 (0.83)	0.0435 (1.10)	1.65	58	4.0-5.5	20/box 80/case
383538	● Pink	20	1.75/ 45	0.0325 (0.83)	0.0435 (1.10)	1.65	51	4.0-5.5	20/box 80/case
383539	● Green	18	1.25/ 32	0.0385 (0.98)	0.0515 (1.31)	1.65	84	4.0-7.0	20/box 80/case
383540	● Green	18	1.75/ 45	0.0385 (0.98)	0.0515 (1.31)	1.65	79	4.0-7.0	20/box 80/case

BD Nexiva™ Diffusics™ Closed IV Catheter System



Catalog no.	Color	Gauge	Catheter length (in)	Extension tube ID (mm)	Gravity flow rate (mL/min)	Maximum power injector pressure limit setting (PSI)	Maximum power injector flow rate at 325 PSI (mL/sec) for contrast media viscosity 27.5 cP (mPa-s)*	Packaging
383691	● Yellow	24	0.75/ 19	1.22	21	325	3.00	20/box 80/case
383692	● Blue	22	1.00/ 25	1.65	45	325	6.50	20/box 80/case
383693	● Pink	20	1.00/ 25	1.65	68	325	10.0	20/box 80/case
383694	● Pink	20	1.25/ 32	1.65	64	325	10.0	20/box 80/case
383695	● Green	18	1.25/ 32	1.65	90	325	15.0 mL/sec	20/box 80/case

*The catheter system has been tested at the listed flow rates, however, due to variations in add-on devices, tubing, contrast media temperature, contrast media viscosity and pressure limit settings these flow rates may not be achievable.

Please consult product labels, inserts and instructions for use for indications, contraindications, warnings and precautions.

Becton Dickinson Sweden AB
 Marieviksgatan 25
 117 43 Stockholm
 Sweden

Becton Dickinson Denmark A/S
 Lyskær 3E, 1. sal
 2730 Herlev
 Denmark

Becton Dickinson Norway AS
 Gjerdrums vei 8
 0484 Oslo
 Norway

Becton Dickinson Finland Oy
 Äyritie 18
 01510 Vantaa
 Finland



bd.com

BD, the BD Logo, Diffusics, MaxZero, Nexiva and Q-Syte are trademarks of Becton, Dickinson and Company or its affiliates.
 © 2023 BD. All rights reserved. BD-85002 (03/23) BD-79505

