

Technical Data Sheet

Company Name:	Local company name	Telephone
	Cook Österreich GMBH	+43 17 956 7121
Cook Medical Europe LTD.	Cook Belgium BV/SL	+32 27001633
Address:	Cook Danmark ApS	+45 38 48 76 07
	Cook Medical Distribution Partners	+353 61 239 240
O' Halloran Road	Cook Medical Finland Oy	+358 9725 19996
National Technology Park	Cook Deutschland GMBH	+49 69 500 72 804
Limerick, Ireland	Cook Medical Hungary KFT.	+36 1 777 9199
	Cook Medical Europe LTD.	+353 61 239 252
	Cook Nederland BV	+31 202013367
Web:	Cook Norge AS	+47 23162968
www.cookmedical.com	Cook Medical SP. Z.O.O.	+48 22 3060159
	Cook Sweden AB	+46 85 876 94 68
	Cook Switzerland GmbH	+41 44 800 96 09
	Cook UK LIMITED	+44 207 3654 183

Vigilance Contact:

CMEU Regulatory Affairs -Post Market

E-mail:

PostMarketEUAuthorisedRep@cookmedical.com

Tel: +353 61 33 44 40

Product:

Connecting tube with drainage bag connector

Product Code:

CTU14.0-30-ST CTU14.0-50-ST

DIMENSIONS

CTU14.0-30-ST

- Size 14Fr
- Length 30cm


• Connections: Stopcock

CTU14.0-50-ST

• Size - 14Fr

• Length – 50cm

• Connections: Stopcock



Components:

Connecting tube
Drainage bag connector

Description:

The one way stopcock allows control of drainage flow during drainage pouch exchange

Intended use:

Used for connection of external drainage catheter to drainage pouch.

Contraindications:

None Known

Precautions:

None Known

Recommended storage instructions:Store in a dry location, away from temperature extremes

Patient contacting	N/A
materials	

CE Marked	Yes
Class	I
Notified body number	0123
Manufacturer	Cook Incorporated

Single use	Yes
Sterilized	Yes
Address of sterilization	6300 North Matthews Drive
	Ellettsville, IN 47429 USA
Method of sterilization	ETO
Shelf life	Up to 36 months
Biological origin	No
Latex	No
Biocompatibility	Yes
Quantity supplied	1