

# Qufora<sup>®</sup> Leg Bags | Single Chamber | Drainable

Q-Flow Valve or T-Valve | Adjustable or Fixed Tube Length | Single Packed

# Discreet. Comfortable. Easy-to-use.

The sterile **Qufora® leg bags** are designed to be worn discreetly under clothing and can be attached to either a urinary catheter or a urisheath. The Qufora® leg bags come with rounded corners for added comfort and are easily fixed to the leg using the Qufora® leg straps.

# **PRODUCT BENEFITS**

- The **single chamber** bags come with different tube lengths or volumes.
- The Qufora<sup>®</sup> leg bags have a soft adjustable or fixed inlet tube. The adjustable tube makes it possible to personalize the tube's length. One sterile leg bags also come with a separate connector/cap.
- A non-return valve ensures that urine will not flow back into the tube, reducing the risk of infection.
- The Qufora<sup>®</sup> leg bags have a universal conical inlet connector (X-mas tree), making it easy to attach a urinary catheter or urisheath to the bag.
- The front has a scale for quick estimation of the urine volume in the bag and the back has a skin friendly soft backing for increased comfort, even in warm conditions.
- The patented Q-Flow valve is easy, ergonomic and safe to operate with one hand also for people with reduced finger dexterity. When turning it 90°, the valve is fully open. At 180° the valve is still fully open and now placed in a safe position without any risk of closing unexpectedly during the night.
- The Q-Flow valve has a silicone tube for easy connection to a urine bag for the night, decreasing the risk of infection when changing the bag.
- The Q-Flow valve provides for a free-flow drainage of urine in a hygienic and safe manner.
- The T-Valve, with a white "stopper", is easy and safe to operate with one hand. It has a silicone tube for easy connection to a urine bag for the night, decreasing the risk of infection when changing the bag.
- The sterile Qufora<sup>®</sup> leg bags come with a **needle-free sample port** for easy urine sampling. The sample port is specially designed for easy cleaning.
- The Qufora<sup>®</sup> leg bags are single packed with the REF number, lot number and expiration year/month printed on the protective bag.

## Qufora® Accessories (not included):

The Qufora<sup>®</sup> leg straps are designed for the Qufora<sup>®</sup> leg bags. The straps (with two buttons and silicone anti-slip) are adjustable for individual needs and the bags include **holes** for the straps.

Medical Device
350 ml 600 ml 750 ml 800 ml
Non-sterile
Sterile

Data sheet



# Qufora® Leg Bags | Single Chamber | Drainable

Q-Flow Valve or T-Valve | Adjustable or Fixed Tube Length | Single Packed

# Product description

- Qufora<sup>®</sup> leg bags are available in several versions with different bag volumes or lengths of the inlet tube.
- The Qufora<sup>®</sup> leg bags are single-use medical devices.
  All materials are 100% latex- and DEHP-free.

# **QUFORA® LEG BAG FEATURES**

- Single chamber
- Adjustable or fixed inlet tube length with a frosted surface (inner diameter: 5.5mm) (11451503, 11451803, 11401803: inner diameter: 5.0mm)
- Separate or fixed connector with protective cap
- Non-return valve
- Conical inlet connector (X-mas tree)
- Non-woven soft backing
- Q-Flow valve with soft silicone outlet tube and free flow (outlet tube: inner diameter: 8.0mm)
- T-Valve with white "stopper" and soft silicone outlet tube (outlet tube: inner diameter: 8.0mm)

PE

PP/PE

Silicone

- Translucent and white films
- Scale on the front
- Single packed
- Eyelets for straps

## MATERIALS

#### Films

- Inlet tube
- Connector
- Connector (sample port)
- Connector cap
- Non-return valve
- Soft backing
- Q-Flow valve
- T-Valve/white "stopper"
- Outlet tube
- Single bag (product)
- Single bag (product)
- Poly bags

# PACKAGING

•	Retail box Shipper box IFU	Cardboard Corrugated cardboard One (1) per retail box or inner poly bag
•	Shelf-life Storage	Five (5) years Dry in temperatures between 0°C and 35°C; protected from direct sunlight; humidity level max 80%; not to be exposed to heavy pressure

PVC (plasticizer: DOTP/DEHT)

PVC (plasticizer: DOTP/DEHT)

PVC (plasticizer: DOTP/DEHT)

LDPE (non-sterile products)

LDPE (inner and outer bag)

HDPE & PE/PET (sterile products)

Non-woven PET

LLDPE/POM;

PP (bags with adjustable tube length)

Inlet tube: PVC (plasticizer: DOTP/DEHT)

Inlet tube: PVC (plasticizer: DOTP/DEHT)

PVC; Port: TPE (bags with fixed tube length)





# QUFORA® LEG BAGS - SINGLE CHAMBER - DRAINABLE NON-STERILE – Q-FLOW VALVE - ADJUSTABLE OR FIXED TUBE LENGTH - SINGLE PACKED

REF	Volume	Tube length	Connector	Soft backing	Quantity poly bag	Quantity shipper box	GTIN-14 <sup>5</sup>
11061901	350 ml	30 cm <sup>2</sup>	X-mas	Yes	104	100	15710187005850
11223802 <sup>1</sup>	600 ml	10 cm <sup>3</sup>	X-mas	No	25	100	15710187001203
11221902	600 ml	30 cm <sup>3</sup>	X-mas	Yes	25	100	15710187001234
11221502	600 ml	50 cm <sup>2</sup>	X-mas	Yes	25	100	15710187004181
11421502	800 ml	50 cm <sup>2</sup>	X-mas	Yes	25	100	15710187004167

<sup>1</sup>Slanting outlet tube; <sup>2</sup>Adjustable tube length; <sup>3</sup>Fixed tube length; <sup>4</sup>Retail box; <sup>5</sup>GTIN barcode number for retail box or poly bag

#### Non-sterile

# QUFORA® LEG BAGS - SINGLE CHAMBER - DRAINABLE - NON-DRAINABLE NON-STERILE - T-VALVE - ADJUSTABLE OR FIXED TUBE LENGTH - SINGLE PACKED

REF	Volume	Tube length	White "stopper"	Soft backing	Quantity poly bag	Quantity shipper box	GTIN-14 <sup>5</sup>
11251502	600 ml	50 cm <sup>3</sup>	Yes	Yes	25	100	15710187001173
11251902	600 ml	30 cm⁴	Yes	Yes	25	100	15710187001142
11259802 <sup>1</sup>	600 ml	10 cm <sup>4</sup>	Yes	No	25	100	15710187001111
11451503	750 ml	50 cm⁴	Yes	No	25	100	15710187004631
11451803	750 ml	10 cm <sup>4</sup>	Yes	No	25	100	15710187004662
11401803 <sup>2</sup>	750 ml	10 cm <sup>4</sup>	No	No	25	100	15710187004693

All bags with X-mas tree connector

<sup>1</sup>Slanting outlet tube; <sup>2</sup>Non-drainable bag with tear mark; <sup>3</sup>Adjustable tube length (loose connector/cap); <sup>4</sup>Fixed tube length; <sup>5</sup>GTIN barcode number for poly bag

### Sterile

# Needle-free sample port connector

# QUFORA® LEG BAGS – SINGLE CHAMBER - DRAINABLE STERILE - Q-FLOW VALVE OR T-VALVE - ADJUSTABLE OR FIXED TUBE LENGTH - SINGLE PACKED

REF	Volume	Tube length	Needle free sample port	Soft backing	Quantity retail box	Quantity shipper box	GTIN-14 <sup>6</sup>
12061901 <sup>1</sup>	350 ml	30 cm <sup>3</sup>	Yes⁵	Yes	10	100	15710187007038
12251902 <sup>2</sup>	600 ml	30 cm⁴	No <sup>5</sup>	Yes	10	100	15710187003078

<sup>1</sup>Q-Flow Valve; <sup>2</sup>T-Valve with white "stopper"; <sup>3</sup>Fixed length and adjustable option with separate connector/cap; <sup>4</sup>Fixed tube length; <sup>5</sup>X-mas tree connector; <sup>6</sup>GTIN barcode number for retail box

Non-sterile

# **QUFORA® LEG STRAPS**

REF	Description	Width / length	Quantity retail box	Quantity shipper box	GTIN-141
21100301	Velcro <sup>®</sup> leg straps 1 pair with silicone anti-slip and buttons	3 cm / 60 cm	10 pairs	100 pairs	15710187001715

Sterile

## **QUFORA® Q-FLOW VALVE**

REF	Description	Quantity retail box	Quantity shipper box	GTIN-141
21104201	Q-Flow valve with 5 cm soft silicone tube	30	300	15710187004754

## **QUFORA® T-VALVE**

н	REF	Description	Quantity retail box	Quantity shipper box	GTIN-14 <sup>1</sup>
	21104101	T-Valve with 5 cm soft silicone tube	30	300	15710187001807

## **QUFORA® LEVER VALVE**

	REF	Description	Quantity retail box	Quantity shipper box	GTIN-14 <sup>1</sup>
	21104301	Lever valve with 5 cm soft silicone tube	30	300	15710187001838

<sup>1</sup>GTIN barcode number for retail box

# **PRODUCT CLASSIFICATION**

 Medical Devices conforming to MDR 2017/745 in Class I/Class I sterile (Is)

# **APPLIED STANDARDS**

- ISO 13485 Quality Management System
- ISO 8669-2 Urine collection bags
  - Part 2: Requirements and test methods
- GMDN Urinebags: Non-sterile: 58924; Sterile: 58923 Straps: 38861; Valves (sterile): 31975

# STERILIZATION

Gas sterilized with ethylene oxide (EtO)



160720

MBH-International A/S Gydevang 28-30 DK-3450 Allerød Denmark

Importer

Tel: +45 88 30 05 92 info@qufora.com www.qufora.com