

# SANDEL® TIME OUT® SKIN MARKER

HEALTHCARE SAFETY DEVICES

## PRODUCT DESCRIPTION

<b>Intended Use</b>	The Gentian violet ink is used to mark an anatomical site for surgery.
<b>Material</b>	Polyester (PET) Polypropylene (PP) Gentian Violet Ink Ruler: paper
<b>Weight</b>	~0.28 ounce (~8 grams)
<b>Overall Dimensions (inch)</b>	• 5.5" L x 0.6" H (14 L x 1.5 H cm) • Ruler (1019-R) : 6" (~ 15 cm)

## PACKAGING AND STORAGE

	Packaging	Total/case
1019	25/Bx 4Bx/Cs	100
1019-R		

## MANUFACTURING AND SAFETY STANDARDS

<b>Lot Traceability</b>	Yes
<b>Re-Usable</b>	Disposable Product, Single Use Only
<b>Cleaning</b>	N/A
<b>Sterilisation Method</b>	Gamma Radiation
<b>Re-Sterilisation</b>	N/A
<b>CE Marked</b>	CE Classification: Class I sterile
<b>Patent Number</b>	US D538,851



## PRODUCT NUMBER

<b>Sterile (shelf life of 3 years)</b>	
1019	1019-R

➤ **Contact your Ansell representative for ordering or more information.**

**North America**  
Tel: (866) 764-3327  
Email: sandel-cs@ansell.com  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email: insidesalesus@ansell.com

**Latin America & Caribbean**  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email: luvas.medicas@ansell.com

**Europe, Middle East & Africa**  
Tel: +32 (0)2 528 74 00  
Email: info@ansell.eu

**Asia Pacific**  
Tel: +603 8310 6688  
Email: apac.medical@ansell.com  
Japan Tel: +813 5549 8151  
Japan Email:  
info.medical.jp@ansell.com

**Australia & New Zealand**  
Tel: +61 3 9270 7270  
Email: protection@ap.ansell.com