

# Blood collection devices and accessories

## Safety blood collection needles

### BD Vacutainer® Eclipse™ Blood Collection Needle

The BD Vacutainer® Eclipse™ safety needle for venous blood sampling has a fully integrated safety shield, which once activated, protects against needlestick injuries. This safety shield is an integral part of the needle and its orientation aligns with the needle bevel. The safety mechanism is designed for single-handed activation. The fully integrated safety shield engages over the needle with an audible click, irreversibly locking with a triple closure mechanism.



#### Further information

BD recommends for all needles without a pre-attached holder that a BD Vacutainer® One-Use Holder is used.

Clinical and technical information is available on request.



### BD Vacutainer® Eclipse™ Blood Collection Needles



Cat. no.	Size	Needle length	Colour code	Box/Case
368609	21 G (0.8 mm)	32 mm		48/480
368610	22 G (0.7 mm)	32 mm		48/480

### BD Vacutainer® Eclipse™ Blood Collection Needle with Pre-Attached Holder

With this safety needle, the holder is already pre-attached, so it is not necessary to manually assemble the needle and holder. This ready-for-use blood sample needle and holder is individually packaged in a sterile blister.



### BD Vacutainer® Eclipse™ Blood Collection Needles with Pre-Attached Holder

Cat. no.	Size	Needle length	Colour code	Box
368650	21 G (0.8 mm)	32 mm		100
368651	22 G (0.7 mm)	32 mm		100