

evercare® Visitor gown with buttons

REF: 2712-01, 2714-01, 2716-01

evercare® Visitor gown can be used as visitor and patient gown. The gown is made of the same material as Clean Air suits which means low linting and low bacterial penetration to support infection control. Comfortable and warming cuffs in the arms that are in full length. Closed with buttons to allow for easy dress on and dress off. Comes in 3 sizes. Single use and Unisex.

Key Features:

- ▶ Material properties meet EN13795 for Clean Air Suit
- Visitor gown can be used as well as a patient gown
- Made of SMS nonwoven for better protection
 - » low linting
 - » low bacterial penetration
- Comfortable cuffs in the arm sleeves
- Userfriendly buttons



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2712-01	Visitor gown with buttons, blue	L, length 120 cm	5/20
2714-01	Visitor gown with buttons, blue	XL, length 128 cm	5/20
2716-01	Visitor gown with buttons, blue	XXL, length 140 cm	5/20

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices













REACH regulation:

The product(s) covered by this data sheet does(do) not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation:

EN ISO 10993-5, EN ISO 10993-10

Compliance with legal requirements:

Medical Devices Directive 93/42/EEC, Class I

Compliance with product standards:

N/A

Quality standards:

EN ISO 13485

Label information and packaging compliance:

EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristic	
Body, sleeve, collar	SSMMS PP nonwoven, antistatic, blue	
Cuffs	PES, white	
Poppers	POM, pvc-free, white	
Neck label	99% nylon & 1% calcium carbonate, white	
Sewing thread	PES, white	

Warnings and Recommendations

Warnings:	The product is single use. If reused performance of product may deteriorate, cross contamination may occur.	
Storage Recommenda- tion	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.	
Disposal Recommendation	Dispose in accordance with local regulations.	
Shelf Life	5 Years	

Packing Information

Package	Material	Size (L x W x H) mm	
Single Pack	N/A	N/A	
Dispenser Pack	PE bag	400 x 300 x 360	
Transport Carton	Carton, TCF or ECF quality	380 x 297 x 325	
Packing Tape	Transparent, PVC-free	N/A	

Barcodes

REF	Single pack	Dispenser	Transport Caton
2712-01	N/A	6 438129 323405	6 438129 523409
2714-01	N/A	6 438129 323412	6 438129 523416
2716-01	N/A	6 438129 323429	6 438129 523423



OneMed Group Oy Address: Metsäläntie 20 Postal Address: P.O. Box 10 FI-00321 Helsinki, Finland Phone +358 20 786 6800 www.onemed.com Domicile Helsinki Business ID 2039640-1