

APPILIMED



Zone industrielle Route de Pra de Plan 1 CH-1618 CHATEL-ST-DENIS Switzerland TEL FAX E-mail Website : +41 (0) 21 948 92 74 : +41 (0) 21 948 90 36 : info@applimed.ch : www.applimed.ch

## **PRODUCT DATA**

I-Product Description & Classification							
I/1	Designation	Metal Applicator					
I/2	General Description	Applicator with Aluminum Stick, 150x0.9 mm, Cotton head					
I/3	GMDN N°	33954					
	CND N°	M010103					
I/4	Reference Number	1801183					
I/5	Class						
		According to Annex VIII of the Regulation 2017/745 (MDR), Rules 4 and 5					
I/6	Notify Body	CE					
I/7	Composition Raw	100% cotton (head) and Aluminum Stick					
	Material						
I/8	Sterile	🗌 Yes 🖾 No					
		lf, yes					
		ETO Gamma Steam					
I/9	First Packaging	Number of pieces in Packaging:         100 piece per pouch					
I/10	Picture						
		the second secon					
		and the state of a state					
		and the second s					
		a second s					

II- Description of use and operating of device						
II/1	Intended use/indication:	A device intended used to apply a substance (eg medication) to the Ear/Nose/throat (ENT) treatment and to care the superficial wound for absorption, barrier or compressing Alone, single use and transient use (<60min)				
II/2	User:	Medical staff or professional				
II/3	Contraindication of use	Not use, if the device is dirty or visual defects If packaging is damaged or open do not use.				
II/4	Special requirement for use	Should not be use with Alcohol and Octenisept				
II/5	Handling and storage conditions:	Standard conditions or atmospheric conditions				

III- Supplementary information						
III/1	Compatibility with others	n/a				
	medical devices (in the set)					
III/2	Device contents animal	NO				
	substances					
III/3	Device contents medicinal	NO				
	substances					
111/4	Device comes into contact with substances	Applimed has demonstrated compatibility with the following substances: - Water - Peroxide water (3%) - Betadine - Eosine - Merfen				
III/5	Used with Accessories	NO				
III/6	Used in combination with	NO				
	other Medical Device					







Zone industrielle Route de Pra de Plan 1 CH-1618 CHATEL-ST-DENIS Switzerland TEL FAX E-mail Website : +41 (0) 21 948 92 74 : +41 (0) 21 948 90 36 : info@applimed.ch : www.applimed.ch

	IV- MANUFACTURING								
IV/1	Manufacturing Operations	<ul> <li>□ peel pack</li> <li>➢ pouches/bags</li> <li>□ boxes</li> <li>⊠ cartons</li> <li>□ cleaning</li> <li>□ sterilization</li> </ul>	100 pieces per pouch 100 boxes per carton						
IV/2		Example of a la <b>REF</b> 1801183 [ Code <b>18-1183</b> <b>10.000</b> <b>Applicateur aluminium, 150x(</b> Applicator aluminium, tête coton Applicator aluminium, tête coton Applicator aluminium, baumwollkopf Aplicador de aluminio, cabeza de algodón Applikator Aluminium, bomuldshoved Applicator e in alluminio, testa di cotone Applicator alumininia, puuvilla kärki Applicator aluminium, katoenen tip <b>CEE MD</b> Applimed SA Pra de Plan 1 <b>EC REP</b>	Lot 22.04.01 2023-04 Stk/pieces 0.9mm algodão Viena AS 0 Stenløse, Danmark						
IV/3	Place of sterilisation:	n/a							